

Design Guide for selected DEDICAM[®] restorations and attachments with exocad[®] DentalCAD

November 2023

M-1733-PRT-EN-INT-BHCL-00-112023

Content



P. 4

Edge design of thin tapered margin line on abutments



P. 9

Design of threaded holes M1.4 on one-piece titanium abutments



P. 21

Design of one-piece zirconia abutments for CERALOG® Hexalobe implants




P. 134

Design of healing abutments



P. 52

Design of anti-rotation protections on abutments



P. 69

Design of screw channels with the aim of an attachment



P. 95

Design of primary parts for telescopic crowns



P. 156

Design of angulated screw channels

Inhalt



P. 31

Design eines MK1 Attachments an Brücken und Kronenblöcke



P. 79

Recommendation on cross-section design for bars



P. 84

Adding a MK1 cantilever attachment to a bar



P. 101

Attaching a Preci-Vertex® with interlock and circumference to crowns and bridges





P. 150

Design of printed models

Edge design of thinly tapered margin line on abutments

Edge design of thinly tapered margin line on abutments



Problem:


Thin margin line can negatively influence the milling result (partial chipping at the edges) and leads to redesign and remake of the structure as a **possible cause of delivery delays**.

Edge design of thinly tapered margin line on abutments

Note: specific edge reinforcement

Abutment design

- 1 → Activate Button „Advanced“
- 2 → Profile borders: Set minimum height 0.1mm (often default height is 0.2mm which may be used as well)




Edge design of thinly tapered margin line on abutments

Note: specific edge reinforcement


Edge reinforcement in sectional view

- Margin line shows a specific edge reinforcement of approx. 0.1 – 0.2mm and can be milled without loss of cervical contour.




Edge design of thinly tapered margin line on abutments

Edge reinforcement = 0.1mm



Edge reinforcement = 0.2mm



Design of threaded holes M1.4 on one-piece titanium abutments for horizontal screwed crowns, fixed with the “Bredent screw”


Design of threaded holes M1.4 on one-piece titanium abutments

Design abutment:

Note:

Wax up or virtual tooth makes it easier to create individual tooth

Minimum wall thickness for coping or crown where the threaded hole is planned: **0.9mm**



Design of threaded holes M1.4 on one-piece titanium abutments


Design abutment and create horizontal screw connection:

Abutment design finished

Wizard – Abutment Design:

- Top
- Bottoms
- Advanced


Press „Next“



Design of threaded holes M1.4 on one-piece titanium abutments

Design abutment and create horizontal screw connection:

Abutment design finished



Wizard – Free-forming: Att.

Check mark „Remove“

Select „Library“

Select attachment:

DEDICAM M1.4 (for Bredent screw)

Design of threaded holes M1.4 on one-piece titanium abutments

Design abutment and create horizontal screw connection:

Align attachment position:

View direction on planned horizontal screw connection

Important note:


Attachment is not displayed authentically. This is necessary and offers qualitative advantages

Selected attachment:

DEDICAM M1.4 (for Bredent screw)

Insertion direction

View and Surface




Design of threaded holes M1.4 on one-piece titanium abutments

Design abutment and create horizontal screw connection :

Positioning of the attachment:

Check mark

- Move
- Rotate
- Position (1) meets functional and aesthetic aspects




Design of threaded holes M1.4 on one-piece titanium abutments

Design abutment and create horizontal screw connection:

Attachment fine adjustment

Pay attention to the groove on the attachment

- Angle (2)
- Depth(3)



Design of threaded holes M1.4 on one-piece titanium abutments


Design abutment and create horizontal screw connection:

Attachment fine adjustment:

View direction on planned horizontal screw connection

Note:

Wax up or virtual tooth makes it easier to create individual tooth.



Design of threaded holes M1.4 on one-piece titanium abutments


Design abutment and create horizontal screw connection:

DEDICAM Attachment „M1.4 for Bredent Screw“ :

Grooves for the independent use of Bredent Screw M1.4 (according to IFU)

Grooves on the DEDICAM Attachment:

- Groove must be enclosed from abutment (1)
- Conical part of screw head 0.3mm enclosed from abutment (2)
- Groove as border until the screw might be reduced. max. 1.3mm (3)
- Total length of screw (4)



Design of threaded holes M1.4 on one-piece titanium abutments


Design abutment and create horizontal screw connection:

Control of correct Positioning from DEDICAM Attachment:

- Show screw channel (1)
- Threaded hole with at least 1mm depth into the abutment (2)
- Show min. thickness and attached parts (3)

Attention: Threaded hole attachment must be placed above abutment min. thickness and abutment screw.

If possible threaded hole might not be placed into abutment screw channel.



Design of threaded holes M1.4 on one-piece titanium abutments

Abutment with horizontal screw connection for Bredent screw M1.4

Finished abutment with the horizontal screw design

Note:
File-splitting is not supported

exocad

Design of one-piece zirconia abutments for CERALOG® Hexalobe implants



Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Order creation:

1. Select tooth position
2. Select indication

Recommendation:

Select „Coping“ or „Offset coping“



Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Order form:

3. Select material for coping or offset-coping


4. *Options & Parameters*

4.a Define type of scan

4.b *Implant-based?:* „On Custom Abutment“

4.c Abutment material: „ZrO2 CERALOG One-Piece Abutment“

5. Press „OK“ to confirm




Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Select CAD library:

6. CAD-Library: „CERALOG® Hexalobe one-piece abutment DEDICAM“

Attention: Ensure you took the correct library


7. Press „Best fit matching“ and confirm with „Next“



Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Wizard Generate Abutment Bottoms:

1. Visualize and respect outer blank limit for design. Hint: Blank size differs from Titanium blank
2. Under „Advanced“ set profile borders value „Height“ to 0.11mm
3. Set „Radius“ always to 0mm (exocad versions provided by AmannGirrbach have mostly 0.2mm on default)




Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Wizard Generate Abutment Bottoms:

4. Never check mark „Allow below implant“
5. Check mark „Free-form“ to finalize bottom

Recommendation: Respect outer blank limit – keep blank half transparent during design




Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Wizard Abutment Design above margin „TOP“

6. Abutment design above margin
7. Do not violate outer blank or abutment angulation of 30°
8. Drag „green arrow“ or set another insertion axis to correct

Recommendation: Respect outer blank limit – keep blank half transparent during design




Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Blank visualization informs about maximum outer shell, angulation limit and minimum thickness

- 9. Screw channel
- 10. Angulation limit of 30°
- 11. Minimum thickness around abutment screw

Please note: Designs which do not meet the requirements will be rejected.




Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Wizard Free-forming:
„Smooth/Flattern“

12. Finalize abutment surface
with „Smooth/Flattern“

13. Design finalized → Press
„Next“

Note: Avoid sharp edges and
corners on the abutment



Design of one-piece zirconia abutments for CERALOG® Hexalobe implants

Abutment design: Final design is always saved incl. screw channel

14. Screw channel

Edge rounding on screw channel


Screw channel edges might be slightly adapted on production site if needed to prevent cracks.

15. Design finished

16. Press „Next“ and send design incl. *constructioninfo.xml* to Camlog

IMPORTANT

Visualization of the hexalobe-connection is distorted and for safety reasons only millable by Camlog



Adding a MK1 attachment to a bridge or crown block

Adding a MK1 attachment to a bridge or crown block

Note:


In order to position attachments to fixed bridges or crown blocks or to cut them by the gingiva. Note the explanations with the example designs when creating the order.



For the design of attachments it is necessary to use the DEDICAM[®] CAD library.

Adding a MK1 attachment to a bridge or crown block

Job definition Multitide mode



■ Coping

Actions

- Scan
- Design
- Manufacture
- Model Creator
- dentalshare
- Copy
- Print

exocad

6136

Tooth shades: A3
Scan mode: One stone model only

Example: Blocked copings on tooth 43 + 44 with MK1 distal on tooth 44



Adding a MK1 attachment to a bridge or crown block

The screenshot shows the 'Tooth 43' configuration window. On the left, a category menu lists various dental components. The 'Coping' option is highlighted with a red box. The main area displays a 'Material' list with 'CoCr Typ4 attachment combination' and 'Ti6Al4V attachment combination' highlighted in red boxes. To the right, the 'Options & parameters' panel shows settings for scanning and printing, with a red box around the 'OK' button at the bottom right.

← Tooth 43 Material configuration (local): Default

Crowns/Copings

- Anatomic crown
- Coping**
- Pressed crown
- Offset coping
- Provisional crown
- Preform crown

Pontics

- Anatomic pontic
- Reduced pontic
- Pressed pontic
- Provisional pontic

Inlays, onlays and veneers

- Inlay/Onlay
- Offset inlay
- Veneer

Digital copy milling

- Anatomic waxup
- Reduced waxup
- Pontic waxup

Primary units

- Bar pillar
- Bar segment
- Attachment
- Telescopic crown

Bite splint

- Bite splint
- Bite splint (missing tooth)

Residual dentition

- Antagonist
- Adjacent tooth
- Missing tooth

Material

5-Axis / Laser / 3D Print

- CoCr Typ4 attachment combination**
- CoCr Typ4 Primary Part Telescope
- Ti6Al4V attachment combination**
- Ti6Al4V Primary Part Telescope
- IPS e.max ZrCAD MO/LT 4-16 Elements

Options & parameters

Scan a pre-op model? No

Scan gingiva separately? No

Implant based? No implant

Shrink mm

Minimal thickness mm

Gap thickness of cement mm

▶ ADVANCED PARAMETERS

Clear **OK**




Order creation: Tooth 43

- *Coping*
- Material: *CoCr Typ4 attachment combination*
or
- *Ti6Al4V attachment combination*
- *OK to confirm*

Adding a MK1 attachment to a bridge or crown block

Job definition Multitide mode



Actions

- Scan
- Design
- Manufacture
- Model Creator
- dentalshare
- Copy
- Print

exocad

6136

Tooth shades: A3
Scan mode: One stone model only



Order creation: Tooth 44

- *Primary telescope*

Note:

MK1 Attachment will be placed here

Adding a MK1 attachment to a bridge or crown block

← Tooth 44 Material configuration (local): Default

Crowns/Copings

- Anatomic crown
- Coping
- Pressed crown
- Offset coping
- Provisional crown
- Preform crown

Pontics

- Anatomic pontic
- Reduced pontic
- Pressed pontic
- Provisional pontic

Inlays, onlays and veneers

- Inlay/Onlay
- Offset inlay
- Veneer

Digital copy milling

- Anatomic waxup
- Reduced waxup
- Pontic waxup

Primary units

- Bar pillar
- Bar segment
- Attachment
- Telescopic crown**

Bite splint

- Bite splint
- Bite splint (missing tooth)

Residual dentition

- Antagonist
- Adjacent tooth
- Missing tooth

Material

5-Axis / Laser / 3D Print

- CoCr Typ4 attachment combination**
- CoCr Typ4 Primary Part Telescope
- Ti6Al4V attachment combination**
- Ti6Al4V Primary Part Telescope

Options & parameters

Scan a pre-op model? No

Scan gingiva separately? No

Implant based? No implant

Minimal thickness 0.5 mm

Gap thickness of cement 0.06 mm

▶ ADVANCED PARAMETERS

Clear OK

Order creation: Tooth 44

- *Coping*
- *Material: CoCr Typ4 attachment combination*
or
- *Ti6Al4V attachment combination*
- *OK to confirm*

Adding a MK1 attachment to a bridge or crown block

Fitting parameter 43 und 44: check values

Note: Parameter should be identical on all stumps → press “Next”

Important: „Don't block out undercuts“ may not be marked


The screenshot displays the exocad software interface for the 'Wizard' of 'Crown Bottoms'. The interface is divided into several panels:

- Left Panel:** Contains 'CEMENT GAP' settings. The 'No cement gap' option is selected. A 'Gap' of 0.06 mm is set. Below, 'SELECT ZONES BY DISTANCE' is set to 'From margin' with a value of 1 mm. A red box highlights these settings.
- Middle Panel:** 'CROWN BORDER PARAMETERS' with five settings: 1. Horizontal (0.1 mm), 2. Angled (0.1 mm), 3. Angle (65°), 4. Vertical (0 mm), 5. Below margin (0 mm). A red box highlights the 'BORDER' tab.
- Right Panel:** 'UNDERCUTS' settings. The 'Don't block out undercuts' checkbox is checked. A red box highlights this checkbox. Below are 'MILLING' settings: 'Anticipate milling' (checked), 'Diameter' (1.2 mm), 'Bullnose/flat tool' (unchecked), and 'Tool tip flat percentage' (80%).
- Bottom Panel:** Navigation buttons: 'Back' and 'Next' (highlighted with a red box) for each of the three panels. A tip below the buttons reads: 'Click "Next" to load default tooth models'.

The central 3D model shows a tooth with two yellow attachments. Green arrows point to the attachment areas, and a red arrow points to the 'Don't block out undercuts' checkbox. The exocad logo is visible in the bottom right corner.


Adding a MK1 attachment to a bridge or crown block

Adapt anatomy design to clinical situation: to reduce full anatomy → press “Next”



Adding a MK1 attachment to a bridge or crown block


Reduced anatomy design on tooth 43: for “Insertion axis telescope” → press “Next”



Adding a MK1 attachment to a bridge or crown block

Define telescope insertion axis: View direction = insertion direction


MK1 Attachment follows defined insertion axis



Adding a MK1 attachment to a bridge or crown block


Primary telescope: edit parallel surfaces –

Flat distal surface remaining for MK1 Attachment




Adding a MK1 attachment to a bridge or crown block

Primary telescope: Add control points to edit distal surface for the MK1 Attachment



Adding a MK1 attachment to a bridge or crown block

Primary telescope: Add control point to edit distal surface for the MK1 Attachment



Adding a MK1 attachment to a bridge or crown block

Primary telescope: Add control point to edit distal surface for the MK1 Attachment

Occlusal extension of added control points

Control points may be located very close to each other

Wizard Primary Telescope

TELESCOPE PARAMETERS

Angle

Diam (C)

STICK (BOTTOM)

Crown Gingiva None

Visualize areas with friction

ADANCED ANGLE CONTROL...

INSERTION AXIS

Undo Redo

Back Next

Speichern

Experte

Werkzeuge

Abstände zeigen

TruSmile

Neue Ansicht

Ansicht 1

Ansicht 2

Ansicht 3

Ansicht 4

Ansicht 5

exocad

9675/64

Adding a MK1 attachment to a bridge or crown block

Primary telescope: Add Gripper to edit distal surface for the MK1 Attachment

Drag in distal direction indicated by arrows

Hint: visualize full anatomy for perfect editing of the connection surface to the MK1 Attachment

exocad

5675/64

Tree View:

- EN/AUSBLENDEN
- Kiefer-Scandaten
- Alle anderen Teile
- Anatomie
 - Anatom. Pontic (45)
 - Vollkronen (43)
 - Vollkronen (44)
- Geschlumpft
- Teleskope
- ZÄHNE
 - 43: Anat. Kappchen
 - 44: Primärteleskop
 - 45: Attachment
- VERSTECKT

Wizard Primary Telescope:

TELESCOPE PARAMETERS

Angle: 0°

Diam (C): 1.2 mm

STICK (BOTTOM): Crown, Gingiva, None

Visualize areas with friction

ADVANCED ANGLE CONTROL...

INSERTION AXIS

Undo Redo

Back Next

Toolbar:

- Speichern
- Experte
- Werkzeuge
- Abstände zeigen
- Tru
- TruSmile
- Neue Ansicht
- Ansicht 1
- Ansicht 2
- Ansicht 3
- Ansicht 4
- Ansicht 5
- Navigation icons

Adding a MK1 attachment to a bridge or crown block

Primary telescope: edit veneering surface → afterwards press “Next”

Clearly defined metal/veneering material surface demarcation

Remaining control points not used for flat MK1 surface will be adapted to suit coping design

Wizard Primary Telescope

TELESCOPE PARAMETERS

Angle: 0

Diam (C): 1.2 mm

STICK (BOTTOM): Crown Gingiva None

Visualize areas with friction

Undo Redo

Back Next

Speichern

Experte

Werkzeuge

Abstände zeigen

TruSmile

Neue Ansicht

Ansicht 1

Ansicht 2

Ansicht 3

Ansicht 4

Ansicht 5

exocad

5675/64

Adding a MK1 attachment to a bridge or crown block

Telescope free-form: Mark DEDICAM MK1 Attachment from the library – select DEDICAM MK1 “cut to gingiva”

1. Mark + Add

2.

Attention:

1. Library: DEDICAM MK1 Attachment
2. Typ: DEDICAM MK1 – „cut to gingiva“

Important: select valid quadrant 1 + 3 or 2 + 4


exocad

Adding a MK1 attachment to a bridge or crown block

Telescope free-form: Attachment – DEDICAM MK1 “cut on gingiva” positioning


3. Place MK1 Attachment by „Move” and „Rotate”

4. Attention: „Scale” must be disabled. If not, current software version is unsuitable - update is highly recommended.



Adding a MK1 attachment to a bridge or crown block

Telescope free-form: Attachment – DEDICAM MK1 “cut on gingiva” adapt to gingiva → press “Next”



Adding a MK1 attachment to a bridge or crown block

Frame design 43 and 44 blocked through connector → press “Next”

7. Connector shape = Shape definition

8. Connector Free = connector shape follows „Gripper“ position

Wizard connectors

SHAPE | FREE

CROSS SECTION

Set connector cross section

Set connector height and width

Min. Area mm²

Visualize connector thickness

SHAPES

Apply cross section / shape change

Move by mouse • Deform with Shift key
Move docking point with Ctrl key

Back Next

Click 'Next' to merge parts to restoration.

Wizard connectors

SHAPE | FREE

CONNECTOR CONTROL POINTS

Move point
Ctrl+click to add new control point
Move in/out
Move freely
Shift+click: Move all

CONNECTOR THICKNESS

Very thick 0% 50% 100% Just okay Too thin

Undo Redo


Back Next

Click 'Next' to merge parts to restoration.

4.49
3.22
9.92

Adding a MK1 attachment to a bridge or crown block

Finalize and save the design: minor adjustments possible through free-form tool → press “Next”



Design of anti-rotation protection on abutments


Design of anti-rotation protection on abutments

Advantage: Efficient and easy to use

1. Go to „TOP“ on Wizard Abutment Design

Note:

Number of Control points is dependend on the abutment design



Design of anti-rotation protection on abutments

Advantage: Efficient and easy to use

2. Add new control point in the middle of two existing control points

Recommendation:


Add control point on the approximal side

Procedure:

- Ctrl + left mouse button
- New control point added
- Positioning between two existing control points

Note:

Distance between the existing control points is about 2.0 to max. 3.0mm




Design of anti-rotation protection on abutments

Advantage: Efficient and easy to use

3. Select from „Tools“ → „Measurement tool“:
Measure the distance from preparation margin to orange colored spline (Mouspointer – left mouse button)

Recommendation:

Shoulder width from abutment shoulder:
shoulder: lower control point
respectively orange dashed line
approx. 0.8mm




Design of anti-rotation protection on abutments

Advantage: Efficient and easy to use

4. Extend shoulder width to approx. 1.3mm

Procedure:

Drag control point's blue arrow from approx. 0.8mm to approx. 1.3mm inwards




Design of anti-rotation protection on abutments

Advantage: Efficient and easy to use

5. Select from „Tools“ → „Measurement tool“ to proof the shoulder width on anti-rotation
6. Measure the distance from preparation margin to orange colored spline (Mouspointer – left mouse button)

Recommendation:

Shoulder width from abutment shoulder as an anti-rotation protection: approx. 1.3mm




Design of anti-rotation protection on abutments

Option File-Splitting

After finalization of the abutment design, the CAD software offers to design a crown (File-splitting)

Note:
DEDICAM allows file-splitting only for single tooth and bridge restoration up to three units




Design of anti-rotation protection on abutments

Option File-Splitting

Notes:


- Check parameter e.g.
 - Cement gap
 - Drill diameter
- Check material related drill radius on DEDICAM Software-Parameterchart



Design of anti-rotation protection on abutments

Option File-Splitting

Coping / crown design



Design of anti-rotation protection on abutments

Advantage: Efficient and easy to use


Suitable anti-rotation of coping / crown follows abutment design

→ Cement gap is even

Note:

The abutment shoulder has no cement gap approx. 1mm from preparation margin.

All other surfaces have approx. 60 to 80 µm cement gap.



Design of anti-rotation protection on abutments

Alternative to the before mentioned anti-rotation protection


After the initial design of the abutment the anti-rotational is installed in the „Free-Forming – Attachments“ mode.

The following selection is available:

- Attachments: Subtract
- Library:
 - Parametric design (Rotation)
 - BoltDesign

Note:


- Parameter adaption required
 - 1** → Diameter: minimum 1.5mm
 - 2** → Height: minimum 5 - 6mm (longer as upper abutment part)
 - 3** → Diameter (Head): reduce e.g. 1.7mm (easier positioning)



Design of anti-rotation protection on abutments

Alternative to the before mentioned anti-rotation protection

- Position the abutment on occlusal view such as lateral abutment flanks are visible to define the insertion direction of the anti-rotational
- Insertion direction **1**
→ View
- Place the attachment on the desired spot **2**



Design of anti-rotation protection on abutments

Alternative to the before mentioned anti-rotation protection


- Place the attachment on the desired spot **2**
- 2D cross section – slightly above the margin: Check correct and millable positioning **3**

Important note:

- Place the attachment (Bolt design \varnothing min. 1.5mm) less than the maximum diameter into the abutment

Diameter = 1.5mm



maximum diameter on abutment =
1.43mm



Design of anti-rotation protection on abutments

Alternative to the before mentioned anti-rotation protection

Verify the correct position of the attachment before applying it to the design **4**



Design of anti-rotation protection on abutments

Alternative to the before mentioned anti-rotation protection



After subtraction sharp edges must be smoothed in „Free-Forming“; „Free-Smooth/Flattern“ **5**

→ Activate button „Smooth/Flattern“ **6**

Suggestion:

- Put „Strength“ on a low level **7**
- Put „Brush size“ on a low level **8**
- All sharp edges will be smoothed **8** with this tool

Important note: Due to the drills used in the production „the supplied anti-rotation protection on the structure may deviate minimally.



Design of anti-rotation protection on abutments


Alternative to the before mentioned anti-rotation protection

Check the anti-rotational width after smoothing by using 2D cross section again

Width must be **approx. 1.7mm**

Important note:

Due to the drills used in the production the supplied anti-rotation protection on the structure may deviate minimally.



Design of anti-rotation protection on abutments


This anti-rotational protection is suitable for a precisely fitting coping / crown whether designed in file-splitting or via second design.

Note:

Orders in file-splitting are only processed with anti-rotation protection ensuring that the milling and the fit of the coping / crown can be guaranteed.

Under consideration of the minimum wall thickness this anti-rotational protection is recommended for the following abutment types:

- Meso structures for titanium bases CAD/CAM, for crowns
- One-piece abutments made of titanium
- One-piece abutments for CERALOG Hexalobe implants



Design of screw channels with the aim of an attachment


Design of screw channels with the aim of an attachment

Restrictions:

- Implant parts scanned as dies
- Not available for IPS e.max CAD
- Required software version: 2.2 Valletta or higher


Design of screw channels with the aim of an attachment

Important: Attachment selection and application must be done at last design step after switching to expert mode!






Design of screw channels with the aim of an attachment

- 1 After calling up expert mode restoration has to be merged and saved first
- 2 Now functionality «Free-form merged reconstruction» is available



Design of screw channels with the aim of an attachment



- 1 Select tab «Att.» and «Remove»
- 2 «OK» to confirm possible warning message
- 3 Two sizes of screw channels are offered
- 4 Important: Mark «Allow any changes»



Design of screw channels with the aim of an attachment

- 1 Align the attachment by its tip looking into the cavity
- 2 Do **not** use «Scale» tool
- 3 Check attachment to implant axis orientation by sectional tool tool


Important: For trouble-free integration of the screw ensure the attachments axis is aligned with the implant axis



Design of screw channels with the aim of an attachment

Divergences between implant axis and screw channel axis may prevent the screw from receiving its thread.



Camlog will not carry out any verification about functionality prior and post to manufacturing.



Design of screw channels with the aim of an attachment



Attention: Attachment should not touch crown margin!

If necessary adapt attachment length and / or axis






Design of screw channels with the aim of an attachment

Correctly placed attachment to create a screw channel





Design of screw channels with the aim of an attachment

- 1 To apply correct placed attachment press «Apply» and «OK»
- 2 Confirm free-form changes again
- 3 Added crew channel on final restoration




Recommendation on cross-section designs for bars

Recommendation on cross-section designs for bars



Note:

Preci-Horix cross-section 4.6mm^2 plus cross section below Preci-Horix geometry



Measure of bar height and width below Preci-Horix geometry

e.g.


$$2.45 \times 0.67 = 1.64\text{mm}^2$$

$$\text{Preci-Horix} + 4.60\text{mm}^2$$

$$\text{Total} = 6.24\text{mm}^2$$



For bars with cross-section-reducing attachments (Preci-Horix), it should be noted that the reduced bar cross-section must be compensated by a larger width or height.

Recommendation on cross-section designs for bars



Note:

Preci-Horix cross-section 4.6mm^2 plus cross section below Preci-Horix geometry



Measure of bar height and width below Preci-Horix geometry

e.g.


$$2.42 \times 1.84 = 4.45\text{mm}^2$$

$$\text{Preci-Horix} + 4.60\text{mm}^2$$

$$\text{Total} = 9.05\text{mm}^2$$


For bars with cross-section-reducing attachments (Preci-Horix), it should be noted that the reduced bar cross-section must be compensated by a larger width or height.

Recommendation on cross-section designs for bars



Recommended bar width of min. 2.2mm ensures manufacturing of suprastructures

Recommendation on cross-section designs for bars



It is recommended not to go below a cross-section of approx. 8 to 9mm².

Note: The cross-section is smaller on bars for prefabricated bar matrices (Micro/Macro Dolder).


Adding a MK1 cantilever attachment to a bar

Adding a MK1 cantilever attachment to a bar

Bar type: DEDICAM
Primary Bar 2° tapered 1

Cantilever bar width: min. 2.8mm

Note:
It is the recommendation of the MK1 attachment manufacturer to use a 2° tapered bar type.




Adding a MK1 cantilever attachment to a bar

Select DEDICAM MK1 attachment according its placing region 1. +.3. Q. or 2.+ .4. Q.

1


Note:

In case the attachment shall be cut to the gingiva, select MK1 attachment with „cut-to-gingiva“ naming.




Adding a MK1 cantilever attachment to a bar

Push the MK1 attachment into the bar profile. **1**




Adding a MK1 cantilever attachment to a bar

Do not push the MK1 attachment further into the bar profile than the marked line indicates. **1**



Adding a MK1 cantilever attachment to a bar

Change to side view for positioning of the MK1 attachment in terms of height. Therefore, use the green arrow. **1**




Adding a MK1 cantilever attachment to a bar

Attention:

According to the IFU, the maximum extension of the bar including attachment is 10mm measured from the outer diameter of the implant.

Use the 2D cross section to verify the length.




Adding a MK1 cantilever attachment to a bar

The ideal solution is to place the functional part of the MK1 attachment in the interdental space of the denture teeth. 1

Note:


In order to ensure operability by the patient it is recommended to place the functional part of the MK1 attachment not further distal than the 2nd premolar.



Adding a MK1 cantilever attachment to a bar

After final positioning of the MK1 attachments activate the button „Snap to Gingiva“ and apply by clicking „Next➔“


1



1


Adding a MK1 cantilever attachment to a bar

Do not use any other tools from the „Free-Form Bar“ to finalize the bar design.



Adding a MK1 cantilever attachment to a bar

Send the design via
Dentalshare or our eService
to Camlog.



Design of a primary part of a telescopic crown

Design of a primary part of a telescopic crown

Stump parameter, internal fit for primary components

Gap

Cement Gap 0.06

Additional spacing

Axial 0.02

Radial 0.00

Border

Horizontal 0.10

Angled 0.10

Angel 65°

Vertical 0.00

Below margin 0.00

Undercuts


Angle 0.00

Size 0.10

Milling

Anticipate milling → check mark

Diameter 1.2



Design of a primary part of a telescopic crown

Stump parameter, internal fit for primary components

- CoCr Typ4 Primary Part Telescope
- Ti6Al4V Primary Part Telescope

Angle: 0.00° - 6.00°

- Value can be changed: telescope = 0° / double crown 2 – 6° (Note: use same value per jaw)

Minimum thickness: 0.50mm

- Value should not be changed if possible in order to ensure that there is enough material thickness even after corrections.

Border: H: 0.1mm; A: 0.1mm; A: 65°; V: 0.0mm

- For accuracy of fit use optimised CAD parameter.

Milling: Ø 1.2mm

- „Anticipate milling“ → check mark
Diameter 1.2mm

Design of a primary part of a telescopic crown

Order creation


Select **1.** Material configuration (local)

2. DEDICAM 2018-06-20 or higher and

click „OK“ **3.** to confirm

4. Indication „Telescope crown“

5. Material „CoCr Typ4 Primary Part Telescope“ or „Ti6Al4V Primary Part Telescope“



Design of a primary part of a telescopic crown

Design primary components as follows:

1. Select „Uniform insertion direction for jaw“

2. Telescope angle = 0° or double crown angle 2-6°

3. \emptyset bur = 1.2mm

4. Height function surface = Height ~5.0mm

5. abutment shoulder = epi- or supragingival





Design of a primary part of a telescopic crown

Finish of the primary components design

Right quadrant = primary telescope

Left quadrant = primary double crown 2°

Note: parallel and conical walled telescopes must never be mixed.



Attaching a Preci-Vertix® with interlock and circumference to crowns and bridges

Attaching a Preci-Vertex® with interlock and circumference

Note:

In order to position attachments to fixed bridges or crown blocks or to cut them by the gingiva. Note the explanations with the example designs when creating the order.




For the design of attachments it is necessary to use the DEDICAM® CAD library.

Attaching a Preci-Vertex® with interlock and circumference

Job definition

Multidie mode



Tooth shades

A2

Scan mode

One stone model only

Actions



- Scan
- Design
- Manufacture
- Model Creator
- dentalshare
- Copy
- Print

exocad

Example: Copings blocked on tooth 43 + 44 with a Preci-Vertex distal to 44 and an Interlock between 43 + 44



Attaching a Preci-Vertex® with interlock and circumference



Order creation: Tooth 43

- Primary telescope
- Material: CoCr type 4 or Ti6Al4V attachment combination
- Press OK to confirm

Attaching a Preci-Vertex® with interlock and circumference

Definition



Aktion

- Scannen
- CAD
- CAM
- Model Creator
- dentalshare
- Kopieren

exocad

5675



Order creation: Tooth 44

- Primary telescope

Note:

Preci-Vertex attachment and Interlock will be placed here.

Attaching a Preci-Vertex® with interlock and circumference

← Tooth 44 Material configuration (local): Default

Crowns/Copings

- Anatomic crown
- Coping
- Pressed crown
- Offset coping
- Provisional crown
- Preform crown

Pontics

- Anatomic pontic
- Reduced pontic
- Pressed pontic
- Provisional pontic

Inlays, onlays and veneers

- Inlay/Onlay
- Offset inlay
- Veneer

Digital copy milling

- Anatomic waxup
- Reduced waxup
- Pontic waxup

Primary units

- Bar pillar
- Bar segment
- Attachment
- Telescopic crown

Bite splint

- Bite splint
- Bite splint (missing tooth)

Residual dentition

- Antagonist
- Adjacent tooth
- Missing tooth

Material

Subaxis / Laser / 3D Print

- CoCr Typ4 attachment combination
- CoCr Typ4 Primary Part Telescope
- Ti6Al4V attachment combination
- Ti6Al4V Primary Part Telescope

Options & parameters

Scan a pre-op model? No

Scan gingiva separately? No

Implant based? No implant

Minimal thickness 0.4 mm

Gap thickness of cement 0.06 mm

▶ ADVANCED PARAMETERS

Order creation: Tooth 44

- Primary telescope
- Material : CoCr Typ 4 or Ti6Al4V attachment combination
- Press OK to confirm


Clear OK



Attaching a Preci-Vertex® with interlock and circumference

Job definition

Multidie mode



Telescopic crown

Tooth shades
A2

Scan mode
One stone model only

Click to connect

Actions

- Scan
- Design
- Manufacture
- Model Creator
- dentalshare
- Copy
- Print

exocad




Order creation: Tooth 43 + 44

- Block selected teeth

Attaching a Preci-Vertex® with interlock and circumference

Scan orientation: View direction = Insertion direction

Direction followed by Preci-Vertex circumference and Interlock



Attaching a Preci-Vertex® with interlock and circumference

Fitting parameter 43 und 44: check values

Note: Parameter should be identical on all stumps → press “Next”

Important: „block out undercuts“ may not be marked

The screenshot displays the 'Wizard Crown Bottoms' dialog box in the exocad software. The 'GAP' tab is active, showing parameters for 'CEMENT GAP' (set to 0.06 mm) and 'CROWN BORDER PARAMETERS' (1. Horizontal: 0.1 mm, 2. Angled: 0.1 mm, 3. Angle: 65°, 4. Vertical: 0 mm, 5. Below margin: 0 mm). The 'Don't block out undercuts' checkbox is checked. The 'Next' button is highlighted with a red box. A 3D model of a dental arch is visible in the background.

Wizard Crown Bottoms

GAP BORDER UNDERCUTS

CEMENT GAP

No cement gap

Gap 0.06 mm

Add another zone...

SELECT ZONES BY DISTANCE

From margin 1 mm

PARAMETER EXPLANATION

1. Horizontal 0.1 mm

2. Angled 0.1 mm

3. Angle 65°

4. Vertical 0 mm

5. Below margin 0 mm

Don't block out undercuts

Don't block out zone near prepline:

Size 0.1 mm

MILLING

Anticipate milling

Diameter 1.2 mm

Bullnose/flat tool

Tool tip flat percentage 80 %

Show undercuts Apply

Back Next

Back Next

Back Next


exocad

Add custom view

6136/64

Attaching a Preci-Vertex® with interlock and circumference

Adapt anatomy design to clinical situation → press “Next”




Attaching a Preci-Vertex® with interlock and circumference

Telescope insertion direction 43 and 44 → press “Next”

The screenshot displays the exocad software interface. The main view shows a 3D dental model of a maxillary arch with two green arrows indicating the insertion direction for teeth 43 and 44, both at a 2.6° angle. A wizard dialog titled 'Telescope Insertion Direction' is open on the left. It shows a tooth chart with teeth 43 and 44 highlighted in blue. The wizard has a 'Next' button highlighted with a red box. The 'Back' button is also visible. The exocad logo is in the bottom right corner of the software interface.

Attaching a Preci-Vertex® with interlock and circumference

Full anatomy design 43 + 44 will be reduced to primary crowns: further “Telescope” design
→ press “Next”



Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: edit parallel surfaces – modification as coping. Add control point if needed.


Added control points

- Add control point – Mousepointer on the spline
- ctrl + left mouse button
- Control point will be added

Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: edit parallel surfaces – modification as coping. Add "Gripper" if needed.

– Narrow interdental space between 43 and 44



Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: palatal surface preparation for shoulder

➤ Add control point – Mousepointer on the spline
➤ ctrl + left mouse button
➤ Control point will be added

Added control point

Wizard Primary Telescope

TELESCOPE PARAMETERS

Angle

Diam (C)

STICK (BOTTOM)

Crown Gingiva None

Visualize areas with friction

ADVANCED ANGLE CONTROL...

INSERTION AXIS

Undo Redo

Back Next

Save Expert Tools Instances TruSmile


Add custom view

0136/54

Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Shoulder adapted – closed interdental space between 43 and 44

Hint: visualize full anatomy for perfect editing



Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Add Gripper for editing distal surface of the Preci-Vertex attachment

➤ Add control point – Mousepointer on the spline
➤ ctrl + left mouse button
➤ Control point will be added

Information: Control points may be located very close to each other

Wizard Primary Telescope

TELESCOPE PARAMETERS

Angle

Diam. (C) mm

STICK (BOTTOM)

Visualize areas with friction

➤ ADVANCED ANGLE CONTROL...

➤ INSERTION AXIS

Undo Redo

Back Next

exocad

Save
Expert
Tools
Show distances
TruSmile

Add custom view

01/10/24

Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Add Gripper for editing distal surface of the Preci-Vertex attachment

visualize full anatomy for perfect editing

Occlusal extension of added control point

Wizard Primary Telescope

TELESCOPE PARAMETERS

Angle

Diam (C)

STICK (BOTTOM)

Crown Gingiva None

Visualize areas with friction

ADVANCED ANGLE CONTROL...

INSERTION AXIS

Undo Redo

Back Next

Save

Expert

Tools

Show distances

TruSmile

exocad

0136/04

Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Add Gripper for editing distal surface of the Preci-Vertex attachment and shoulder

The screenshot displays the exocad software interface for designing a primary telescope. The main 3D view shows a yellow and green dental model with several control points (blue spheres) and a red arrow pointing to a conical icon. A text box on the right says "Bring control point to distal".

- Adapt shoulder
- Occlusal area – mark icon (conical) between control points
- abutment shoulder / chamfer must be higher than the gingiva

Transform control point into conical icon by mouse click
Note: Redo for each control point

3.14 [0.0°]
0.93 [32.7°]

Bring control point to distal

Wizard Primary Telescope

TELESCOPE PARAMETERS

Angle

Diam (Q)

STICK (BOTTOM)

Crown Gingiva None

Visualize areas with friction

ADVANCED ANGLE CONTROL...

INSERTION AXIS

Undo Redo

Back Next

Save Expert Tools Show distances Tru Smile


exocad

Add custom view

Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Edit vestibular. distal surface and shoulder → press “next”

- Edit vestibular surface – bring control point to cervical
- Adapt incisal edge on full anatomy



Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Free-Form Telescope

➤ Add/Remove and Smooth/Flatten possible

Note: Telescope surfaces may not be free-formed

Wizard Free-Form Telescope ...

FREE ATT.

Add/Remove Smooth/Flatten

BRUSH

Strength (use Ctrl + mouse wheel) < >

Brush size (Shift + mouse wheel) < >

Type: [Icon] [Icon] [Icon]

Undo Redo

Back Next

exocad

6136/04

Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Add attachment

- Select attachment
- Library: DEDICAM-fits_for_Preci-Vertex
- Typ: DEDICAM-fits_for_Preci-Vertex length (4, 6 + 8mm)

Note:
Select attachment named DEDICAM

The screenshot displays the 'Free-Form Telescope' wizard in the exocad software. The interface is divided into three main sections. The top section shows the 'ATT.' tab selected, with 'Add' and 'Remove' buttons. The 'Library' dropdown is set to 'DEDICAM-fits_for_Preci-Vertex' and the 'Type' dropdown is set to 'DEDICAM-fits_for_Preci-Vertex_20150925'. The middle section shows a list of attachments, with 'DEDICAM-Kleberretention' selected. The bottom section shows the 'DEDICAM-fits_for_Preci-Vertex®-4mm-025-Lstl' attachment selected. The background shows a 3D model of a dental arch with a grey telescope attached. The exocad logo is in the bottom right corner.


Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Add attachment

- Place attachment on distal surface
- Lower max. the Preci-Vertex radius into distal surface
- Mark cut to gingiva

Note: Max. radius lowered on distal surface

Attention: Correction needed. otherwise Preci-Vertex matrice is unplaceable






Attaching a Preci-Vertex® with interlock and circumference

Primary telescope: Add attachment - correct position

- Preci-Vertex occlusal orientated
- Preci-Vertex must penetrate Gingiva
- Press „Apply“

Note:
Max. radius lowered on distal surface

Note:
Preci-Vertex must penetrate Gingiva



Attaching a Preci-Vertex® with interlock and circumference

Select interlock according approximal space

▼ EIN/AUSBLENDEN ?

- ☑ Kiefer-Scandaten
- ☑ Alle anderen Teile
- ☑ Anatomie
- ☑ Teleskope
- ☑ Angehängte Teile
- ☑ ZÄHNE

VERTECKT Anzeigen

➤ Select attachment

➤ Library: DEDICAM Interlock

➤ Typ: DEDICAM Interlock 1.5 + 3.0mm diameter possible

Wizard Free-Form Telescope ...

FREE ATT.

+ Add - Remove

Library: DEDICAM-Interlock_201

Type: DEDICAM-Interlock_20150810

DEDICAM-M14-forBredentScrew

DEDICAM-MK1-bolt_20150825

DEDICAM-Screw-channel 2018-04-06

generic

microtec_TK1

Rhein83

Si-tec

Sterngold

Wizard Free-Form Telescope ...

FREE ATT.

+ Add - Remove

Library: DEDICAM-Interlock_201

Type: DEDICAM-Interlock_1.5

DEDICAM-Interlock_1.5.stl

DEDICAM-Interlock_3.0.stl

Top View Surface

☑ Cut on gingiva 0 mm

☐ Allow any changes

Apply Unload

Undo Redo

Back Next

Note:
Attachment cut to Gingiva

exocad

Add custom view

0130/64


Attaching a Preci-Vertex® with interlock and circumference

Select interlock according approximal space

- Library: DEDICAM Interlock
- Typ: DEDICAM Interlock 1.5 + 3.0mm diameter possible
- Mark „Remove“

Attention:

Mark Remove; attachment color switches from YELLOW to RED




The screenshot displays the exocad software interface. On the left, a sidebar shows a tree view with options like 'Ein/Ausblenden', 'Kiefer-Scandaten', 'Alle anderen Teile', 'Anatomie', 'Teleskope', and 'Angehängte Teile'. The main workspace shows a 3D model of a dental arch with a yellow interlock being attached to a tooth. A red interlock is also visible. A wizard window titled 'Wizard Free-Form Telescope ...' is open, showing the 'Remove' button highlighted in red. The wizard window also shows the library 'DEDICAM-Interlock_201' and the type 'DEDICAM-Interlock 1.5'. The bottom right corner of the interface shows the 'exocad' logo and a navigation panel with 'Neue Ansicht', 'Ansicht 1', 'Ansicht 2', 'Ansicht 3', and 'Ansicht 4'.

Attaching a Preci-Vertex® with interlock and circumference

Vertical and horizontal interlock position

- Library: DEDICAM Interlock
- Type: DEDICAM Interlock 1.5 + 3.0mm diameter possible
- Mark „delete“



Positioning:

Vertical – identical to step/chamfer

Horizontal – balanced approximal space (Omegaform)


A screenshot of the CAD software interface. The 'Wizard Free-Form Telescope' dialog box is open, showing options for 'Add' and 'Remove' components. The 'Library' is set to 'DEDICAM-Interlock_201' and the 'Type' is 'DEDICAM-Interlock_1.5'. The 'INSERTION DIRECTION' section shows 'Top', 'View', and 'Surface' options. The 'Cut on gingiva' checkbox is checked, and the 'Allow any changes' checkbox is unchecked. The 'Apply' and 'Unload' buttons are visible. The 'Back' and 'Next' navigation buttons are also present.

A screenshot of the CAD software interface showing the 'Save' button and other tool icons. The 'Save' button is highlighted. Other icons include 'Expert', 'Tools', 'Show distances', 'TruSmile', and 'Add custom view'. The 'cad' logo is visible at the bottom.

Attaching a Preci-Vertex® with interlock and circumference

Vertical and horizontal interlock position

- After correct positioning → mark „Apply“
- **Attention:** Interlock will be removed only on one side
- Go to „Expert Mode“ to cut the other side too



Positioning:

Vertical – identical to step/chamfer

Horizontal – balanced approximal space (Omega shape)

Attaching a Preci-Vertex® with interlock and circumference

Vertical and horizontal interlock position after one side cut → go to Expert Mode

- Interlock must be repositioned
- Precise copy of removed interlock position recommended
- After correct placement → Mark “Apply”
- **Attention:** If necessary process needs to be redone → mark “Undo”

Note: Interlock will be removed only on one side
Redo placement on Expert Mode

The screenshot displays a CAD software interface for dental model assembly. The main workspace shows a 3D model of a dental arch with a grey interlock component and a red telescope component. A red arrow points to the interlock, and another red arrow points to the telescope. The interface includes a top-left menu with options like 'Kiefer-Scandaten', 'Alle anderen Teile', 'Anatomie', 'Teleskope', and 'Angehängte Teile'. A bottom-left panel titled 'Free-Form Telescope ...' contains controls for 'ATT.' (highlighted in red), '+Add', and 'X Remove' (highlighted in red). Below these are 'Move', 'Rotate', and 'Scale' buttons. The 'INSERTION DIRECTION' section has 'Top', 'View', and 'Surface' options. A checkbox 'Cut on gingiva' is checked, and 'Allow any changes' is unchecked. At the bottom of this panel are 'Apply', 'Undo' (highlighted in red), and 'Redo' buttons. The bottom toolbar features various icons, with 'Free-Form Telescope ...' (highlighted in red) and 'Delete Reconstructions' (highlighted in red) being prominent. A right-side toolbar includes 'Save', 'Wizard' (highlighted in red), 'Tools', 'Show distances', and 'TruSmile'. The bottom right corner shows 'Add custom view' and '6136/64'.

Attaching a Preci-Vertex® with interlock and circumference

Vertical and horizontal interlock position after cutting (omega shape) on Expert Mode. Now continue with connector.

- After correct positioning → “Apply” attachment
- „OK” to confirm and switch to connector design

The screenshot displays the DEDICAM software interface. On the left, a tree view shows the model structure with 'Teleskope' expanded. The main workspace shows a 3D dental model with a red telescope being positioned. The 'Free-Form Telescope' dialog box is open, showing options for 'Add', 'Remove', 'Move', 'Rotate', and 'Scale'. The 'Apply' button is highlighted in red. The 'OK' button is also highlighted in red. The 'Wizard' button in the top right toolbar is highlighted in red. The bottom toolbar contains various tools, with 'Free-Form Telescope' highlighted in red.


Attaching a Preci-Vertex® with interlock and circumference

Connector design → Interlock must stay untouched → otherwise redesign is required

- Select form / cross section → connector “free-forming”
- Shape connector by control points → Mark „Apply”

Connector: Form
Select value cross section
select basic shape
Apply changes

Connector „Free”
Connector design



Attaching a Preci-Vertex® with interlock and circumference

Connector design → Interlock must stay untouched → otherwise redesign is required

- Do not place the connector into the interlock → Mark “Apply”
- Press „OK“ to finalize design

The screenshot displays a CAD software interface for dental design. A 3D model of a tooth and its connector is shown. The connector is a cylindrical piece with a dashed line indicating its placement. Dimensions are shown: 4.47 (vertical), 2.09 (horizontal), and 7.54 (circumference). The 'Connectors' panel on the left shows 'SHAPE' set to 'FREE' and 'CONNECTOR CONTROL POINTS' with a green arrow pointing to a point on the connector. Below this, 'CONNECTOR THICKNESS' is shown with a color scale from 'Very thick' (blue) to 'Very thin' (red). The 'Wizard' button is highlighted with a red box. The 'Note' box contains the text: 'Note: Check interlock – Connector should not „fill“ interlock space'. The 'OK' button in the 'Connectors' panel is also highlighted with a red box. The software interface includes a top menu bar, a left sidebar with 'EB/AUSBLENDEN' and 'VERSTECKT' options, and a bottom toolbar with various tools like 'Insertion Direction', 'Bottoms', 'Copy Tooth', 'Load Custom Tooth Model', 'Shrink Telescopers Anatomy', 'Virtual Waxup Bottom', 'Free-form', 'Telescopers', 'Free-Form Telescope ...', 'Connector', 'Merge and save restorations', 'Model Alignment', and 'Delete Reconstructions'. The 'Connector' tool is highlighted with a red box. The right sidebar shows 'Speichern', 'Wizard', 'Werkzeuge', 'Abstände zeigen', and 'Tru Smile'.

Attaching a Preci-Vertex® with interlock and circumference

Connector design → Interlock must stay untouched → otherwise redesign is required

save design

→ stl file incl. *Constructioninfo* needs to be provided to production site

The screenshot shows the exocad software interface. On the left, a tree view shows a hierarchy of parts: Kiefer-Scandaten, Alle anderen Teile, Anatomie, Teleskope, Verbinder, and Fertige Teile. The 'Fertige Teile' folder is expanded. In the center, a 3D model of a dental restoration is shown. On the right, a vertical toolbar contains icons for Save, Expert, Tools, Show distances, and TruSmile. The 'Expert' icon is highlighted with a red box. At the bottom left, a 'Wizard Merge and save restorations' dialog is open, showing 'I'm done' as the selected next step. A note is overlaid on the 3D model.

Note:
It is not required to free-form finalized parts

exocad

Design of healing abutments

Design of healing abutments

General notes / prerequisites for custom healing abutments:

Custom healing abutments made of PEEK can be additionally ordered - without redesign - as individual impression posts for open or closed impression taking within the same order.

Availability (July 2022):

	CAMLOG®	CONOLOG®	CERALOG® ¹	iSy® on Implant shoulder	BioHorizons® ²	Further implant systems
Titanium healing abutments <small>¹ ZrO₂ for CERALOG</small>	✓	✓	✓	✓	✓	✓
PEEK healing abutments	✓	✓	✓	✓	✓	✗
PEEK impression posts <small>² BioHorizons: For open tray technique only</small>	✓	✓	✓	✓	✓	✗

Design of healing abutments

Order form:

1. Define tooth position
2. Select prosthetics
«Anatomic crown»
3. Scan mode:
e.g. Digital impression scan

Patient Multi-die

← Tooth 22 Material configuration (local): Default

Crowns and copings

- Anatomic crown
- Coping
- Pressed crown
- Eggshell crown (Provisional)
- Offset coping
- Mockup

Pontics

- Anatomic pontic
- Reduced pontic
- Pressed pontic
- Eggshell pontic (Provisional)

Inlays, onlays and veneers

- Inlay/Onlay
- Offset inlay
- Veneer

Digital copy milling

- Anatomic waxup
- Reduced waxup
- Pontic waxup

Removables and appliances

- Full denture
- Partial denture
- Bite splint
- Primary telescopic crown
- Secondary telescopic crown
- Attachment

Bars

- Bar pillar
- Bar segment
- Offset substructure

Legend:
■ Anatomic crown
■ Adjacent tooth
■ Antagonist

Workflow icons:
Scan
Design
Model Creator
Manufacture
dentalshare
Copy
Print
Settings...
About...


exocad

Tooth shade: A1
Scan mode: Digital impression scan

Design of healing abutments

Order form:





1. Select «Custom Abutment» from «Options and Parameters»
2. Select «PEEK Gingiva former» from Abutment material; Color «white»



Design of healing abutments

CAD library selection for CAMLOG, CONELOG, CERALOG, iSy and BioHorizons:

- 1 Example: DEDICAM CAMLOG library
- 2 Select your desired set of healing abutment / impression post
- 3 Implant diameter
- 4 Regular for CAMLOG: \varnothing 3.3 – 6.0
- 5 Platformswitched for CAMLOG: \varnothing 3.8PS – 6.0PS



Design of healing abutments




After import, alignment, cropping and data matching of the scan data, the healing abutment design starts.

1 Due to the selection of the anatomy (crown, bridge) a draft of the basic restoration shape is defined on the scan.

2 Options are „Copy“ or „Mirror“ of the corresponding tooth

3 as „Virtual pre-op“ or „Direct copy“





4 Tooth placement tools for positioning and scaling



Design of healing abutments

The additional anatomy (crown, bridge) supports the design of the healing abutments. Healing abutments are the basis, the foundation, for the prosthetic restoration.

If implant positions serve for bridges **always** create a bridge „virtual Wax-up“



Design of healing abutments



The anatomy supports the circular design of the healing abutments.

This is created comparable to the design of an abutment.

Note:




1 If needed add „Gripper“ to design an anatomical shape.
The "abutment shoulder" runs at the level of the gingiva.

2 The distance between the healing abutments and the neighboring tooth should be at least 1.0 - 1.5mm.



Design of healing abutments

- 1 After defining the outer contour and the subgingival design, the abutment shoulder is reduced to the value 0.1mm
- 2 Drag all „Gripper“ as far as possible to the abutment shoulder
- 3 Push the upper part of the gingiva former downwards with the green arrow
- 4 The distance between the healing abutment and the neighboring tooth should be at least 1.0 - 1.5mm



Design of healing abutments

Rounding and reducing the height of the healing abutments.



The healing abutment should protrude approx. 1.5mm above gingival level.

Reason: Soft tissue should not "grow" over the healing abutment during healing.

Depending on its size, the healing abutment displaces a corresponding volume of soft tissue; therefore, the height should not exceed the displaced gingiva.

1 Use the free-forming tools „Add/Remove“ and „Smooth/Flatten“ to shape the upper part

2 The minimum height of healing abutment must be considered




Design of healing abutments

The subgingival design should be discussed with the dentist.

Recommendations:

Lower third almost follows the implant diameter.

Upper two-thirds: concave / convex design to the cross section of the planned prosthetics.




Design of healing abutments

Minimum geometry

1


The inner blank limit regarding the screw seating cannot be ignored.



Design of healing abutments

Maximum geometry

- 1 To control the design in height and diameter the “Ti-Base/Preform” can be displayed
- 2 The maximum height of the healing abutment from implant shoulder is 7mm
- 3 The maximum diameter of the healing abutment is 9.9mm
- 4 The screw channel is always central in the maximum geometry
5. **Recommendation:** for large shapes, remove the transparency of the preform in order to better recognize the parts that stand out from the geometry



Design of healing abutments

Since individual healing abutments have oval, square / rectangular or triangular shapes, it is helpful to facilitate the correct insertion of the healing abutments for the dentist by means of marking / notching.

To apply a marker the following are suitable methods:

1 Free-Forming


2 Add/Remove

3 Remove brush size / medium strength

Note: additional smoothing is not necessary

4 **Recommendation:** A position of the marking/groove corresponding to the vestibular cam facilitates the correct insertion of the healing abutment and impression post for the dentist.


Always create only one marking/groove - no double grooves.



Design of healing abutments


Since individual healing abutments have oval, square/rectangular or triangular shapes, it is helpful for the practitioner to make additional dot-like markings when there are several "shape-matched" healing abutments.

1. Quadrant = one identification




2. Quadrant = two identifications

without identification



3. Quadrant = dot identification






4. Quadrant = line identifications

Design of healing abutments

After finishing the design, healing abutment stl-file is created by clicking "I`m done". **1**

Send design data via dentalshare or eService to DEDICAM production



Design of printed models

Design of printed models

General notes / requirements for printed DEDICAM models*:

- Activated Model Creator function on exocad license dongle
- Software version 2.2 Valletta or higher
- Software version 2.4 Plovdiv or higher allows to choose between several analogs in the Model Creator. In earlier versions, only one analog provider per implant family can be stored for system reasons. I.e., for CAMLOG® and CONELOG® either the original analogs from Camlog or the DIM analogs from NT.
- Software version 3.0 Galway allows an individual gap dimension to the analog. Therefore, a separate library is required. Please contact the DEDICAM Technical Service to use this function.
- Gingiva masks designed with version 3.0 Galway are produceable by DEDICAM print model partner *
- Checking and printing of the design data transmitted to Camlog is done by Innovation MediTech GmbH (Dreve)*

* This service might not be available in your country

Design of printed models


Overview of implant analogs available in the DEDICAM CAD libraries for printed models

	CAMLOG®	CONOLOG®	CERALOG®	iSy®	BioHorizons®	Further implant-systems
Original analogs from Camlog	✓	✓	✓	✓	✗	✗
DIM Analogs from NT-Trading	✓	✓	✗	✗	✓	✓

Design of printed models

From software version 2.4 Plovdiv it is possible to choose between several analogs in the Model Creator. In earlier versions, only one analog provider per implant family can be stored for system reasons. I.e. for CAMLOG® and CONELOG® either the original analogs from Camlog or the DIM analogs from NT.

- 1 Option original implant analog from Camlog
- 2 Option DIM analog from NT



For the selected model analog, multiple model analogs are available.

Select analog to use for tooth: 11

CONELOG - Camlog Analog
CAMLOG Biotechnologies GmbH

CONELOG - NT DIM Analog
Camlog Biotechnologies GmbH




Design of printed models

New features in version 3.0 Galway:

1 «Stump» settings

2 Value on «Shaft taper» generates conical stump shape




Design of printed models

New features in version 3.0 Galway:

«Settings»

- Default value «Horizontal distance» of 0,05mm must be reduced between **0.00** and **0.02mm** to ensure the proper fit of Camlog analogs for printed models.

Attention: Setting parameters for the model analogs are only available after library update by the DEDICAM Technical Service! Inhouse-manufacturers must determine the value adjusted to their printer themselves.



Design of an angulated screw channel


Design of an angulated screw channel

General notes:

The Titanium-base CAD-CAM free is suitable for a design with an angulated screw channel up to a maximum of 25° for CAMLOG and CONELOG. Exception: CONELOG® GH 2.0 mm → up to 15°.











The Titanium-base CAD-CAM free is contraindicated for bridges. The angulation of the screw channel results in the abutment screw being trapped in the Titanium-base after the restoration has been bonded.

To tighten the abutment screw with the angulated screw channel, the "Ballpoint" screwdriver Art. no. J5319.050x must always be used.



Design of an angulated screw channel

The overview table shows the dependencies between the prosthetic components and a design with or without an angulated screw channel

	CAMLOG	Screw	CONELOG	Screw	Screwdriver
Angulated screw channel design (Ti-base CAD/CAM free)	 K2247.xxxx (short) K2265.xxxx (long)	 J4005.1601 or J4005.2001	 C2247.xxxx (short) C2265.xxxx (long)	 C4015.1601 or C4015.2001	 J5319.050x
Straight screw channel design (Ti-base CAD/CAM)	 K2244.xxxx	 J4005.1601 or J4005.2001	 C2242.xxxx	 C4015.1601 or C4015.2001	 J5317.0502

Design of an angulated screw channel

General notes:

The BioHorizons Hybrid titanium base is suitable for an angulated screw channel design of up to a maximum of 15°.








The BioHorizons Hybrid titanium base with an angulated screw channel is contraindicated for bridges.

With an angulated screw channel, the yellow anodized abutment screw, and the "Precision Angled" screwdriver. Art. no. BZ5334.2014 must be used.




Design of an angulated screw channel

The overview table shows the dependencies between the BioHorizons components and a design with an angulated screw channel


	Titanium base hexed	Titanium base non-hexed	Screw	Screwdriver
Straight screw channel design	 <p>All diameters</p>	 <p>All diameters</p>	 <p>BZ4015.0010 (PXMUAS)</p>	 <p>BZ5305.5020 (135-351)</p>
Angulated screw channel design	 <p>All diameters</p>	<p>X not possible</p>	 <p>BZ4022.0010 (PXPAS)</p>	 <p>BZ5334.2014 (PADM14)</p>

Design of an angulated screw channel

1. Select the implant system CAMLOG or CONELOG from the «Select implant analog» box



2. Selection of the diameter and the height of the chimney, and if available, the gingival height




Example: CAMLOG 3.8 Ti-base free long K2265.3848


Type and Ø Chimney height Art.-Nr.

Design of an angulated screw channel

1. Select the desired scan body and perform Best-Fit-Matching




2. The CAD/CAM base free can be rotated to the preferred orientation using the arrow



Design of an angulated screw channel

In the "Generate abutment bases" step under "Advanced", the angulated screw channel function appears. Either by clicking with the mouse **1.** or by dragging **2.** the green dot in the desired direction, the individual inclination of the screw channel can be set.


Note: The inclination of the screw channel causes the screw cannot be removed once the restoration is bonded to the titanium base.



Design of an angulated screw channel

In addition to the angulation of the screw channel, its positioning to the anti-rotation surfaces must be considered.

If the screw channel is positioned too far to the left or right and thus covers the titanium base's anti-rotation surface, the restoration's anti-rotation function will be lost. There is no control on the production side.



Disclaimer

The information provided does not qualify the viewer to adopt or implement the product in a clinical setting. For proper use of the product(s), please refer to the relevant instructions for use (IFU) and work instructions.

BioHorizons® Implants and Multi-unit abutments are manufactured by BioHorizons.

CAMLOG®, CERALOG®, CONELOG® and iSy® products, in addition to all custom-made DEDICAM® products, are manufactured by: ALTATEC GmbH | Maybacherstr. 5 | 71299 Wimsheim | Germany

CAMLOG®, CERALOG®, COMFOUR®, CONELOG®, DEDICAM®, iSy® and Tube-in-Tube® are registered trademarks of CAMLOG Biotechnologies GmbH.

BioHorizons® is a registered trademark of BioHorizons.

The referenced trademarks may, however, not be registered in all markets. Any mentioned third-party trademarks are the property of the respective trademark owner. CAMLOG Biotechnologies GmbH or BioHorizons respectively may market products only upon clearance and approval by the respective national authorities. Therefore, not all products may be available in your country.

Media ID: M-1733-PRT-EN-INT-BHCL-00-112023

© CAMLOG Biotechnologies GmbH, 2022. All rights reserved.