

# PROSTHETIC REHABILITATION WITH BIMAXILLARY IMPLANT-SUPPORTED FIXED DENTAL PROSTHESES USING THE CAMLOG® VARIO SR SYSTEM

a perfect fit™



Prof. a.c. Marco Finotti  
Università di Padova, Italy



Graduated at the University of Padua with a degree in Medicine and Surgery, specializing in Orthognathodontics, member of the "Medical Association" and the "Odontoiatric Association" of Milan, training in "Professional development program of bone biology" at Harvard and in "Continuing, education in surgical periodontal therapy" at Gothenburg, Professor of Periodontal Surgery and Implant-prosthetics at the University Popolare of Caserta for the 1996-97 academic year and at the University Popolare of Milan for the 1997-98 academic year, Adjunct Professor at the University of Padua since 1998-99, self-employed in Padua and Milan, carrying out exclusively surgical activity, speaker at National and International courses and conventions, author of many publications in the field of implantology, has participated to an important experimental research on the effects of orthodontic loads on implants, receiving an award for the best presentation at the SIDO 2000 conference, active member of S.I.O. (Italian Osteointegration Society), member of S.I.D.O. (Italian Orthodontic Society), Fouding Member of the Piezosurgery Academy.

## IMPLANTS USED

Tooth	18	17	16	15	14	13	12	11	21	22	23	24	25	26	27	28
Impl. type				SL		SL	SL			SL	SL		SL			
Impl. length				13		13	13			13	13		11			
Impl. Ø				4.3		4.3	4.3			4.3	4.3		4.3			
Impl. surface				PP		PP	PP			PP	PP		PP			

Tooth	48	47	46	45	44	43	42	41	31	32	33	34	35	36	37	38
Impl. type					SL	SL	SL			SL	SL	SL				
Impl. length					13	13	13			13	13	13				
Impl. Ø					3.8	3.8	3.8			3.8	3.8	3.8				
Impl. surface					PP	PP	PP			PP	PP	PP				

Impl. type: ROOT-LINE (RL)/SCREW-LINE (SL) Impl. surface: Promote® (P)/Promote® Plus (PP)

## PROSTHETICS

- standard
  - platform switching
  - removable
  - fixed
  - crown
  - bridge
  - cement-retained
  - screw-retained
  - partially edentulous
  - fully edentulous
  - other
- 
- Universal abutment
  - Esthomic® abutment
  - Telescope abutment
  - Gold-plastic abutment
  - Ceramic abutment
  - CAD/CAM abutment on titanium base
  - Temporary abutment
  - Logfit® abutment
  - Locator® abutment
  - Bar abutment
  - Ball abutment
  - Vario SR abutment
  - other

## INFORMATION ON PATIENT AND TREATMENT

The patient, a 63 year old smoker (15 cigarettes a day for 40 years) presented himself with a severe periodontal compromised situation of the residual dentition. In the upper jaw he has a fixed restoration from 15 to 24 with a porcelain to metal fused bridge presenting a mobility of grad 3 based on the now irreversibly damaged periodontal situation. In the lower jaw despite the lingual splint we find mobility between 33 to 43, the dental elements 34 and 35 present a second grade mobility and tooth 48 is mesially inclined. The patient reports a disturbed TMJ function with audible noise during mouth opening and closure, besides pain which is diffusing into the facial muscles and paravertebral region.

After an appropriate and adequate clinical and radiographic evaluation (panoramic x-ray, CT scan) the choice was to extract the remaining teeth in both jaws with a provisional rehabilitation by means of two removable dentures to :

- Functionally re-educate the patient
- Recover the horizontal and vertical dimensions
- Resolve the mal coordination condyle meniscal eliminating the muscular and articulation disturbances through possible modifications of the removable dentures

During gnathologic therapy the patient was treated surgically with the insertion of 6 implants per jaw to support, after successful osseointegration, an adapted fixed restoration to recover the exact mastication with correct occlusal tables. The gnathological situation resolved, during the phase of osseointegration of the implants the occlusal parameters are reproduced with prosthetic guides which duplicate perfectly the provisional restorations transferring thus the individual therapeutic data. The impression for the realization of the provisional dentures is created by underlining the prosthetic guide with silicon. The precision and passivity of the final prosthetic suprastructures are achieved using splinted plaster transfers. The provisional resin dentures with a metal reinforcement are screw-retained on the abutments with a torque of 25 Ncm.

### Initial situation




Fig. 1: Radiographic evaluation of the initial situation.




Fig. 2: Panoramic radiograph after tooth extraction.



Fig. 3: Removable temporary dentures.

### Impression taking



Fig. 4: Soft tissue healing after tooth extraction in the mandible.



Fig. 5: Impression with the temporary dentures to create the radiographic and surgical templates.




Fig. 6: Lab-putty impression to duplicate.



Fig. 7: Upper surgical stent.



Fig. 8: Upper prosthetic guides.



Fig. 9: Upper and lower prosthetic guides.

### Implant placement



Fig. 10: Upper and lower surgical stents.




Fig. 11: Flap preparation with a Quantum Molecular Resonance Bistoury.



Fig. 12: Alveoplasty with piezosurgery.



Fig. 13: Window preparation for sinus grafting.



Fig. 14: Careful implant insertion in the area of the sinus floor elevation.



Fig. 15: 6x CAMLOG® SL 4.3 mm implants in the maxilla.

### Bone Augmentation procedure



Fig. 16: Healing caps and sutures with an apically repositioned flap.



Fig. 17: 6x CAMLOG® SL 3.8 mm implants in the mandible.



Fig. 18: Autogenous bone and Bio-Oss®.



Fig. 19: Submerged healing with matrix and continuous sutures in the mandible.



Fig. 20: Radiographic control.



Fig. 21: Healing after the second stage surgery in the mandible.



Fig. 22: Try-in of prosthetic guides.



Fig. 23: Prosthetic guides with implant position references.



Fig. 24: Inserted Vario SR abutments.

### Impression taking on abutment level



Fig. 25: Open tray impression in the maxilla.



Fig. 26: Silicon rebased prosthetic guide with the transfer copings attached with resin to it in the maxilla.



Fig. 27: Silicon rebased prosthetic guide with the transfer copings attached with resin to it in the mandible.



Fig. 28: Master cast with Vario SR analogs.



Fig. 29: Vario SR impression copings splinted with plaster for a passiv-fit/precision try-in.



Fig. 30: Vario SR impression copings splinted with plaster for a passiv-fit/precision try-in.



Fig. 31: Verification of passivity and precision of both impressions.



Fig. 32: Upper and lower master cast with tooth set-up in wax.



Fig. 33: Try-in of the tooth set-up, verification of esthetics.



Fig. 34: Test of the teeth over the metal reinforcement on the model.



Fig. 35: Basal view of the final denture.



Fig. 36: Softtissue healing in the maxilla.



Fig. 37: Screw retained temporary in the maxilla.



Fig. 38: Softtissue healing in the mandible.



Fig. 39: Screw retained provisional in the mandible.

## CONCLUSIONS

The use and positioning of Vario SR Abutments helps to stabilize the soft tissues at gingiva level due to the fact that the abutments will not be removed anymore once the provisional restoration is inserted.

The passivation, first of the provisional and then the final restoration on the implants is supported by the conical characteristics of the Vario SR Abutments. This prevents any tension, which, if present, might be responsible for a prosthetic failure due to overload, even with correctly positioned implants. In addition, perfect passivity insures the preservation of the aesthetics, short and long terms stability. The convenient removability of the screw retained prosthesis, if necessary, is possible without interfering with the gingival seal.

The stability is granted by the screw retaining system: screws with adequate dimensions combined with a large number of threads and appropriate depth avoiding thus any friction or tension of the screw joint, granting therefor the stability of the fixation.

The possibility to remove the rehabilitation will provide inspection over time, either at prosthetic level or at the implant-prosthetic emergence. This applies as well to the professional hygiene recalls.

This type of rehabilitation, on 6 implants, in accordance to the international literature shows great restorative efficacy in edentulous patients, in single or both jaws. The correct positioning of the implants will help to preserve the marginal bone and simple routine maintenance of the restoration for the patient to his greatest satisfaction.

## Final restoration in situ

---



Fig. 40: Buccal view final restoration.



Fig. 41: Detailed view of anterior teeth and artificial gums.



Fig. 42: Radiographic evaluation after insertion of the final restoration.

## REFERENCES

- SLOT W, GERRY M, RAGHOEBAR GM, VISSINK A, HENNY JA, MEIJER HJ.** Maxillary overdentures supported by anteriorly or posteriorly placed implants opposed by a natural dentition in the mandible: a 1-year prospective case series study. *J Clin Periodontol* **2013** Mar; 40(3): 303–10
- GULIÈ F, RAGHOEBAR GM, TER MEULEN JW, VISSINK A, MEIJER HJ.** Mandibular overdentures supported by 6-mm dental implants: a 1-year prospective cohort study. *Clin Implant Dent Relat Res* **2012** May; 14 Suppl 1: e59–66. doi: 10.1111/j.1708-8208.2011.00358.x. [epub ahead of print]
- ANTOUN H, BELMON P, CHERFANE P, SITBON JM.** Immediate loading of four or six implants in completely edentulous patients. *Int J Periodontics Restorative Dent* **2012** Feb; 32(1): e1–9
- RAGHOEBAR GM, MEIJER HJ, STELLINGSMA K, VISSINK A.** Addressing the atrophied mandible: a proposal for a treatment approach involving endosseous implants. *Int J Oral Maxillofac Implants* **2011**; 26: 607–17
- ANDRIOTELLI M, ATT W, STRUB JR.** Prosthodontic complications with implant overdentures: systematic literature review. *Int J Prosthodont* **2010**; 23: 195–203
- OHKUBO C, BAEK KW.** Does the presence of antagonistic remaining teeth affect implant overdenture success? A systematic review. *J Oral Rehabil* **2010**; 37: 306–12
- SLOT W, RAGHOEBAR GM, VISSINK A, HUDDLESTON SLATER JJ, MEIJER HJ.** A systematic review of implant-supported maxillary overdentures after a mean observation period of at least 1 year. *J Clin Periodontol* **2010**; 37: 98–110
- AGLIETTA M, SICILIANO VI, ZWAHLEN M, BRAGGER U, PJETURSSON BE, LANG NP, ET AL.** A systematic review of the survival and complication rates of implant supported fixed dental prostheses with cantilever extensions after an observation period of at least 5 years. *Clinical Oral Implants Research* **2009**; 20: 441–51
- CARLSSON GE.** Dental occlusion: modern concepts and their application in implants prosthodontics. *Odontology* **2009**; 97: 8–17
- YOUNIS L, TAHER A, ABU-HASSAN MI, ET AL.** Evaluation of bone healing following immediate and delayed dental implant placement. *J Contemp Dent Pract* **2009**; 10: 35
- SADOWSKY SJ.** Treatment consideration for maxillary implant overdentures: systematic review. *J Prosthet Dent* **2007**; 97: 340–48
- STEINEBRUNNER L, WOLFART S, BÖSSMANN K, KERN M.** In vitro evaluation of bacterial leakage along the implant-abutment interface of different implant systems. *Int J Oral Maxillofac Implants*. **2005** Nov–Dec; 20(6): 875–81
- STELLINGSMA K, SLAGTER AP, STEGENGA B, RAGHOEBAR GM, MEIJER HJ.** Masticatory function in patients with an extremely resorbed mandible restored with mandibular implant-retained overdentures: comparison of three types of treatment protocols. *J Oral Rehabil* **2005**; 32: 403–10
- YAMAMOTO H, MAEDA Y, SOGO M, MURAKAMI S.** Minimum number of implants for occlusal support in edentulous mandible. *J Rest Dent* **2004**; (abstract 0821)
- JACOBS R.** Preoperative radiologic planning of implants in compromised patients. *Periodontology* 2000; 33: 12, **2003**
- VISSER A, RAGHOEBAR GM, MEIJER HJ, VAN T HOF MA, STEGENGA B, VISSINK A.** A randomized prospective clinical trial on the effectiveness of three treatment modalities for patients with lower denture problems. A 10 year follow-up study on patient satisfaction. *Int J Oral Maxillofac Surg* **2003**; 32: 498–503
- DULA K, MINI R, VAN DER STELT P, ET AL.** The radiographic assessment of implant patients: Decision-making criteria. *Int J Oral Maxillofac Implants* **2001**; 16: 80
- SCHWAETZ-ARAD D, GULAYEV N, CHAUSHU G.** Immediate versus non-immediate implantation for full-arch fixed reconstruction following extraction of all residual teeth: a retrospective comparative study. *J Periodontol* **2000**; 71: 923
- STANFORD CM.** Biomechanical and functional behavior of implants. *Adv Dent Res* **1999**; 13: 88
- HOBKIRK JA, HAVTHOULAS TK.** The influence of mandibular deformation, implant numbers, and loading position on detected forces in abutments supporting fixed implant superstructures. *J Prosthet Dent* **1998**; 80: 169
- CHIAPASCO M, GATTI C, ROSSI E, ET AL.** Implant-retained mandibular overdentures with immediate loading: Results on 226 consecutive cases. *Clin Oral Implants Res* **1997**; 8: 48
- TARNOW DP, EMTIAZ S, CLASSI A.** Immediate loading of threaded implants at stage 1 surgery in edentulous arches. Ten consecutive case reports with 1 to 5-year data. *Int J Oral Maxillofac Implants* **1997**; 12: 319
- ASSIF D, MARSHAK B, SCHMIDT A.** Accuracy of implant impression techniques. *Int J Oral Maxillofac Implants* **1996**; 11: 216–22
- CARR AB, BRUNSKI JB, HURLEY E.** Effects of fabrication, finishing and polishing procedures on preload in prostheses using conventional 'gold' and plastic cylinders. *Int J Oral Maxillofac Implants* **1996**; 11: 589–98

**BRÅNEMARK PI, SVENSSON B, VAN STEENBERGHE D.**

Ten-year survival rates of fixed prostheses on four or six implants ad modum Brånemark in full edentulism. Clin Oral Implants Res **1995** Dec; 6(4): 227–31

**RANGERT B, KROGH PH, LANGER B, ET AL.**

Bending overload and implant fracture: A retrospective clinical analysis. Int J Oral Maxillofac Implants **1995**; 10: 326

**SPECTOR MR, DONOVAN TE, NICHOLLS JL.**

An evaluation of impression techniques for osseointegrated implants. J Prosthet Dent **1990**; 63: 444–47

**CONTACT**

Prof. a c. Finotti Marco

Via Sorio 12

I-35141 Padova

Tel: 0039/049723413

Fax: 0039/049723302

e-mail: [fi-med@studiofinotti.it](mailto:fi-med@studiofinotti.it)

[www.centromedicovesalio.it](http://www.centromedicovesalio.it)

**HEADQUARTERS**

CAMLOG Biotechnologies AG | Margarethenstrasse 38 | CH-4053 Basel | Switzerland

Tel +41 61 565 41 00 | Fax +41 61 565 41 01 | [info@camlog.com](mailto:info@camlog.com) | [www.camlog.com](http://www.camlog.com)

**camlog**