NovoMatrix™ Reconstructive Tissue Matrix

Instructions
For Use

Table of Contents

Sommaire / Inhaltsverzeichnis / Indice / Índice / İçindekiler / Índice

Instructions for Use English	3–6
Instrucciones de uso Castellano; Español	7–10
Gebrauchsanweisung Deutsch	11–14
Notice d'utilisation Français.....	15–18
Istruzioni per l'uso Italiano.....	19–22
Kullanım Talimatları Türkçe.....	23–26
Instruções de Utilização Português.....	27–30

Device Description

NovoMatrix™ Reconstructive Tissue Matrix (NovoMatrix RTM) is a collagen matrix that is derived from porcine dermis and then processed and preserved in a patented aqueous phosphate buffered solution containing matrix stabilizers. NovoMatrix RTM consists of a sterilized sheet of processed porcine dermis provided in prescribed sizes and packaged in a double pouch.

Use of NovoMatrix RTM provides for a strong and biocompatible implant and will incorporate into the patient's tissue with associated cell and microvascular ingrowth.

Indications For Use

NovoMatrix RTM is indicated for:

- Localized gingival augmentation to increase keratinized tissue (KT) around teeth
- Alveolar ridge reconstruction for prosthetic treatment
- Guided tissue regeneration procedures in recession defects for root coverage

Contraindications

- NovoMatrix RTM should not be placed where active infection or inflammation exists. Before placement, the physician should be confident that any active or recent infection has been properly treated.
- Patients with known allergy to collagen should not be treated with NovoMatrix RTM.
- NovoMatrix RTM is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.
- Polysorbate 20 is a component of the aqueous phosphate buffered solution and therefore NovoMatrix RTM should not be used in patients with a known sensitivity to this material.

Warnings

- Do not reuse or resterilize. This product is only for single patient one-time use.
- Do not use if the package is opened or damaged.
- Do not use if seal is broken or compromised.
- Do not use if the temperature monitoring device located on the product carton does not display "OK".
- After use, handle and dispose of all unused product and packaging in accordance with accepted medical practices and applicable national and regional environmental laws on disposal of packaging and biological materials.
- NovoMatrix RTM cannot be reused once it has been removed from the packaging and/or is in contact with a patient without increased risk of patient-to-patient contamination and subsequent infection.

Precautions

- Discard NovoMatrix RTM if handling has caused possible damage or contamination.
- Discard NovoMatrix RTM if it is past the use-by-date of the product [as indicated as year (4 digits), month (2 digits), day (2 digits)].

- Ensure that NovoMatrix RTM is put into a sterile basin and covered with room temperature sterile saline or room temperature sterile lactated Ringer's solution for a minimum of 2 minutes prior to implantation.
- Place NovoMatrix RTM in maximum possible contact with healthy, well-vascularized tissue to promote cell ingrowth and tissue remodeling.
- NovoMatrix RTM should be hydrated and moist when the package is opened. If NovoMatrix RTM is dry, do not use.
- Performance of NovoMatrix RTM may be impacted by patients' conditions such as uncontrolled metabolic diseases and immune disorders or radiotherapy in the oral cavity.

Adverse Events

Potential adverse events which may result from placement of an implant or graft include, but are not limited to, the following: wound or systemic infection, inflammation, soft tissue dehiscence, hypersensitive allergic or other immune response.

Special instructions for use in Periodontology

A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene.

Storage

- NovoMatrix RTM is a sterile medical device that should be stored in a clean, dry location at room temperature.
- Refer to the temperature monitor included to ensure that product has been stored within its temperature limits. Only use the product if the included temperature monitor displays "OK" on the screen. If screen displays anything other than "OK," do not use the product.*
- NovoMatrix RTM should be stored in its original packaging.
- The use-by-date of the product is indicated as year (4 digits), month (2 digits), day (2 digits).

Sterilization:

This product is sterilized by electron beam irradiation.

Instructions For Preparing NovoMatrix RTM For Surgical Use

These instructions are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Clinicians should be familiar with surgical procedures and techniques involving periodontal procedures before using NovoMatrix RTM.

Required Materials

- Sterile forceps.
- Soaking fluid: room temperature sterile saline or room temperature sterile lactated Ringer's solution.
- One sterile basin per piece of NovoMatrix RTM.

Preparation Instructions

1. Open the carton and remove the pouch.
2. Peel open the outer pouch and remove the inner pouch using aseptic technique. The inner pouch is sterile and may be placed directly into the sterile field.
3. Open the inner pouch carefully and aseptically remove the matrix. Always use sterile gloves or forceps when handling NovoMatrix RTM.
4. Place NovoMatrix RTM in a sterile basin covered by room temperature sterile saline or room temperature sterile lactated Ringer's solution for a minimum of 2 minutes.
5. Keep NovoMatrix RTM submerged in the room temperature sterile solution until ready for implantation. NovoMatrix RTM can be retained in sterile solution for a maximum of 4 hours.


Implantation Instructions

1. Prepare the surgical site using standard techniques.
2. NovoMatrix RTM may be folded, trimmed or cut as required to fit the surgical site using aseptic technique, ensuring allowance for overlap.
3. Transfer NovoMatrix RTM to the surgical site using sterile gloves or forceps.
4. Suture NovoMatrix RTM into place. The device should be sutured tension-free to the surrounding tissue and should be covered by the flap.
5. Complete the standard surgical procedure.
6. Discard any unused portions of NovoMatrix RTM as per institutional procedures.

Post-Operative Care

Proper post-operative care should be practiced to minimize contamination to the surgical wound.

Definitions

	USE BY DATE		INNER POUCH		OUTER POUCH
	CONTENTS		CORRECT SETTING		*INCORRECT SETTING
	DO NOT USE IF PACKAGE IS DAMAGED		OPEN HERE		CONSULT INSTRUCTIONS FOR USE
	DO NOT RE-STERILIZE		STERILIZED USING IRRADIATION		TEMPERATURE LIMITATIONS
	FOR SINGLE USE ONLY		MANUFACTURED BY		AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

Customer Support

*Call customer support on day of receipt for immediate return authorization.

For product complaints and potential adverse events, please contact your local Sales Representative or Customer Support.

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NOVOMATRIX™ is a trademark of BioHorizons
Allergan® and its design are trademarks of Allergan, Inc.
© 2019 Allergan. All rights reserved.
LifeCell Corporation
Allergan Pharmaceuticals International Limited