
















SCREW OVERVIEW

For the CAMLOG[®] Implant system



ABUTMENT AND PROSTHETIC SCREWS INTRAORAL USE

Implant-Abutment connection






	Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	Ø 6.0 mm		
	M 1.6			M 2.0		Tightening torque	
CAMLOG® Abutment screw							
 <p>Temporary Abutments PEEK, incl. PS</p> <p>Scanbody</p> <p>ScanPost for Sirona Scanbody</p>						tightening by hand*	
 <p>Temporary Abutments titanium, crown and bridge</p>							
 <p>Esthomic® Abutments, incl. PS</p>						20 Ncm	
 <p>Universal Abutment, incl. PS</p> <p>Telescope Abutment</p> <p>Gold-plastic Abutment</p> <p>Logfit® Abutment</p>	 <p>10.5 mm</p> <p>J4005.1601</p>	 <p>10.5 mm</p> <p>J4005.2001</p>					
 <p>Ceramic Abutment</p>							
 <p>Titanium bases CAD/CAM, crown and bridge</p>							
 <p>Vario SR Abutments, 20° and 30° angled</p>							
CAMLOG® Vario SR Abutment screws							
 <p>Vario SR Abutment, straight</p>	 <p>11.9 mm</p> <p>J4007.1600</p>	 <p>11.9 mm</p> <p>J4007.2000</p>			20 Ncm		
CAMLOG® Abutment screws with reduced head, light blue anodized							
 <p>COMFOUR™</p> <p>Bar Abutments, 17° and 30° angled</p>	 <p>9.5 mm</p> <p>J4004.1601</p>	 <p>9.5 mm</p> <p>J4004.2001</p>			20 Ncm		

PS: Platform Switching

*Optionally für temporary abutments titanium: tightening torque after successful healing phase 20 Ncm













ABUTMENT AND PROSTHETIC SCREWS INTRAORAL USE

Abutment-Prosthetic connection

		Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	Ø 6.0 mm		
		M 1.6			M 2.0		Tightening torque	
Prosthetic screws for bar abutments, light blue anodized								
 <p>COMFOUR™ Bar Abutments, 17° and 30° angled</p>	3.6 mm	 J4012.1601		3.8 mm	 J4012.2001		15 Ncm	
	Vario SR Prosthetic screw, yellow anodized							
 <p>Vario SR Abutments, straight, 20° and 30° angled</p>	4 mm	 J4005.2004						15 Ncm









AUXILIARY SCREWS INTRA- AND EXTRAORAL USE

Abutment-Prosthetic connection

		Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	Ø 6.0 mm		
		M 1.6			M 2.0		Tightening torque	
Prosthetic screws for bar abutments, light blue anodized								
 <p>Scanning cap for bar abutments</p>	3.6 mm	 J4012.1601		3.8 mm	 J4012.2001		tightening by hand	
	Screws for bar abutments, for impression taking open tray and for soldering, light blue anodized							
 <p>COMFOUR™ Bar abutments, straight, 17° and 30° angled</p>	12 mm	 J4012.1610		12.2 mm	 J4012.2010		tightening by hand	
	17 mm	 J4012.1615		17.2 mm	 J4012.2015			
	22 mm	 J4012.1620		22.2 mm	 J4012.2020			
	Plastic screws for bar abutment, as fixation and bonding aid, beige							
	29 mm	 J4009.1627		29.2 mm	 J4009.2027		tightening by hand	

LAB SCREWS EXTRAORAL

Lab analog-Abutment connection











	Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	Ø 6.0 mm	
	M 1.6			M 2.0		Tightening torque
CAMLOG® Lab screws*, brown anodized						
 <p>Temporary Abutments PEEK, incl. PS</p> <p>Scanbody</p> <p>ScanPost for Sirona Scanbody</p>						tightening by hand
 <p>Temporary Abutments titanium, crown and bridge</p>						
 <p>Esthomic® Abutments, incl. PS</p>						
 <p>Universal Abutment, incl. PS</p> <p>Telescope Abutment</p> <p>Gold-plastic Abutment</p>	<p>10.5 mm</p>  <p>J4006.1601</p>			<p>10.5 mm</p>  <p>J4006.2001</p>		
 <p>Ceramic Abutment</p>						
 <p>Titanium bases CAD/CAM, crown and bridge</p>						
 <p>Vario SR Abutments, 20° and 30° angled</p>						
CAMLOG® Bonding aids, white (enclosed with the titanium bases CAD/CAM)						
 <p>Titanium bases CAD/CAM, crown and bridge</p>	<p>27.5 mm</p>  <p>J4009.1600</p>			<p>27.5 mm</p>  <p>J4009.2000</p>		tightening by hand
CAMLOG® Vario SR Lab screws*, brown anodized						
 <p>Vario SR Abutment, straight</p>	<p>11.9 mm</p>  <p>J4008.1600</p>			<p>11.9 mm</p>  <p>J4008.2000</p>		tightening by hand
CAMLOG® Lab screws with reduced head*, light blue partially anodized						
 <p>COMFOUR™</p> <p>Bar Abutments, 17° and 30° angled</p>	<p>9.5 mm</p>  <p>J4004.1600</p>			<p>9.5 mm</p>  <p>J4004.2000</p>		tightening by hand

PS: Platform Switching

*Note: Lab screws may not be used on patients.

LAB SCREWS EXTRAORAL

Abutment-Prosthetic connection

		Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	Ø 6.0 mm	
		M 1.6			M 2.0		Tightening torque
Lab prosthetic screws for bar abutments*, brown anodized							
	Scanning cap for bar abutments						
	COMFOUR™ Bar Abutments, straight, 17° and 30° angled	3.6 mm  J4013.1601			3.8 mm  J4013.2001		tightening by hand
	Bar lab analog for bar abutments						
Vario SR Prosthetic screw, yellow anodized							
	Vario SR Abutments, straight, 20° and 30° angled				4 mm  J4005.2004		tightening by hand
	Vario SR Analog						
Prosthetic screw for bar abutments*, for fabrication of the wax up on the bar sleeve for titanium bonding base, Passive-Fit, on the bar lab analog							
	Titanium bonding base for bar abutments and bar sleeve for titanium bonding base, burn-out, Passive-Fit	5.5 mm  J4005.1602			5.5 mm  J4005.2002		tightening by hand

*Note: Lab screws may not be used on patients.



HEADQUARTERS

CAMLOG Biotechnologies AG | Margarethenstr. 38 | 4053 Basel | Switzerland
Phone +41 61 565 41 00 | Fax +41 61 565 41 01 | info@camlog.com | www.camlog.com

Manufacturer CAMLOG®/CONOLOG® Products: ALTATEC GmbH, Maybachstr. 5, 71299 Wimsheim, Germany

