


CONELOG®
SYSTEM

PRODUCT CATALOG 2014

Valid from May 2014

a perfect fit™



camlog

ADDED VALUE THROUGH COMPETENCE IN IMPLANTOLOGY

A partnership with CAMLOG encompasses the exceptional price-performance ratio of the CONELOG® Products and a custom-built service program. CAMLOG's specific expertise clearly encompasses more than the product, its safety and ease of use. With all our implant expertise, we are fully committed to our customers. Rediscover the implant expertise of CAMLOG in your daily routine time and time again, and turn it into your own added value.

We are happy to help.

FOR PERFECTIONISTS

FOR MINIMALISTS

FOR FANS



THE CONELOG® IMPLANT SYSTEM SYSTEM INFORMATION	2	CONELOG® PROSTHETICS	37
NEW PRODUCTS	2	PROSTHETICS OVERVIEW	38
SURGERY	11	TORQUE WRENCH SETTINGS	42
PLANNING	12	IMPRESSION TAKING, BITE REGISTRATION, FABRICATION OF THE PLASTER MODEL	44
X-RAY PLANNING FOILS AND TRANSFER PICTURES	12	IMPRESSION TAKING, OPEN/CLOSED TRAY	44
CT-PLANNING	13	BITE REGISTRATION	44
IMPLANT OVERVIEW	14	FABRICATION OF THE PLASTER MODEL	45
CONELOG® SCREW-LINE IMPLANTS	15	CROWN AND BRIDGE RESTORATIONS	45
SURGERY SET CAMLOG®/CONELOG® SCREW-LINE	16	CONELOG® TEMPORARY ABUTMENT	45
SURGICAL INSTRUMENTS SCREW-LINE	17	CONELOG® ESTHOMIC® ABUTMENTS	45
OSTEOTOMY SET SCREW-LINE	20	CONELOG® UNIVERSAL ABUTMENT	47
GUIDE SYSTEM	24	CONELOG® GOLD-PLASTIC ABUTMENT	47
CONELOG® HEALING CAPS	30	CONELOG® LOGFIT® ABUTMENTS	47
ALTAPIN SET	31	ACCESSORIES FOR CONELOG® ABUTMENTS	48
INSTRUMENTS	32	CROWN, BRIDGE AND HYBRID RESTORATIONS	49
SURGICAL INSTRUMENTS	32	CONELOG® TITANIUM BASES CAD/CAM	49
SURGICAL-PROSTHETIC INSTRUMENTS	34	ACCESSORIES FOR CAD/CAM PROSTHETICS	50
PROSTHETIC INSTRUMENTS	35	CONELOG® VARIO SR ABUTMENTS	51
		ACCESSORIES FOR VARIO SR	53
		HYBRID RESTORATIONS	54
		CONELOG® BAR ABUTMENTS	54
		CONELOG® BALL ABUTMENTS	56
		CONELOG® LOCATOR® ABUTMENTS	58
		CONELOG® UNIVERSAL ABUTMENT	61
		CONELOG® TELESCOPE ABUTMENT	61
		ACCESSORIES FOR CONELOG® ABUTMENTS	61
		INSTRUMENTS FOR DENTAL TECHNICIANS	62
		CONELOG® SELECTION ABUTMENTS	63
		AUXILIARY ARTICLES	65
		INFORMATION	69
		INDEX	73
		ALPHABETICAL	74
		ARTICLE NUMBER	81
		ORDER FORM	87

THE CONELOG® IMPLANT SYSTEM

SYSTEM INFORMATION

GENERAL SYSTEM INFORMATION ABOUT THE CAMLOG® AND CONELOG® IMPLANT SYSTEM

The CONELOG® Implant System is based on the many years of clinical and laboratory experience with the CAMLOG® Implant System. The CONELOG® Implant System is a user-friendly, consistently prosthesis-oriented implant system.

All CAMLOG® and CONELOG® Products are continuously fabricated based on the latest state of technology. The CAMLOG® and CONELOG® Implant System was progressively enhanced through in-house research and development in collaboration with hospitals, universities and dental technicians. Today, it represents the state of the art in implant dentistry.

The CAMLOG® Implant System is well documented scientifically. Numerous studies based on various parameters, e.g. implant surface, time of implantation and/or implant loading, primary stability, connection design or type of superstructure, support this. The long-term outcomes for the CAMLOG® Implant System are persuasive.

NEW PRODUCTS

GUIDE SYSTEM

The CAMLOG® Guide System for template-guided implant bed preparation and insertion has been extended by the CONELOG® SCREW-LINE Implants and is now called Guide System.

Guide System CONELOG® SCREW-LINE Implants are available with a 3.3 mm diameter in the lengths 9, 11, 13 and 16 mm as well as with 3.8 and 4.3 mm diameters in the lengths 7, 9, 11, 13 and 16 mm.

Newly available are guided pilot drills and a driver with ISO shaft for Guide System implants. The form drills were fitted with a four-edged cutting geometry with chip grooves to provide greater smoothness of operation. The guiding sleeves were shortened by 1 mm for minimally invasive implantation at greater gingival thickness, provided with opposing bevels for use in narrow tooth gaps and fitted with a marking to give better groove alignment. The template drills are fitted newly with a spiral-shaped cutting edge and a longer shaft. All drills, gingiva punches, insertion posts and guiding sleeves are color-coded throughout.

The new Guide System comprises 85 items.

LOCATOR®

A new replacement male is available for the Locator® Anchoring system for extended range without retention. The replacement male is colored gray and serves to adapt the cumulated retention force of a Locator® retained prosthesis on implants with diameters 3.8/4.3/5.0 mm.

CONELOG® SCANPOSTS FOR SIRONA SCANBODY

In conjunction with the Sirona Scanbodies the new CONELOG® ScanPosts serve for the digital recording of the implant position or laboratory analog position in relation to the remaining teeth and the soft tissue.

CONELOG® TITANIUM BASE CAD/CAM

In future, the Sirona inCoris ZI meso blocks can also be used with the CONELOG® Implant System. The newly modified CONELOG® Titanium base CAD/CAM, including the dark purple, anodized abutment screw and a bonding aid, will be available in the gingival heights 0.8 mm and 2.0 mm from September 2014 with a new article number.

THE CONELOG® IMPLANT SYSTEM

SYSTEM INFORMATION


CONELOG® SCREW-LINE IMPLANTS

CONELOG® SCREW-LINE Implants are endosseous implants, available in different lengths and diameters. They are placed surgically in the maxillary and/or mandibular bone and serve as anchors for functional and esthetic oral rehabilitations in partial or fully edentulous patients. The prosthetic treatment is then completed with single crowns, bridges or full dentures, which are fastened to the corresponding CONELOG® Components with the CONELOG® SCREW-LINE Implants. Different implant sizes can be used in the same arch.

The CONELOG® SCREW-LINE Implant represents a conical self-tapping screw implant in its geometry and is available with the Promote® surface.

Implantation

CONELOG® SCREW-LINE Implants are intended as immediate or delayed immediate implants. The selected healing technique can be either submerged or transgingival. The taper of the outer implant body of 3° – 9° (depending on length and diameter) allows easy insertion due to self-centering.



CONELOG® SCREW-LINE
Implant, Promote® plus

THE CONELOG® IMPLANT SYSTEM


SYSTEM INFORMATION

Promote® surface

The CONELOG® SCREW-LINE Implant is available with the Promote® surface. The abrasive-blasted, acid-etched Promote® surface extends over the entire body of the implant to the acid-etched 45° implant shoulder.

The Promote® surface has proven its worth as a surface for anchoring dental implants in the bone and has shown good scientific results in osteohistology and in pull-out trials.


The results suggest that the Promote® surface leads to rapid osseointegration of the CONELOG® Implants.




Promote® surface

CONELOG® SCREW-LINE IMPLANT INNER CONFIGURATION

CONELOG® SCREW-LINE Implants are equipped with a cone (7.5°) and three grooves in the inner configuration for positioning CONELOG® Abutments. The CONELOG® Abutments are apical with a cone and three cams, and lock into the conical connection and the three grooves of the implant. The CONELOG® Abutment does not cover the implant shoulder. A CONELOG® Abutment screw is used to fix CONELOG® Abutments in the CONELOG® SCREW-LINE Implant with a defined torque.



Conical CONELOG® Implant abutment connection



THE CONELOG® IMPLANT SYSTEM

SYSTEM INFORMATION


For optimal positioning of the abutments in the implant, they should be aligned in the bone so that one of the three grooves points vestibularly. The drivers include markings on the outside that correspond to the three grooves of the CONELOG® Implant inner configuration.




CONELOG® IMPLANT/BAR ABUTMENT/BALL ABUTMENT CONNECTION

Different CONELOG® Abutments, CONELOG® Bar abutments, CONELOG® Ball abutments, and CONELOG® Locator® Abutments with several geometries for anchoring of implant-retained full prostheses are available for the CONELOG® Implant System. These abutments differ from each other in two different connection types of their apical regions.

CONELOG® Bar, Ball and Locator® Abutments are equipped with a thread in the apical region which engages with the inner thread of the CONELOG® Implant. These abutments are screwed into the CONELOG® Implant using the corresponding drivers at a defined torque.




Due to the design of the screw connection, CONELOG® Bar, Ball, and Locator® Abutments have no cams.



Example: CONELOG® Ball abutment (\varnothing 4.3 mm) in a CONELOG® SCREW-LINE Implant

PRODUCTION PRECISION

The inner and outer geometry of the CONELOG® Implants and abutments are rotary machined for the most part. The tolerances can therefore be kept very low. The result is excellent part precision without impacting the material structure. The CONELOG® Implant abutment connection ensures a very precise, stable and rotation-resistant connection to the CONELOG® Prosthetic components.



CONELOG® SCREW-LINE Implant




CONELOG® Abutment

CONELOG® COMPONENTS

For the CONELOG® SCREW-LINE Implants, some CONELOG® Components are available such as CONELOG® Cover screws, CONELOG® Healing caps, CONELOG® Impression posts and CONELOG® Prosthetic components.



CONELOG® Cover screw
CONELOG® Healing cap bottleneck, wide body, cylindrical



CONELOG® Impression post
open and closed tray

THE CONELOG® IMPLANT SYSTEM

SYSTEM INFORMATION

CONELOG® PROSTHETIC COMPONENTS

CONELOG® SCREW-LINE Implants can be provided with a pallet of wide-ranging, anatomically appropriate prosthetic components. Esthetic solutions are easily achieved. All CONELOG® Abutments are color-coded by implant diameter.

CONELOG® Temporary abutment


The CONELOG® Temporary abutment (titanium alloy) can be used for immediate implantation or after opening for temporary restoration.



CONELOG® Esthomic® Abutments

The CONELOG® Esthomic® Abutments include an anatomical shoulder and are available in various gingival heights.


- CONELOG® Esthomic® Abutment, straight
- CONELOG® Esthomic® Abutment, 15° and 20° angled, in type A and B
- CONELOG® Esthomic® Abutment, inset



Angled CONELOG® Esthomic® Abutments, type A and B

Angled abutments type A and B are available in the CONELOG® Esthomic® Line of abutments. Type A has a cam opposing the angle direction. In contrast to type A, type B has a 60° offset cam. Both types make six different rotation positions possible to achieve an optimal prosthetic axis alignment.

CONELOG® Abutment cam alignments



CONELOG® Gold-plastic abutment

The CONELOG® Gold-plastic abutment can be used to fabricate single-crowns, individual implant abutments (mesostructures) for cementable bridge restorations and primary abutments for bridging implant axis divergences in the telescopic crown technique using the cast-on technique.



THE CONELOG® IMPLANT SYSTEM

SYSTEM INFORMATION


CONELOG® Titanium base CAD/CAM

The CONELOG® Titanium base CAD/CAM enables restorations on CONELOG® Implants with individual and highly precise two-piece abutments made from zirconium oxide, titanium or cast alloy. The base is available in the gingival heights 0.8 and 2.0 mm and is used as a bonding base for individual, implant-supported reconstructions such as mesostructures as well as crown, bridge and double-crown restorations.



CONELOG® Logfit® Abutments

The Logfit® Prosthetic System enables the fabrication of cement-retained crown and bridge restorations. The Logfit® Prosthetic System consists of color-coded CONELOG® Logfit® Abutments with two gingival heights (1.0 and 2.5 mm), Logfit® Impression caps, Logfit® Analogs and Logfit® Burn-out plastic copings with and without rotation protection for the manufacture of cast crown and bridge reconstructions.



CONELOG® Vario SR abutments

CONELOG® Vario SR abutments are available in straight, in 20° and 30° angled versions. These abutments can compensate for large implant axis divergences. Prefabricated Vario SR prosthetic components make fabrication of occlusally screw-retained crown, bridge and bar restorations possible. Burn-out plastic copings can be used to fabricate casted crown and bridge restorations. The titanium caps have a retention surface on the outside and are designed for temporary or final bridge restorations made of plastic. Titanium bases for bars are available for laser-welded bar restorations.



IMPORTANT NOTE

Because of the conical inner configuration of the CONELOG® SCREW-LINE Implants, they are only compatible with CONELOG® Components.

CONELOG® Universal and Telescope abutment

The CONELOG® Universal and Telescope abutments can be used for individually fabricated cementable crown and bridge restorations and for telescopic crown restorations. The abutments are made of titanium alloy and can be custom trimmed.



CONELOG® DISCONNECTOR FOR CONELOG® ABUTMENTS

The CONELOG® Disconnector for CONELOG® Abutments is used to remove or decouple CONELOG® Abutments from the CONELOG® Implants/Lab analogs. If the abutment does not come loose, the torque wrench (locked setting) can be placed on the disconnector and the abutment can be loosened by turning clockwise.



CONELOG® Bar-, Ball- and Locator® Abutments

For the CONELOG® Implant System, CONELOG® Bar, Ball, and Locator® Abutments are available in different gingival heights for anchoring of implant-retained full prostheses. For CONELOG® Bar abutments, there are several bases available for burn-out, cast-on, soldering, and laser-welding techniques. CONELOG® Bar-, Ball- and Locator® Abutments are perfectly suited for movable anchorings of full prostheses.



CONELOG® Bar, Ball and Locator® Abutment

THE CONELOG® IMPLANT SYSTEM

SYSTEM INFORMATION

GENERAL SAFETY INSTRUCTIONS AND WARNINGS

The descriptions in this product catalog are not sufficient to allow immediate use of the CONELOG® Implant System. Instruction by a surgeon experienced in using the CONELOG® Implant System is strongly recommended.

An improper procedure during use of the CONELOG® Products can lead to failures, implant damage, bone loss or unsatisfactory esthetic results. The products may be used only by dentists, surgeons and dental technicians trained in the implant system. Use of these products requires specialized knowledge and experience in implant dentistry. Detailed information on the choice of suitable implants, prosthetic components, treatment planning and application is contained in the user information. Your local representative will be happy to provide you with information.

Consulting and technical advice on the use of our products is provided orally, in writing, by electronic media and/or by demonstration. The provided information represents the current state of science and technology known to us at the point in time that the product is placed on the market. This does not exempt the user from the responsibility to personally test the product for suitability for the intended purpose, indication and procedure. Handling and use of the product take place outside of our control and is the direct responsibility of the user. All liability for damages resulting from such use is waived. ALTATEC GmbH/CAMLOG Biotechnologies AG does not warrant nor provide replacement services when non-system components are used.

The CONELOG® Implant System is part of a comprehensive treatment concept and must be used only with the pertinent original parts and tools according to the suggestions and instructions for use provided by the manufacturer. All components of the CONELOG® Implant System are matched precisely to one another. The use of third party components can affect the function and safety of the system.

Instruments and system components are designed for specific implant and prosthetic lines and implant diameters. Their use with any other implant or prosthetic product lines or a different diameter can lead to the mechanical failure of system components, tissue injury, or esthetically unsatisfactory results. Because of this, some implant and prosthetic product lines have some dedicated components and tools. Pay attention to the color markings when choosing instruments for the required implant diameter.

Products intended for single use must not be reused because safe preparation and/or functional safety cannot be ensured.

Within the bounds of our sales and delivery terms, we guarantee the flawless quality of our products.

Because of the small sizes involved, it may happen that a product is swallowed and/or aspirated. Aspiration can lead to dyspnea and in the worst case to asphyxiation. For this reason, products must be secured in general against aspiration and swallowing during intraoral use.

Where indications are listed for a particular product, it should be noted that any indications that are not listed are in fact contraindicated.

Sterile packaged products must be kept dry, out of direct sunlight and at room temperature. The packaging must be checked before opening for damage and the expiration date and should only be opened immediately before the products are used.

Not all products are available in all countries.

PACKAGING UNITS

If not otherwise specified, the packaging unit is 1 each.

EXPLANATION OF SYMBOLS/CHARACTERS ON LABELS/PRODUCT LABELS AND PACKAGE INSERTS

STERILE	R	Sterilization using irradiation
		Non-sterile
		Caution, consult accompanying documents
		Use by
		Do not reuse
REF		Article number
LOT		Lot number
		Manufacturer
		Temperature limitation
Rx only		Caution: US Federal law restricts this device to sale by or on the order of a dentist or physician.

COLOR-CODING OF THE SURGICAL AND PROSTHETICAL CONELOG® PRODUCTS

Color	Diameter
gray	3.3 mm
yellow	3.8 mm
red	4.3 mm
blue	5.0 mm

THE CONELOG® IMPLANT SYSTEM

SYSTEM INFORMATION

Exterior packaging:

Sealed, color-coded, folding carton with product label


Secondary packaging:

Sealed

Implant internal packaging (primary packaging):

Sealed, color-coded

Example product label of exterior implant packaging



SURGERY	11
PLANNING	12
X-RAY PLANNING FOILS AND TRANSFER PICTURES	12
CT-PLANNING	13
IMPLANT OVERVIEW	14
CONELOG® SCREW-LINE IMPLANTS	15
SURGERY SET CAMLOG®/CONELOG® SCREW-LINE	16
SURGICAL INSTRUMENTS SCREW-LINE	17
OSTEOTOMY SET SCREW-LINE	20
GUIDE SYSTEM	24
CONELOG® HEALING CAPS	30
ALTAPIN SET	31
INSTRUMENTS	32
SURGICAL INSTRUMENTS	32
SURGICAL-PROSTHETIC INSTRUMENTS	34
PROSTHETIC INSTRUMENTS	35

SURGERY


X-RAY PLANNING FOILS / X-RAY TRANSFER PICTURES

Art. No.

Article


∅
in mm

C5300.9010




**X-Ray Planning foil 1.25:1,
for CONELOG® SCREW-LINE Implants,
Magnification 25%**

C5300.9011



**X-Ray Planning foil 1.4:1,
for CONELOG® SCREW-LINE Implants,
Magnification 40%**

C5300.9080



C5300.9081

3.3

C5300.9082

3.8

C5300.9083

4.3

5.0

SURGERY
CT-PLANNING

Art. No.	Article	Material
A2002.2000	 CT-tube for drill Ø 2.0 mm, corrugated tubing Length 4.0/10.0 mm internal diameter 2.1 mm external diameter 2.5 mm (10 units)	Titanium alloy
A2222.2200	 CT-tube for drill Ø 2.2 mm, corrugated tubing Length 4.0/10.0 mm internal diameter 2.3 mm external diameter 2.7 mm (10 units)	Titanium alloy
A2050.2600	 Drill for CT-tube (for A2002.2000), Ø 2.6 mm	Stainless steel
A2050.2800	 Drill for CT-tube (for A2222.2200), Ø 2.8 mm	Stainless steel


SURGERY

IMPLANT OVERVIEW


CONELOG® SCREW-LINE IMPLANTS, PROMOTE® PLUS

	\varnothing 3.3 mm	\varnothing 3.8 mm	\varnothing 4.3 mm	\varnothing 5.0 mm
Length 7 mm		 C1062.3807	 C1062.4307	 C1062.5007
Length 9 mm	 C1062.3309	 C1062.3809	 C1062.4309	 C1062.5009
Length 11 mm	 C1062.3311	 C1062.3811	 C1062.4311	 C1062.5011
Length 13 mm	 C1062.3313	 C1062.3813	 C1062.4313	 C1062.5013
Length 16 mm	 C1062.3316	 C1062.3816	 C1062.4316	 C1062.5016
\varnothing apical (mean value)		2.7 mm	3.5 mm	3.9 mm
				4.6 mm

Note: The implant length is the distance from the apical curve to the occlusal plane surface of the implant.



SURGERY
CONELOG® SCREW-LINE IMPLANTS

Art. No.	Article	∅ in mm	Length in mm	Material
C1062.3309			9	Titanium Grade 4
C1062.3311	CONELOG® SCREW-LINE Implant, Promote® plus incl. CONELOG® Insertion post and CONELOG® Implant cover screw	3.3	11	
C1062.3313			13	
C1062.3316			16	
C1062.3807			7	
C1062.3809			9	
C1062.3811		3.8	11	
C1062.3813			13	
C1062.3816			16	
C1062.4307			7	
C1062.4309			9	
C1062.4311		4.3	11	
C1062.4313			13	
C1062.4316			16	
C1062.5007			7	
C1062.5009			9	
C1062.5011		5.0	11	
C1062.5013			13	
C1062.5016			16	

All CONELOG® Implants are packed sterile.

Surface: Promote®

Important notes

- Preparation of the implant bed for CONELOG® SCREW-LINE Implants and for CAMLOG® SCREW-LINE Implants is identical.
- CONELOG® SCREW-LINE Implants can only be supplied with CONELOG® Prosthetic components.

SURGERY

SURGERY SET CAMLOG®/CONELOG® SCREW-LINE

Art. No.

Article

J5300.0057

Surgery Set CAMLOG®/CONELOG® SCREW-LINE

contains all necessary color-code ordered surgical instruments, incl. torque wrench and holding key for insertion post



Note: The preparation of the implant bed for CAMLOG® SCREW-LINE Implants and for CONELOG® SCREW-LINE Implants is identical.

SURGERY

SURGICAL INSTRUMENTS SCREW-LINE

Art. No.	Article	\varnothing in mm	Length in mm	Material
J5050.2300*	 Round bur resterilizable	2.3	—	Stainless steel
J5050.3500	 Round bur resterilizable	3.5	—	Stainless steel
J5051.2003	 Pilot drill without coil, resterilizable	2.0	—	Stainless steel
J5051.2000*	 Pilot drill SCREW-LINE reduced coil, resterilizable	2.0	—	Stainless steel
J5051.2800*	 Pre-drill SCREW-LINE resterilizable	1.7-2.8	—	Stainless steel
J5015.0007*	 Depth stop SCREW-LINE	—	7	Stainless steel
J5015.0009*	 for pilot drill (J5051.2000) and pre-drill (J5051.2800) with reduced coil, resterilizable	—	9	Stainless steel
J5015.0011*	 —	—	11	Stainless steel
J5015.0013*	 —	—	13	Stainless steel

* The instruments listed above are also included in the surgery set CAMLOG®/CONELOG® SCREW-LINE.

SURGERY

SURGICAL INSTRUMENTS SCREW-LINE

Art. No.	Article	\varnothing in mm	Length in mm	Material
J5052.3309			9	
J5052.3311			11	
J5052.3313		3.3	13	
J5052.3316			16	
J5052.3807			7	
J5052.3809			9	
J5052.3811		3.8	11	
J5052.3813			13	
J5052.3816			16	
J5052.4307			7	
J5052.4309			9	
J5052.4311		4.3	11	
J5052.4313			13	
J5052.4316			16	
J5052.5007			7	
J5052.5009			9	
J5052.5011		5.0	11	
J5052.5013			13	
J5052.5016			16	
J5015.3300		3.3	–	Titanium alloy
J5015.3800		3.8	–	
J5015.4300		4.3	–	
J5015.5000		5.0	–	

The instruments listed above are also included in the surgery set CAMLOG®/CONELOG® SCREW-LINE.

SURGERY
SURGICAL INSTRUMENTS SCREW-LINE

Art. No.	Article	\varnothing in mm	Length in mm	Material
J5053.3316	 Form drill SCREW-LINE, Cortical bone, resterilizable	3.3	–	Stainless steel
J5053.3816	 Form drill SCREW-LINE, Cortical bone, resterilizable	3.8	–	
J5053.4316	 Form drill SCREW-LINE, Cortical bone, resterilizable	4.3	–	
J5053.5016	 Form drill SCREW-LINE, Cortical bone, resterilizable	5.0	–	
J5054.3309	 Tap SCREW-LINE, with hexagon, resterilizable	3.3	–	Stainless steel
J5054.3809	 Tap SCREW-LINE, with hexagon, resterilizable	3.8	–	
J5054.4309	 Tap SCREW-LINE, with hexagon, resterilizable	4.3	–	
J5054.5009	 Tap SCREW-LINE, with hexagon, resterilizable	5.0	–	

The instruments listed above are also included in the surgery set CAMLOC®/CONELOG® SCREW-LINE.

SURGERY

OSTEOTOMY SET SCREW-LINE

Art. No.	Article	\varnothing in mm	Material
J5418.0020			
			
Osteotomy Set CAMLOG®/CONELOG® SCREW-LINE, straight convex			
Content:			
J5417.2800	Pre-Osteotome SCREW-LINE, straight convex	1.7-2.8	
J5418.3300	Osteotome SCREW-LINE, straight convex	3.3	
J5418.3800	Osteotome SCREW-LINE, straight convex	3.8	
J5418.4300	Osteotome SCREW-LINE, straight convex	4.3	
J5418.5000	Osteotome SCREW-LINE, straight convex	5.0	
J5418.6000	Osteotome SCREW-LINE, straight convex	6.0*	
* for CAMLOG® SCREW-LINE Implant \varnothing 6.0 mm only.			
J5417.2800	Pre-Osteotome SCREW-LINE straight convex	1.7-2.8	Stainless steel
J5418.3300	Osteotome SCREW-LINE straight convex	3.3	Stainless steel
J5418.3800		3.8	
J5418.4300		4.3	
J5418.5000		5.0	

SURGERY

OSTEOTOMY SET SCREW-LINE

Art. No.	Article	\varnothing in mm	Material
J5418.0030			
Osteotomy Set CAMLOG®/CONELOG® SCREW-LINE, angled convex			
Content:			
J5417.2800	Pre-Osteotome SCREW-LINE, straight convex	1.7-2.8	
J5418.3310	Osteotome SCREW-LINE, angled convex	3.3	
J5418.3810	Osteotome SCREW-LINE, angled convex	3.8	
J5418.4310	Osteotome SCREW-LINE, angled convex	4.3	
J5418.5010	Osteotome SCREW-LINE, angled convex	5.0	
J5418.6010	Osteotome SCREW-LINE, angled convex	6.0*	
* for CAMLOG® SCREW-LINE Implant \varnothing 6.0 mm only.			
J5417.2800		1.7-2.8	Stainless steel
J5418.3310		3.3	Stainless steel
J5418.3810		3.8	
J5418.4310		4.3	
J5418.5010		5.0	


SURGERY

OSTEOTOMY SET SCREW-LINE

Art. No.	Article	\varnothing in mm	Material
J5420.0020			
			
Osteotomy Set CAMLOG®/CONELOG® SCREW-LINE, straight concave			
Content:			
J5419.2800	Pre-Osteotome SCREW-LINE, straight concave 1.7-2.8	1.7-2.8	
J5420.3300	Osteotome SCREW-LINE, straight concave 3.3	3.3	
J5420.3800	Osteotome SCREW-LINE, straight concave 3.8	3.8	
J5420.4300	Osteotome SCREW-LINE, straight concave 4.3	4.3	
J5420.5000	Osteotome SCREW-LINE, straight concave 5.0	5.0	
J5420.6000	Osteotome SCREW-LINE, straight concave 6.0*	6.0*	
* for CAMLOG® SCREW-LINE Implant \varnothing 6.0 mm only.			
J5419.2800	Pre-Osteotome SCREW-LINE straight concave	1.7-2.8	Stainless steel
J5420.3300	Osteotome SCREW-LINE straight concave	3.3	Stainless steel
J5420.3800		3.8	
J5420.4300		4.3	
J5420.5000		5.0	

SURGERY

OSTEOTOMY SET SCREW-LINE

Art. No.	Article	\varnothing in mm	Material
J5420.0030			
			
Osteotomy Set CAMLOG®/CONELOG® SCREW-LINE, angled concave			
Content:			
J5419.2800 Pre-Osteotome SCREW-LINE, straight concave 1.7-2.8			
J5420.3310 Osteotome SCREW-LINE, angled concave 3.3			
J5420.3810 Osteotome SCREW-LINE, angled concave 3.8			
J5420.4310 Osteotome SCREW-LINE, angled concave 4.3			
J5420.5010 Osteotome SCREW-LINE, angled concave 5.0			
J5420.6010 Osteotome SCREW-LINE, angled concave 6.0*			
* for CAMLOG® SCREW-LINE Implant \varnothing 6.0 mm only.			

J5419.2800		Pre-Osteotome SCREW-LINE straight concave	1.7-2.8	Stainless steel
J5420.3310		Osteotome SCREW-LINE angled concave	3.3	Stainless steel
J5420.3810			3.8	
J5420.4310			4.3	
J5420.5010			5.0	

SURGERY

GUIDE SYSTEM

TEMPLATE-GUIDED IMPLANT BED PREPARATION AND IMPLANTATION

Art. No.	Article	\varnothing in mm	Length in mm	Material
Guide System Implants				
C1063.3309		3.3	9	Titanium Grade 4
C1063.3311		3.3	11	
C1063.3313		3.3	13	
C1063.3316		3.3	16	
C1063.3807		3.8	7	
C1063.3809		3.8	9	
C1063.3811		3.8	11	
C1063.3813		3.8	13	
C1063.3816		3.8	16	
C1063.4307		4.3	7	
C1063.4309		4.3	9	
C1063.4311		4.3	11	
C1063.4313		4.3	13	
C1063.4316		4.3	16	

All CONELOG® Implants are packaged sterile.

Surface: Promote®

Guide System Pilot Drills

J5043.3309		Guide System Pilot drill set, for \varnothing 3.3 mm, internal irrigation, sterile, Content: Guide System pilot drills \varnothing 2.0 mm, length 5 and 9 mm	3.3	9	Stainless steel
J5043.3311		Guide System Pilot drill set, for \varnothing 3.3 mm, internal irrigation, sterile, Content: Guide System pilot drills \varnothing 2.0 mm, length 5, 9 and 11 mm	3.3	11	
J5043.3313		Guide System Pilot drill set, for \varnothing 3.3 mm, internal irrigation, sterile, Content: Guide System pilot drills \varnothing 2.0 mm, length 5, 9, 11 and 13 mm	3.3	13	
J5044.3316		Guide System Pilot drill, for \varnothing 3.3 mm, internal irrigation, sterile, \varnothing 2.0 mm, length 16 mm	3.3	16	
J5043.4307		Guide System Pilot drill set, for \varnothing 3.8/4.3 mm, internal irrigation, sterile, Content: Guide System pilot drills \varnothing 2.0 mm, length 5 and 7 mm	3.8 4.3	7	
J5043.4309		Guide System Pilot drill set, for \varnothing 3.8/4.3 mm, internal irrigation, sterile, Content: Guide System pilot drills \varnothing 2.0 mm, length 5 and 9 mm	3.8 4.3	9	

Note: All Guide System drills and gingiva punches are intended for single use only!

SURGERY GUIDE SYSTEM

TEMPLATE-GUIDED IMPLANT BED PREPARATION AND IMPLANTATION

Art. No.	Article	\varnothing in mm	Length in mm	Material
Guide System Pilot Drills				
J5043.4311	 Guide System Pilot drill set, for $\varnothing 3.8/4.3$ mm , internal irrigation, sterile, Content: Guide System pilot drills $\varnothing 2.0$ mm, length 5, 9 and 11 mm	3.8 4.3	11	Stainless steel
J5043.4313	 Guide System Pilot drill set, for $\varnothing 3.8/4.3$ mm , internal irrigation, sterile, Content: Guide System pilot drills $\varnothing 2.0$ mm, length 5, 9, 11 and 13 mm	3.8 4.3	13	
J5044.4316	 Guide System Pilot drill for $\varnothing 3.8/4.3$ mm , internal irrigation, sterile, $\varnothing 2.0$ mm, length 16 mm	3.8 4.3	16	

Guide System Surgery Sets and Form Drills

J5045.3309	 Guide System Surgery set, SCREW-LINE, $\varnothing 3.3$ mm , internal irrigation, sterile, Content: Guide System pre-drill $\varnothing 3.3$ mm, length 5 mm and Guide System form drill $\varnothing 3.3$ mm, length 9 mm	3.3	9	Stainless steel
J5045.3311	 Guide System Surgery set, SCREW-LINE, $\varnothing 3.3$ mm , internal irrigation, sterile, Content: Guide System pre-drill $\varnothing 3.3$ mm, length 5 mm and Guide System form drills $\varnothing 3.3$ mm, length 9 and 11 mm	3.3	11	
J5045.3313	 Guide System Surgery set, SCREW-LINE, $\varnothing 3.3$ mm , internal irrigation, sterile, Content: Guide System pre-drill $\varnothing 3.3$ mm, length 5 mm and Guide System form drills $\varnothing 3.3$ mm, length 9, 11 and 13 mm	3.3	13	
J5046.3316	 Guide System Form drill, SCREW-LINE, $\varnothing 3.3$ mm , internal irrigation, sterile, $\varnothing 3.3$ mm, length 16 mm	3.3	16	

Note: All Guide System drills and gingiva punches are intended for single use only!

SURGERY

GUIDE SYSTEM

TEMPLATE-GUIDED IMPLANT BED PREPARATION AND IMPLANTATION

Art. No.	Article	\varnothing in mm	Length in mm	Material
Guide System Surgery Sets and Form Drills				
J5045.3807	 Guide System Surgery set, SCREW-LINE, Ø 3.8 mm, internal irrigation, sterile, Content: Guide System pre-drill Ø 3.8 mm, length 5 mm and Guide System form drills Ø 3.8 mm, length 7 mm	3.8	7	Stainless steel
J5045.3809	 Guide System Surgery set, SCREW-LINE, Ø 3.8 mm, internal irrigation, sterile, Content: Guide System pre-drill Ø 3.8 mm, length 5 mm and Guide System form drills Ø 3.8 mm, length 9 mm	3.8	9	
J5045.3811	 Guide System Surgery set, SCREW-LINE, Ø 3.8 mm, internal irrigation, sterile, Content: Guide System pre-drill Ø 3.8 mm, length 5 mm and Guide System form drills Ø 3.8 mm, length 9 and 11 mm	3.8	11	
J5045.3813	 Guide System Surgery set, SCREW-LINE, Ø 3.8 mm, internal irrigation, sterile, Content: Guide System pre-drill Ø 3.8 mm, length 5 mm and Guide System form drills Ø 3.8 mm, length 9, 11 and 13 mm	3.8	13	
J5046.3816	 Guide System Form drill, SCREW-LINE, Ø 3.8 mm, internal irrigation, sterile, Ø 3.8 mm, length 16 mm	3.8	16	
J5045.4307	 Guide System Surgery set, SCREW-LINE, Ø 4.3 mm, internal irrigation, sterile, Content: Guide System pre-drill Ø 4.3 mm, length 5 mm and Guide System form drills Ø 4.3 mm, length 7 mm	4.3	7	
J5045.4309	 Guide System Surgery set, SCREW-LINE, Ø 4.3 mm, internal irrigation, sterile, Content: Guide System pre-drill Ø 4.3 mm, length 5 mm and Guide System form drills Ø 4.3 mm, length 9 mm	4.3	9	
J5045.4311	 Guide System Surgery set, SCREW-LINE, Ø 4.3 mm, internal irrigation, sterile, Content: Guide System pre-drill Ø 4.3 mm, length 5 mm and Guide System form drills Ø 4.3 mm, length 9 and 11 mm	4.3	11	
J5045.4313	 Guide System Surgery set, SCREW-LINE, Ø 4.3 mm, internal irrigation, sterile, Content: Guide System pre-drill Ø 4.3 mm, length 5 mm and Guide System form drills Ø 4.3 mm, length 9, 11 and 13 mm	4.3	13	
J5046.4316	 Guide System Form drill, SCREW-LINE, Ø 4.3 mm, internal irrigation, sterile, Ø 4.3 mm, length 16 mm	4.3	16	

Note: All Guide System drills and gingiva punches are intended for single use only!

SURGERY GUIDE SYSTEM

TEMPLATE-GUIDED IMPLANT BED PREPARATION AND IMPLANTATION

Art. No.	Article	\varnothing in mm	Length in mm	Material
Guide System Form Drills Cortical Bone				
J5048.3309	 Guide System Form drill, SCREW-LINE, Cortical bone, internal irrigation, sterile	3.3	9	Stainless steel
J5048.3311		3.3	11	
J5048.3313		3.3	13	
J5048.3316		3.3	16	
J5048.3807		3.8	7	
J5048.3809		3.8	9	
J5048.3811		3.8	11	
J5048.3813		3.8	13	
J5048.3816		3.8	16	
J5048.4307		4.3	7	
J5048.4309		4.3	9	
J5048.4311		4.3	11	
J5048.4313		4.3	13	
J5048.4316		4.3	16	

Note: All Guide System drills and gingiva punches are intended for single use only!

SURGERY

GUIDE SYSTEM

TEMPLATE-GUIDED IMPLANT BED PREPARATION AND IMPLANTATION

Art. No.	Article	\varnothing in mm	Length in mm	Material
Guide System Gingiva Punches, sterile				
J5041.3300	 Ø 3.3 J5041.3300	3.3		Stainless steel
J5041.3800	 Ø 3.8 J5041.3800	3.8		
J5041.4300	 Ø 4.3 J5041.4300	4.3		
Guide System Guiding sleeves, height 3.0 mm (2 units)				
J3714.3303		3.3		Titanium alloy
J3714.3803		3.8		
J3714.4303		4.3		
Guide System Insertion Posts				
C2026.3300		3.3		Titanium alloy
C2026.3800		3.8		
C2026.4300		4.3		

Note: All Guide System drills and gingiva punches are intended for single use only!

SURGERY

GUIDE SYSTEM

TEMPLATE-GUIDED IMPLANT BED PREPARATION AND IMPLANTATION

Art. No.	Article	\varnothing in mm	Length in mm	Material
Guide System Instruments				
J3713.3300	 Guide System Template drill, for Guide System Guiding sleeve	3.3		Stainless steel
J3713.4300	 Guide System Seating tool, for Guide System Guiding sleeve	3.8	4.3	
J3716.3300	 Guide System Check-up pin, for Guide System Guiding sleeve	3.3		Stainless steel
J3716.4300	 Guide System Driver, for Guide System Implant \varnothing 3.3/3.8/4.3 mm, manual/wrench	3.8	4.3	
J5301.3300	 Guide System Check-up pin, for Guide System Guiding sleeve	3.3		Stainless steel
J5301.4300	 Guide System Driver, for Guide System Implant \varnothing 3.3/3.8/4.3 mm, with ISO shaft for angled hand piece	3.8	4.3	
J5303.4300	 Guide System Driver, for Guide System Implant \varnothing 3.3/3.8/4.3 mm, manual/wrench	3.3	3.8	Stainless steel
J5304.4300	 Guide System Driver, for Guide System Implant \varnothing 3.3/3.8/4.3 mm, with ISO shaft for angled hand piece	3.3	3.8	

SURGERY

CONELOG® HEALING CAPS

CONELOG® Healing cap, cylindrical, titanium alloy, sterile					
GH	Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	
GH 2.0	C2015.3320	C2015.3820	C2015.4320	C2015.5020	
GH 4.0	C2015.3340	C2015.3840	C2015.4340	C2015.5040	
GH 6.0		C2015.3860*	C2015.4360*	C2015.5060*	

* Suitable for bite registration.



CONELOG® Healing cap, wide body, titanium alloy, sterile					
GH	Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	
GH 4.0	C2014.3340	C2014.3840	C2014.4340	C2014.5040	
GH 6.0		C2014.3860	C2014.4360	C2014.5060	

CONELOG® Healing cap, bottleneck, titanium alloy, sterile					
GH	Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	
GH 4.0	C2011.3340	C2011.3840	C2011.4340	C2011.5040	
GH 6.0		C2011.3860	C2011.4360	C2011.5060	

Ø = implant diameter, GH = gingiva height in mm

SURGERY

ALTAPIN SET

Art. No.	Article	Material
M5600.0110	 ALTAPin set Membrane fixation system, resterilizable Content: M5100.0010 ALTAPin applicator, straight M5100.0050 ALTAPin pricker (holder) M5100.0070 ALTAPin membrane fixator M5200.0055 ALTAPin pricker, insert (2 units) M1000.0050 ALTAPin magazine, 7 titanium pins (4 units)	Plastic/ Titanium alloy/ Stainless steel
M5600.0210	ALTAPin tray (without content)	Plastic
M5100.0010	 ALTAPin applicator, straight incl. activator	Stainless steel
M5100.0030	 ALTAPin applicator, angled 90° incl. activator	Stainless steel
M5200.0010	 ALTAPin applicator, straight, work element incl. activator	Stainless steel
M5100.0050	 ALTAPin pricker	Stainless steel
M5100.0070	 ALTAPin membrane fixator	Stainless steel
M5100.0100	 ALTAPin surgery mallet	Stainless steel
M5500.0050	 ALTAPin single patient drill, ISO shaft	Stainless steel
M5200.0055	 ALTAPin pricker, insert	Stainless steel
M1000.0050	 ALTAPin magazine 7 titanium pins, sterile (1 unit)	Titanium alloy
M1000.0100	 ALTAPin magazine 7 titanium pins, sterile (3 units)	Titanium alloy

SURGERY

SURGICAL INSTRUMENTS

Art. No.	Article	Dimensions in mm	Material
J5300.2028	 Paralleling pin SCREW-LINE, with depth marks	Ø 1.7-2.8/ 2.0	Titanium alloy
J5002.0005	 Drill extension ISO shaft	26.6	Stainless steel
J5322.0010	 Tap adapter, short for tap SCREW-LINE	18.0	Stainless steel
J5322.0011	 Tap adapter, long for tap SCREW-LINE	23.0	Stainless steel
J5300.0007	 Driver for screw implants with ISO shaft for angled hand piece	27.5	Stainless steel
J5300.0008	 Driver, short for screw implants, manual/wrench, with borehole for screwdriver, hex, long	18.0	Stainless steel
J5300.0009	 Driver, long for screw implants, manual/wrench	27.0	Stainless steel

SURGERY

SURGICAL INSTRUMENTS

Art. No.	Article	Dimensions in mm	Material
J5300.0010	Cardanic driver (30°) for screw implants, adjustable length	–	Stainless steel
			
J5300.0030	PickUp instrument holder for carrying implants	–	Stainless steel
			
J5002.0010	Adapter ISO shaft for angled hand piece	21.0	Stainless steel
			
J5302.0010	Holding key for insertion post	–	Stainless steel
			
C5302.3310		3.3	Stainless steel
C5302.4310		3.8	
C5302.4310		4.3	
J5302.3300		3.3	Titanium alloy
J5302.3800		3.8	
J5302.4300		4.3	
J5002.0012	Cleaning needle for drills with internal irrigation		Stainless steel
			
J5002.0020	Cleaning cannula for drills with internal irrigation		Stainless steel
			

SURGERY

SURGICAL-PROSTHETIC INSTRUMENTS

Art. No.	Article	Length in mm	Material
J5330.8600	Prosthetic set Content: J5320.1030 Torque wrench Screwdriver: J5316.0501 hex, short, manual/wrench J5316.0502 hex, long, manual/wrench J5316.0503 hex, long, ISO shaft J5316.0504 hex, short, ISO shaft		Plastic/ Stainless steel
J5330.8500	Prosthetic tray (without content)		Plastic
J5316.0510		Screwdriver hex, extra short, manual/wrench	14.5 Stainless steel
J5316.0501		Screwdriver hex, short, manual/wrench	22.5 Stainless steel
J5316.0502		Screwdriver hex, long, manual/wrench	30.3 Stainless steel
J5316.0504		Screwdriver hex, short, ISO shaft	18.0 Stainless steel
J5316.0503		Screwdriver hex, long, ISO shaft	26.0 Stainless steel
J5320.1030	 Torque wrench with continuous torque adjustment until maximal 30 Ncm	—	Stainless steel
C5300.1601	 CONELOG® Disconnector for CONELOG® Abutments Ø 3.3/3.8/4.3 mm Thread M 1.6		Stainless steel
C5300.2001	 CONELOG® Disconnector for CONELOG® Abutments Ø 5.0 mm Thread M 2.0		Stainless steel

SURGERY
PROSTHETIC INSTRUMENTS

Art. No.	Article	Length in mm	Material
J5300.0011	 Driver for ball abutment, manual/wrench	18.3	Stainless steel
J5315.0005	 Screwdriver Activator, for ball abutment matrix CM Dalbo®-Plus	—	Stainless steel
C5300.0020	 Driver for CONELOG® Bar abutment	18.6	Stainless steel
C5300.0027	 Driver for impression post and healing caps for CONELOG® Bar abutment	19.1	Stainless steel
J2253.0001	 Driver for Locator®, manual/wrench	24.3	Stainless steel
J2253.0002	 Locator® Instrument threepart	83.0	Stainless steel
J2253.0003	 ZEST ANCHORS INC. 800-262-2310 25° 20° 15° 10° 0° 10° 15° 20° 25° Locator® Angle measurement guide	—	Stainless steel
J2253.0004	 Locator® Parallel post (4 units)	—	Nylon

CONELOG® PROSTHETICS

CONTENT

CONELOG® PROSTHETICS	37
PROSTHETICS OVERVIEW	38
TORQUE WRENCH SETTINGS	42
IMPRESSION TAKING, BITE REGISTRATION, FABRICATION OF THE PLASTER MODEL	44
IMPRESSION TAKING, OPEN/CLOSED TRAY	44
BITE REGISTRATION	44
FABRICATION OF THE PLASTER MODEL	45
CROWN AND BRIDGE RESTORATIONS	45
CONELOG® TEMPORARY ABUTMENT	45
CONELOG® ESTHOMIC® ABUTMENTS	45
CONELOG® UNIVERSAL ABUTMENT	47
CONELOG® GOLD-PLASTIC ABUTMENT	47
CONELOG® LOGFIT® ABUTMENTS	47
ACCESSORIES FOR CONELOG® ABUTMENTS	48
CROWN, BRIDGE AND HYBRID RESTORATIONS	49
CONELOG® TITANIUM BASES CAD/CAM	49
ACCESSORIES FOR CAD/CAM PROSTHETICS	50
CONELOG® VARIO SR ABUTMENTS	51
ACCESSORIES FOR VARIO SR	53
HYBRID RESTORATIONS	54
CONELOG® BAR ABUTMENTS	54
CONELOG® BALL ABUTMENTS	56
CONELOG® LOCATOR® ABUTMENTS	58
CONELOG® UNIVERSAL ABUTMENT	61
CONELOG® TELESCOPE ABUTMENT	61
ACCESSORIES FOR CONELOG® ABUTMENTS	61
INSTRUMENTS FOR DENTAL TECHNICIANS	62
CONELOG® SELECTION ABUTMENTS	63

CONELOG® PROSTHETICS

PROSTHETICS OVERVIEW

PROSTHETICS

IMPLANT IMPRESSION TAKING

CONELOG® Impression post,
open tray



	Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm
C2121.3300				
C2121.3800				
C2121.4300				
C2121.5000				

CONELOG® Impression post,
closed tray



C2110.3300				
C2110.3800				
C2110.4300				
C2110.5000				

Impression cap
for impression post, closed tray



J2111.3300				
J2111.3800				
J2111.4300				
J2111.5000				

CONELOG® Bite registration posts



C2140.3300				
C2140.3800				
C2140.4300				
C2140.5000				

FABRICATION OF THE PLASTER MODEL

CONELOG® Lab analogs



C3010.3300				
C3010.3800				
C3010.4300				
C3010.5000				

CONELOG® ABUTMENTS FOR CROWN AND BRIDGE RESTORATIONS

**CONELOG® Temporary
abutments, titanium alloy**



C2239.3300				
C2239.3800				
C2239.4300				
C2239.5000				

**CONELOG® Esthomic®
Abutments, straight**



GH				
1.5-2.5				
3.0-4.5				
–	C2226.3815	C2226.4315	C2226.5015	
–	C2226.3830	C2226.4330	C2226.5030	

**CONELOG® Esthomic®
Abutments,
15° angled, type A**



GH				
1.5-2.5				
3.0-4.5				
–	C2227.3815	C2227.4315	C2227.5015	
–	C2227.3830	C2227.4330	C2227.5030	

**CONELOG® Esthomic®
Abutments,
15° angled, type B**



GH				
1.5-2.5				
3.0-4.5				
–	C2228.3815	C2228.4315	C2228.5015	
–	C2228.3830	C2228.4330	C2228.5030	

GH = gingiva height

CONELOG® PROSTHETICS

PROSTHETICS OVERVIEW

		\varnothing 3.3 mm	\varnothing 3.8 mm	\varnothing 4.3 mm	\varnothing 5.0 mm	
CONELOG® Esthetic® Abutments, 20° angled, type A	GH 1.5 - 2.5 3.0 - 4.5	—	C2231.3815 C2231.3830	C2231.4315 C2231.4330	C2231.5015 C2231.5030	
CONELOG® Esthetic® Abutments, 20° angled, type B	GH 1.5 - 2.5 3.0 - 4.5	—	C2232.3815 C2232.3830	C2232.4315 C2232.4330	C2232.5015 C2232.5030	
CONELOG® Esthetic® Abutments, Inset	GH 2.0 - 3.3	C2235.3320	C2235.3820	C2235.4320	C2235.5020	
CONELOG® Universal abutments		C2211.3300	C2211.3800	C2211.4300	C2211.5000	
CONELOG® Gold-plastic abutments		C2246.3300	C2246.3800	C2246.4300	C2246.5000	
CONELOG® Logfit® Abutments	GH 1.0 GH 2.5	—	C2550.3810 C2550.3825	C2550.4310 C2550.4325	C2550.5010 C2550.5025	
CONELOG® Telescope abutments, for double crown restorations		—	C2212.3800	C2212.4300	C2212.5000	
CONELOG® CAD/CAM PROSTHETICS FOR CROWNS, BRIDGES AND HYBRID RESTORATIONS	CONELOG® Titanium bases CAD/CAM	GH 0.8 GH 2.0	C2244.3308 C2244.3320	C2244.3808 C2244.3820	C2244.4308 C2244.4320	C2244.5008 C2244.5020
	CONELOG® Titanium bases CAD/CAM	GH 0.8 GH 2.0	C2242.3308* C2242.3320*	C2242.3808* C2242.3820*	C2242.4308* C2242.4320*	C2242.5008* C2242.5020*

GH = gingiva height

* Note: Article numbers will be changed from C2244. to C2242.
as of September 2014

CONELOG® PROSTHETICS

PROSTHETICS OVERVIEW

		\varnothing 3.3 mm	\varnothing 3.8 mm	\varnothing 4.3 mm	\varnothing 5.0 mm
CONELOG® VARIO SR ABUTMENTS FOR CROWNS, BRIDGES AND HYBRID RESTORATIONS OCCLUSALLY SCREW-RETAINED	CONELOG® Vario SR Abutments, straight				
			GH 1.0	— C2560.3808 C2560.4308 C2560.5008	
	CONELOG® Vario SR Abutments, 20° angled				
			GH 3.5-1.9 4.0-1.8	— C2561.3800 C2561.4300 — C2561.5000	
	CONELOG® Vario SR Abutments, 30° angled				
			GH 3.5-1.1 4.5-1.3	— C2562.3800 C2562.4300 — C2562.5000	
	Vario SR Impression caps, closed tray (without rotation securing device)				
				— J2565.4300 J2565.6000	
	Vario SR Impression caps, open tray (without rotation securing device)				
				— J2566.4300 J2566.6000	
	Vario SR analog				
				— J2567.4300 J2567.6000	

GH = gingiva height

CONELOG® PROSTHETICS

PROSTHETICS OVERVIEW

		\varnothing 3.3 mm	\varnothing 3.8 mm	\varnothing 4.3 mm	\varnothing 5.0 mm
CONELOG® ABUTMENTS FOR HYBRID RESTORATIONS	CONELOG® Bar abutments				
		GH 1.0 GH 2.5 GH 4.0	C2255.3310 C2255.3325 –	C2255.3810 C2255.3825 C2255.3840	C2255.4310 C2255.4325 C2255.4340
					C2255.5010 C2255.5025 C2255.5040
	CONELOG® Ball abutment sets, incl. male part and matrix CM Dalbo®-Plus, red dublication aid/spacer, stabilizing ring and ball abutment analog				
		GH 1.5 GH 3.0 GH 4.5	C2250.3315 C2250.3330 –	C2250.3815 C2250.3830 C2250.3845	C2250.4315 C2250.4330 C2250.4345
					C2250.5015 C2250.5030 C2250.5045
	CONELOG® Ball abutments, male part				
		GH 1.5 GH 3.0 GH 4.5	C2249.3315 C2249.3330 –	C2249.3815 C2249.3830 C2249.3845	C2249.4315 C2249.4330 C2249.4345
					C2249.5015 C2249.5030 C2249.5045
	CONELOG® Locator® Abutments				
		GH 1.0 GH 2.0 GH 3.0 GH 4.0 GH 5.0	C2253.3310 C2253.3320 C2253.3330 C2253.3340 –	C2253.3810 C2253.3820 C2253.3830 C2253.3840 C2253.3850	C2253.4310 C2253.4320 C2253.4330 C2253.4340 C2253.4350
					C2253.5010 C2253.5020 C2253.5030 C2253.5040 C2253.5050
CAD/CAM PROSTHETICS	CONELOG® Scanbodies				
		Ø 3.8 / 4.3	C2600.3310	C2600.4310	C2600.5010
	CONELOG® ScanPosts				
		\$	C2620.3306	C2620.3806	C2620.4306
					C2620.5006

GH = gingiva height

CONELOG® PROSTHETICS

TORQUE WRENCH SETTINGS

FOR FINAL FIXATION OF SURGICAL COMPONENTS AND SUPRASTRUCTURES

Article



J5320.1030

CONELOG® Implant cover screws



CONELOG® Healing caps

cylindrical
wide body
bottleneck




Vario SR protection caps



CONELOG® Lab screws

C4006.1601
C4016.1601*
(Ø 3.3/3.8/4.3 mm)

C4006.2001
C4016.2001*
(Ø 5.0 mm)



CONELOG® Vario SR lab screws

C4008.1600
(Ø 3.8/4.3 mm)

C4008.2000
(Ø 5.0 mm)



* available from fall 2014

CONELOG®

Temporary abutments
(Optional: Torque after
finished healing phase
20 Ncm)



CONELOG® Abutment screws


C4005.1601
(Ø 3.3/3.8/4.3 mm)

C4005.2001
(Ø 5.0 mm)



C4015.1601**
(Ø 3.3/3.8/4.3 mm)

C4015.2001**
(Ø 5.0 mm)



CONELOG® Vario SR abutment screws

C4007.1600
(Ø 3.8/4.3 mm)

C4007.2000
(Ø 5.0 mm)



** available from September 2014

CONELOG® Esthomic® Abutments, straight



CONELOG® Esthomic® Abutments, angled 15°/20°, type A and B

CONELOG® Esthomic® Abutments, inset

CONELOG® Gold-plastic abutments



CONELOG® Universal abutments

CONELOG® Telescope abutment



CONELOG® Logfit® Abutments

CONELOG® Titanium bases CAD/CAM

(from left to right)



CONELOG® Vario SR abutments straight, 20°/30° angled



Vario SR prosthetic screw

J4005.2004
(Ø 3.8/4.3/5.0 mm)



Vario SR

plastic copings/titanium
caps/bases for bar



Torque
tightened by hand

20 Ncm

15 Ncm

All screws have to be tightened after at least 5 minutes with the corresponding torque!

CONELOG® PROSTHETICS

TORQUE WRENCH SETTINGS

FOR FINAL FIXATION OF SURGICAL COMPONENTS AND SUPRASTRUCTURES

Article	camlog		Torque in mm	\varnothing
CONELOG® Bar abutments		 J5320.1030	C5300.0020 (\varnothing 3.3/3.8/4.3/5.0 mm)	20 Ncm 30 Ncm 4.3 5.0
Healing caps for bar abutments			C5300.0027 (\varnothing 3.3/3.8/4.3/5.0 mm)	tightened by hand
Prosthetic screw			J4005.1602 (\varnothing 3.3/3.8/4.3/5.0 mm)	15 Ncm
Bases for bar abutments Titanium bonding base for bar abutments		 J5316.0510 (hex, extra short) J5316.0501 (hex, short) J5316.0502 (hex, long)		
CONELOG® Ball abutments		 J5300.0011	20 Ncm 30 Ncm 4.3 5.0	3.3
CONELOG® Locator® Abutments		 J2253.0001	20 Ncm 30 Ncm 4.3 5.0	3.3
CONELOG® Scanbodies (\varnothing 3.3/3.8/4.3/5.0 mm)		 J5316.0501 (hex, short) J5316.0502 (hex, long)	tightened by hand	
CONELOG® ScanPosts (\varnothing 3.3/3.8/4.3/5.0 mm)				

All screws have to be tightened after at least 5 minutes with the corresponding torque!

CONELOG® PROSTHETICS

IMPRESSION TAKING, BITE REGISTRATION, FABRICATION OF THE PLASTER MODEL

IMPRESSION TAKING

Art. No.	Article	\varnothing in mm	Dimensions in mm	Material
CONELOG® Impression posts, open tray, incl. fixing screw (The fixing screw can be shortened extra-oral by 3 mm with a screwdriver, hex.)				
C2121.3300		3.3	19.3/22.3	Titanium alloy
C2121.3800		3.8	19.3/22.3	
C2121.4300		4.3	19.3/22.3	
C2121.5000		5.0	19.3/22.3	

CONELOG® Impression posts, closed tray, incl. impression cap, bite registration cap and fixing screw

C2110.3300		3.3	–	Titanium alloy/ POM
C2110.3800		3.8	–	
C2110.4300		4.3	–	
C2110.5000		5.0	–	

Impression caps for Impression posts, closed tray, (5 units)

J2111.3300		3.3	–	POM
J2111.3800		3.8	–	
J2111.4300		4.3	–	
J2111.5000		5.0	–	

BITE REGISTRATION

CONELOG® Bite registration posts, incl. fixing screw and bite registration cap

C2140.3300		3.3	–	Titanium alloy/ POM
C2140.3800		3.8	–	
C2140.4300		4.3	–	
C2140.5000		5.0	–	

Bite registration caps (5 units)

J2112.3300		3.3	–	POM
J2112.3800		3.8	–	
J2112.4300		4.3	–	
J2112.5000		5.0	–	

POM: Polyoxyethylene
H: height in mm

CONELOG® PROSTHETICS

IMPRESSION TAKING, BITE REGISTRATION, FABRICATION OF THE PLASTER MODEL

FABRICATION OF THE PLASTER MODEL

Art. No.	Article	Ø in mm	GH in mm	Material


CONELOG® Lab analogs for CONELOG® Abutments

C3010.3300		3.3	–	Titanium alloy
C3010.3800		3.8	–	
C3010.4300		4.3	–	
C3010.5000		5.0	–	

SINGLE-CROWN RESTORATIONS

TEMPORARY RESTORATIONS

CONELOG® Temporary abutments, preparable, incl. CONELOG® Abutment screw

C2239.3300	 H 11.0	3.3	–	Titanium alloy
C2239.3800		3.8	–	
C2239.4300		4.3	–	
C2239.5000		5.0	–	

CROWN AND BRIDGE RESTORATIONS

CONELOG® ESTHOMIC® ABUTMENTS FOR CEMENTED CROWN AND BRIDGE RESTORATIONS

CONELOG® Esthomic® Abutments, straight, preparable, incl. CONELOG® Abutment screw

C2226.3815	 H 9.7	3.8	1.5-2.5 3.0-4.5	Titanium alloy
C2226.4315		4.3	1.5-2.5 3.0-4.5	
C2226.4330		5.0	1.5-2.5 3.0-4.5	
C2226.5015		5.0	1.5-2.5 3.0-4.5	
C2226.5030				

CONELOG® Esthomic® Abutments, 15° angled, type A, preparable, incl. CONELOG® Abutment screw

C2227.3815	 H 9.5	3.8	1.5-2.5 3.0-4.5	Titanium alloy
C2227.4315		4.3	1.5-2.5 3.0-4.5	
C2227.4330		5.0	1.5-2.5 3.0-4.5	
C2227.5015				
C2227.5030				

GH: gingiva height

POM: Polyoxyethylene

H: height in mm


CONELOG® PROSTHETICS

CROWN AND BRIDGE RESTORATIONS


CONELOG® ESTHOMIC® ABUTMENTS FOR CEMENTED CROWN AND BRIDGE RESTORATIONS

Art. No.	Article	\varnothing in mm	GH in mm	Material


CONELOG® Esthomic® Abutments, 15° angled, type B, preparable, incl. CONELOG® Abutment screw

C2228.3815	H 9.5		3.8	1.5-2.5 3.0-4.5	Titanium alloy
C2228.3830			4.3	1.5-2.5 3.0-4.5	
C2228.4315			5.0	1.5-2.5 3.0-4.5	
C2228.4330					
C2228.5015					
C2228.5030					

CONELOG® Esthomic® Abutments, 20° angled, type A, preparable, incl. CONELOG® Abutment screw

C2231.3815	H 9.5		3.8	1.5-2.5 3.0-4.5	Titanium alloy
C2231.3830			4.3	1.5-2.5 3.0-4.5	
C2231.4315			5.0	1.5-2.5 3.0-4.5	
C2231.4330					
C2231.5015					
C2231.5030					

CONELOG® Esthomic® Abutments, 20° angled, type B, preparable, incl. CONELOG® Abutment screw


C2232.3815	H 9.5		3.8	1.5-2.5 3.0-4.5	Titanium alloy
C2232.3830			4.3	1.5-2.5 3.0-4.5	
C2232.4315			5.0	1.5-2.5 3.0-4.5	
C2232.4330					
C2232.5015					
C2232.5030					

CONELOG® Esthomic® Abutments, Inset, preparable, incl. CONELOG® Abutment screw

C2235.3320	H 9.0		3.3	2.0-3.3	Titanium alloy
C2235.3820			3.8	2.0-3.3	
C2235.4320			4.3	2.0-3.3	
C2235.5020			5.0	2.0-3.3	

GH: gingiva height H: height in mm

Description of CONELOG® Esthomic® Abutments type A and B



CONELOG® PROSTHETICS

CROWN AND BRIDGE RESTORATIONS

CONELOG® UNIVERSAL ABUTMENT FOR CEMENTED CROWN AND BRIDGE RESTORATIONS

Art. No.	Article	Ø in mm	GH in mm	Material
C2211.3300	H 11.0	3.3*	—	Titanium alloy
C2211.3800		3.8	—	
C2211.4300		4.3	—	
C2211.5000		5.0	—	

CONELOG® Universal abutments, preparable, incl. CONELOG® Abutment screw

C2211.3300	H 11.0	3.3*	—	Titanium alloy
C2211.3800		3.8	—	
C2211.4300		4.3	—	
C2211.5000		5.0	—	

CONELOG® GOLD-PLASTIC ABUTMENTS FOR CROWN AND BRIDGE RESTORATIONS

CONELOG® Gold-plastic abutments, cast-on, incl. CONELOG® Abutment screw

C2246.3300	H 11.7	0,31 g**	3.3*	—	Cast-on gold alloy/POM
C2246.3800		0,36 g**	3.8	—	
C2246.4300		0,36 g**	4.3	—	
C2246.5000		0,55 g**	5.0	—	

* only for crown restorations in the region of the upper lateral and lower lateral and central incisors (Ø 3.3 mm not for double crown restorations)

** Noble metal weight approx.

CONELOG® LOGFIT® PROSTHETIC SYSTEM FOR CEMENTED CROWN AND BRIDGE RESTORATIONS

CONELOG® Logfit® Abutments, incl. CONELOG® Abutment screw

C2550.3810	H 6.0	3.8	1.0	Titanium alloy
C2550.3825		2.5		
C2550.4310		4.3	1.0	
C2550.4325			2.5	
C2550.5010		5.0	1.0	
C2550.5025			2.5	

Logfit® Impression caps

J2551.4300	H 12.0	3.8	—	POM
J2551.4300		4.3	—	
J2551.6000		5.0	—	

Logfit® Analogs

J2552.4300		3.8	—	Titanium alloy
J2552.4300		4.3	—	
J2552.6000		5.0	—	

GH: gingiva height

POM: Polyoxymethylene

H: height in mm

CONELOG® PROSTHETICS

CROWN AND BRIDGE RESTORATIONS

CONELOG® LOGFIT® PROSTHETIC SYSTEM FOR CEMENTED CROWN AND BRIDGE RESTORATIONS

Art. No.	Article	\varnothing in mm	GH in mm	Material
Logfit® Plastic copings, for crowns (with rotation securing device), burn-out				
J2553.4302	H 5.8	3.8	–	POM
J2553.4302		4.3	–	
J2553.6002		5.0	–	

Logfit® Plastic copings, for bridges (without rotation securing device), burn-out

J2553.4301	H 5.8	3.8	–	POM
J2553.4301		4.3	–	
J2553.6001		5.0	–	

POM: Polyoxyethylene
H: height in mm

ACCESSORIES FOR CONELOG® ABUTMENTS

CONELOG® Abutment screws for the definitive screwing of abutments into the implant

C4005.1601		CONELOG® Abutment screw, hex	3.3	Thread	Titanium alloy
			3.8	Thread M 1.6	
			4.3		
C4005.2001			5.0	Thread M 2.0	

CONELOG® Lab screws for the fixation of abutments on the working model

C4006.1601		CONELOG® Lab screw, hex, brown anodized	3.3	Thread	Titanium alloy
			3.8	Thread M 1.6	
			4.3		
C4006.2001			5.0	Thread M 2.0	

Note: Lab screws be used on working models only.

CONELOG® PROSTHETICS

CROWN, BRIDGE AND HYBRID RESTORATIONS

CAD/CAM PROSTHETICS

Art. No.	Article	\varnothing in mm	GH in mm	Material
CONELOG® Titanium bases CAD/CAM, bonding base for individual CAD/CAM fabricated dental prosthesis, incl. CONELOG® Abutment screw for Titanium base CAD/CAM and CONELOG® Bonding aid (POM)				
C2244.3308**	C2242.3308**/*	H 4.7	3.3	0.8 2.0
C2244.3320**	C2242.3320**/*	H 4.7	3.8	0.8 2.0
C2244.3808	C2242.3808*		4.3	0.8 2.0
C2244.3820	C2242.3820*		5.0	0.8 2.0
C2244.4308	C2242.4308*			
C2244.4320	C2242.4320*			
C2244.5008	C2242.5008*			
C2244.5020	C2242.5020*			

**only for crown restorations in the region of the upper lateral and lower lateral and central incisors.

* available from September 2014. The new CONELOG® Titanium base CAD/CAM with modified screw channel is compatible with Sirona inCoris ZI meso blocks and replaces the existing Titanium base CAD/CAM with Art. No. C2244.xxxx. The new Titanium base can only be used with the new CONELOG® Abutment screw for Titanium base CAD/CAM, dark purple anodized.

Note: To achieve high user friendliness and precision, the geometries of the CONELOG® Titanium bases CAD/CAM are already stored in some dental CAD systems as a part library.

For more information see www.camlog.com/products/conelog/digital-technology.

CONELOG® Modeling aids for CONELOG® Titanium base CAD/CAM, burn-out, for fabricating mesostructures and crowns

C2244.3302	C2242.3302*	H 11.0	3.3	POM
C2244.3802	C2242.3802*		3.8	
C2244.4302	C2242.4302*		4.3	
C2244.5002	C2242.5002*		5.0	

* available from fall 2014. The new CONELOG® Modeling aid with modified screw channel replaces the existing modeling aid with Art.-No. C2244.xxxx and can only be used with the new CONELOG® Abutment screw for Titanium base CAD/CAM, dark purple anodized.

GH: gingiva height POM: Polyoxymethylene
H: height in mm

CONELOG® PROSTHETICS

CROWN, BRIDGE AND HYBRID RESTORATIONS

CAD/CAM PROSTHETICS

Art. No.	Article	∅ in mm	Dimensions in mm	Material
CONELOG® Scanbodies, for optical, 3-dimensional localization of CONELOG® Implants in the mouth or CONELOG® Lab analogs in the working model, incl. CONELOG® Abutment screw, sterile				
C2600.3310	H 10	3.3	—	PEEK
C2600.4310		3.8	—	
C2600.5010		4.3	—	
		5.0	—	

Note: Some scan systems require the use of system-innate scan abutments.

For information on the compatibility of the CONELOG® Scanbody with suitable dental CAD systems see www.camlog.com/products/conelog/digital-technology.

CONELOG® ScanPosts for Sirona Scanbody, for digital recording of the CONELOG® Implant or lab analog position, incl. CONELOG® Abutment screw

C2620.3306	H 10.2	3.3	—	Titanium alloy
C2620.3806		3.8	—	
C2620.4306		4.3	—	
C2620.5006		5.0	—	

ACCESSORIES FOR CAD/CAM PROSTHETICS

CONELOG® Abutment screws for Titanium base CAD/CAM, hex, for definitive screw connection into the implant

C4005.1601		C4015.1601*		3.3	Thread M 1.6	Titanium alloy
C4005.2001		C4015.2001*		3.8	—	
				4.3	—	
				5.0	Thread M 2.0	

* available from September 2014. The new modified and dark purple anodized CONELOG® Abutment screw for Titanium base CAD/CAM replaces the existing abutment screw for the new Titanium base CAD/CAM.

CONELOG® Lab screws for Titanium base CAD/CAM, hex, for the fixation on the working model

C4006.1601	C4016.1601*		3.3	Thread M 1.6	Titanium alloy
C4006.2001	C4016.2001*		3.8	—	
			4.3	—	
			5.0	Thread M 2.0	

* available from fall 2014. The new modified CONELOG® Lab screw for Titanium base CAD/CAM replaces the existing lab screw for the new Titanium base CAD/CAM.

Note: Lab screws be used on working models only.

GH: gingiva height


PEEK: Polyetheretherketon

H: height in mm

CONELOG® PROSTHETICS

CROWN, BRIDGE AND HYBRID RESTORATIONS

VARIO SR PROSTHETIC COMPONENTS FOR OCCLUSALLY SCREW-RETAINED RESTORATIONS

Art. No.	Article	\varnothing in mm	GH in mm	Material
CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw, sterile				
C2560.3808		3.8	1.0	Titanium alloy
C2560.4308		4.3	1.0	
C2560.5008		5.0	1.0	
CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw, sterile				
C2561.3800		H 6.5	3.8	3.5-1.9
C2561.4300		H 6.5	4.3	3.5-1.9
C2561.5000		H 7.0	5.0	4.0-1.8
CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw, sterile				
C2562.3800		H 6.0	3.8	3.5-1.1
C2562.4300		H 6.0	4.3	3.5-1.1
C2562.5000		H 7.0	5.0	4.5-1.3
Impression caps, open and closed tray				
J2566.4300		Vario SR impression cap, open tray, incl. fixing screw	3.8	—
J2566.6000			4.3	Titanium alloy
J2566.6000			5.0	
J2565.4300		Vario SR impression cap, closed tray	3.8	—
J2565.6000			4.3	Titanium alloy
J2565.6000			5.0	
Fabrication of the plaster model				
J2567.4300		Vario SR analog	3.8	—
J2567.6000			4.3	Titanium alloy
J2567.6000			5.0	

GH: gingiva height
H: height in mm

CONELOG® PROSTHETICS

CROWN, BRIDGE AND HYBRID RESTORATIONS

VARIO SR PROSTHETIC COMPONENTS FOR OCCLUSALLY SCREW-RETAINED RESTORATIONS

Art. No.	Article	Ø in mm	GH in mm	Material
Cast crown, bridge and bar restorations				
J2563.4302	 Vario SR plastic coping, crown, with triple-surface antirotation, burn-out	3.8 4.3	—	POM
J2563.6002		5.0	—	
J2563.4301	 Vario SR plastic coping, bridge, burn-out (2 units)	3.8 4.3	—	POM
J2563.6001		5.0	—	
J2563.4303	 Vario SR plastic coping, crown, with single-surface antirotation, for angled Vario SR abutments, burn-out	3.8 4.3	—	POM
J2563.6003		5.0	—	
Temporary or final bridge restorations				
J2564.4301	 Vario SR titanium cap, bridge, incl. Vario SR prosthetic screw	3.8 4.3	—	Titanium alloy
J2564.6001		5.0	—	
Laser-welded bar restorations				
J2570.4300	 Vario SR base for bar, laser-weldable, incl. Vario SR prosthetic screw	3.8 4.3	—	Titanium Grade 4
J2570.6000		5.0	—	
Accessories for Vario SR				
J2568.4300	 Vario SR protection cap, sterile	3.8 4.3	—	Titanium alloy
J2568.6000		5.0	—	

GH: gingiva height

H: height in mm

POM: Polyoxymethylene

CONELOG® PROSTHETICS

CROWN, BRIDGE AND HYBRID RESTORATIONS

VARIO SR PROSTHETIC COMPONENTS FOR OCCLUSALLY SCREW-RETAINED RESTORATIONS

Art. No.	Article	∅ in mm	Dimensions in mm	Material
Accessories for Vario SR				
J2569.0001		3.8	4.3	5.0
	Vario SR aligning tool 20°, for insertion post, sterile	—	—	Copolycarbonate
J2569.0002		3.8	4.3	5.0
	Vario SR aligning tool 30°, for insertion post, sterile	—	—	Copolycarbonate
C4007.1600		3.8	4.3	Thread M 1.6
	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight	—	—	Titanium alloy
C4007.2000		5.0	—	Thread M 2.0
	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight	—	—	Titanium alloy
C4005.1601		3.8	4.3	Thread M 1.6
	CONELOG® Abutment screw, hex, for CONELOG® Vario SR abutment, 20°/30° angled,	—	—	Titanium alloy
C4005.2001		5.0	—	Thread M 2.0
	CONELOG® Abutment screw, hex, for CONELOG® Vario SR abutment, 20°/30° angled,	—	—	Titanium alloy
C4008.1600		3.8	4.3	Thread M 1.6
	CONELOG® Vario SR lab screw, hex, brown anodized, for CONELOG® Vario SR abutment, straight	—	—	Titanium alloy
C4008.2000		5.0	—	Thread M 2.0
	CONELOG® Vario SR lab screw, hex, brown anodized, for CONELOG® Vario SR abutment, straight	—	—	Titanium alloy
C4006.1601		3.8	4.3	Thread M 1.6
	CONELOG® Lab screw, hex, brown anodized, for CONELOG® Vario SR abutment, 20°/30° angled	—	—	Titanium alloy
C4006.2001		5.0	—	Thread M 2.0
	CONELOG® Lab screw, hex, brown anodized, for CONELOG® Vario SR abutment, 20°/30° angled	—	—	Titanium alloy
J4005.2004		3.8	4.3	5.0
	Vario SR prosthetic screw, hex, yellow anodized	—	—	Titanium alloy

Note: Lab screws be used on working models only.


CONELOG® PROSTHETICS

HYBRID RESTORATIONS

BAR RESTORATIONS

Art. No.	Article	\varnothing in mm	GH in mm	Material
CONELOG® Bar Abutments				
C2255.3310		3.3	1.0	Titanium alloy
C2255.3325			2.5	
C2255.3810			1.0	
C2255.3825		3.8	2.5	
C2255.3840			4.0	
C2255.4310			1.0	
C2255.4325		4.3	2.5	
C2255.4340			4.0	
C2255.5010			1.0	
C2255.5025		5.0	2.5	
C2255.5040			4.0	
Impression post for CONELOG® Bar abutment				
C2130.4300	H	H 6.5	3.3	Titanium alloy
		H 6.5	3.8	
		H 6.5	4.3	
		H 7.0	5.0	
Healing cap for CONELOG® Bar abutment				
C2030.4300	H 3.7	3.3	—	Titanium alloy
		3.8		
		4.3		
		5.0		
Soldering aid/bar-lab analog for CONELOG® Bar abutment				
C3020.4300		3.3	—	Stainless steel
		3.8		
		4.3		
		5.0		
Screw, hex, for bar abutment, screw for bases for bar abutments for impression taking, open tray and soldering				
J4005.1610		3.3	Thread	Stainless steel
		3.8	M 1.6	
		4.3		
		5.0		

GH: gingiva height
H: height in mm

Art. No.	Article	\varnothing in mm	Dimensions	Material
Base for CONELOG® Bar abutment, burn-out				
C2257.3300		3.3 3.8 4.3 5.0	—	POM
Base for CONELOG® Bar abutment, cast-on				
C2263.4300		Noble metal weight approx. 0,48 g	3.3 3.8 4.3 5.0	Cast-on gold alloy/POM
Base for CONELOG® Bar abutment, solderable				
C2258.4300		3.3 3.8 4.3 5.0	—	Solderable gold alloy
Base for CONELOG® Bar abutment, laser-weldable				
C2262.4300		3.3 3.8 4.3 5.0	—	Titanium Grade 4
Titanium bonding base for CONELOG® Bar abutment, Passive-Fit				
C2260.3300		3.3 3.8 4.3 5.0	—	Titanium alloy

POM: Polyoxyethylene
H: height in mm

CONELOG® PROSTHETICS

HYBRID RESTORATIONS

BAR RESTORATIONS

Art. No.	Article	\varnothing in mm	GH in mm	Material
Sleeve for titanium bonding base, burn-out, Passive-Fit				
C2261.3300	H14.0	3.3 3.8 4.3 5.0	—	POM
Holding pin for CONELOG® Bar abutment/CONELOG® Bar lab analog, yellow, (4 units)				
C2256.1600		3.3 3.8 4.3 5.0	Thread M 1.6	POM
Prosthetic screw, for bar abutments, hex, (also for Passive-Fit)				
J4005.1602		3.3 3.8 4.3 5.0	Thread M 1.6	Titanium alloy

BALL ABUTMENT ANCHORING SYSTEM

CONELOG® Ball abutment sets

incl. male part and matrix CM Dalbo®-Plus, red duplication aid/spacer, stabilizing ring and ball abutment analog

C2250.3315		3.3	1.5	Titanium alloy/
C2250.3330			3.0	Titanium Grade 4/
C2250.3815			1.5	Gold alloy/Brass/
C2250.3830			3.0	Plastic
C2250.3845			4.5	
C2250.4315			1.5	
C2250.4330			3.0	
C2250.4345			4.5	
C2250.5015		5.0	1.5	
C2250.5030			3.0	
C2250.5045			4.5	

GH: gingiva height
H: height in mm

POM: Polyoxyethylene

BALL ABUTMENT ANCHORING SYSTEM

Art. No.	Article	\varnothing in mm	GH in mm	Material
CONELOG® Ball abutments, male part, incl. stabilizing ring				
C2249.3315		3.3	1.5	Titanium alloy/
C2249.3330			3.0	Plastic
C2249.3815		3.8	1.5	
C2249.3830			3.0	
C2249.3845			4.5	
C2249.4315		4.3	1.5	
C2249.4330			3.0	
C2249.4345			4.5	
C2249.5015		5.0	1.5	
C2249.5030			3.0	
C2249.5045			4.5	

Matrix CM Dalbo®-Plus, for ball abutment, incl. lamella retention insert

J2250.0005		3.3	–	Titanium Grade 4/
		3.8		Gold alloy
		4.3		
		5.0		

Lamella retention insert, for matrix CM Dalbo®-Plus

J2250.0007		3.3	–	Gold alloy
		3.8		
		4.3		
		5.0		

Ball abutment analogs, incl. stabilizing ring

C3015.3300		3.3	–	Brass/
		3.8		Plastic
		4.3		
C3015.5000		5.0	–	

GH: gingiva height

Dalbo®-Plus is a registered trademark of Cendres + Métaux SA, Biel, Switzerland.

CONELOG® PROSTHETICS

HYBRID RESTORATIONS

LOCATOR® ANCHORING SYSTEM

Art. No.	Article	\varnothing in mm	GH in mm	Material
CONELOG® Locator® Abutments				
C2253.3310			1.0	Titanium alloy/TiN
C2253.3320			2.0	
C2253.3330		3.3	3.0	
C2253.3340			4.0	
C2253.3810			1.0	
C2253.3820			2.0	
C2253.3830		3.8	3.0	
C2253.3840			4.0	
C2253.3850			5.0	
C2253.4310			1.0	
C2253.4320			2.0	
C2253.4330		4.3	3.0	
C2253.4340			4.0	
C2253.4350			5.0	
C2253.5010			1.0	
C2253.5020			2.0	
C2253.5030		5.0	3.0	
C2253.5040			4.0	
C2253.5050			5.0	

Impression taking

J2253.0200		Locator® Impression cap (4 units)	3.3	–	Aluminum/ Polyethylene
			3.8		
			4.3		
			5.0		

Fabrication of the plaster model

J2253.0340		Locator® Analog, (4 units)	3.3	–	Aluminum
			3.8		
			4.3		
			5.0		

GH: gingiva height TiN: Titanium nitride-coated

Art. No.	Article	∅ in mm	Dimensions in mm	Material
Locator® Male processing packages				
J2253.0102		Locator® Male processing packages (2 units), Content per package: 1 Titanium housing with processing replacement male, 1 Block out spacer white, 1 Replacement male clear, 1 Replacement male pink, 1 Replacement male blue	3.3 3.8 4.3 5.0	– Titanium alloy/ Polyethylene/ Teflon/Nylon
J2253.0112		Locator® Male processing packages for extended range* (2 units), Content per package: 1 Titanium housing with processing replacement male, 1 Block out spacer white, 1 Replacement male green, 1 Replacement male orange, 1 Replacement male red	3.8 4.3 5.0	– Titanium alloy/ Polyethylene/ Teflon/Nylon

* nonlicensed for implant diameter 3.3 mm

Locator® Accessories

J2253.0401		Locator® Block out spacer (20 units)	3.3 3.8 4.3 5.0	– Teflon
J2253.0402		Locator® Processing replacement male (4 units)	3.3 3.8 4.3 5.0	– Polyethylene

CONELOG® PROSTHETICS

HYBRID RESTORATIONS

LOCATOR® ANCHORING SYSTEM

Art. No.	Article	∅ in mm	Dimensions in mm	Material
Locator® Replacement males				
J2253.1005		Locator® Replacement male, clear, STRONG, Div.: 0°-10° (4 units)	3.3 3.8 4.3 5.0	– Nylon
J2253.1003		Locator® Replacement male, pink, MEDIUM, Div.: 0°-10° (4 units)	3.3 3.8 4.3 5.0	– Nylon
J2253.1002		Locator® Replacement male, blue, LIGHT, Div.: 0°-10° (4 units)	3.3 3.8 4.3 5.0	– Nylon
Locator® Replacement males for extended range				
J2253.2004		Locator® Replacement male for extended range*, green, STRONG, Div.: 10°-20° (4 units)	3.8 4.3 5.0	– Nylon
J2253.2003		Locator® Replacement male for extended range*, orange, MEDIUM, Div.: 10°-20° (4 units)	3.8 4.3 5.0	– Nylon
J2253.2002		Locator® Replacement male for extended range*, red, LIGHT, Div.: 10°-20° (4 units)	3.8 4.3 5.0	– Nylon
J2253.2000		Locator® Replacement male for extended range*, gray, NO RETENTION, Div.: 0°-20° (4 units)	3.8 4.3 5.0	– Nylon

* nonlicensed for implant diameter 3.3 mm

Manufacturer Locator®:
Zest Anchors LLC
Escondido, CA 92029, USA

Locator® is a registered trademark of Zest Anchors LLC

CONELOG® PROSTHETICS

HYBRID RESTORATIONS

DOUBLE CROWN RESTORATIONS WITH CONELOG® UNIVERSAL AND TELESCOPE ABUTMENTS

Art. No.	Article	\varnothing in mm	Dimensions in mm	Material
CONELOG® Universal abutments for double crown restorations, preparable, incl. CONELOG® Abutment screw				
C2211.3800	H11.0	3.8	—	Titanium alloy
C2211.4300		4.3	—	
C2211.5000		5.0	—	

CONELOG® Telescope abutments for double crown restorations, preparable, incl. CONELOG® Abutment screw

C2212.3800	H12.0	3.8	—	Titanium alloy
C2212.4300		4.3	—	
C2212.5000		5.0	—	

H: height in mm

ACCESSORIES FOR CONELOG® ABUTMENTS

CONELOG® Abutment screws for the definitive screwing of abutments into the implant

C4005.1601		CONELOG® Abutment screw, hex	3.3	Thread M 1.6	Titanium alloy
			3.8		
			4.3		
C4005.2001			5.0	Thread M 2.0	

CONELOG® Lab screws for the fixation of abutments on the working model

C4006.1601		CONELOG® Lab screw, hex, brown anodized	3.3	Thread M 1.6	Titanium alloy
			3.8		
			4.3		
C4006.2001			5.0	Thread M 2.0	

Note: Lab screws be used on working models only.

CONELOG® PROSTHETICS

INSTRUMENTS FOR DENTAL TECHNICIANS

Art. No.	Article	\varnothing in mm	Material
Universal holders			
C3709.0010	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets for \varnothing 3.3/3.8/4.3/5.0 mm, for grinding CONELOG® Abutments	—	Stainless steel/ Titanium alloy
			
J3709.0015	Universal holder	—	Stainless steel
			
CONELOG® Abutment collets , for universal holder, for grinding CONELOG® Abutments			
C3709.3300		3.3	Titanium alloy
C3709.3800		3.8	
C3709.4300		4.3	
C3709.5000		5.0	
Reworking reamers			
J3711.0010	Reworking reamer, for base for bar abutment, burn-out, plane surface/cone seat	3.3 3.8 4.3 5.0	Stainless steel
			
J3711.0020	Reworking reamer, for base for bar abutment, burn-out, screw seat	3.3 3.8 4.3 5.0	Stainless steel
			

CONELOG® PROSTHETICS

CONELOG® SELECTION ABUTMENTS

Art. No.	Article	\varnothing in mm	GH in mm	Material
C8011.1000				
CONELOG® Abutment selection kit, content (2 units each):				
CONELOG® Esthomic® Selection abutment, straight	3.8	1.5-2.5	POM	
CONELOG® Esthomic® Selection abutment, straight	3.8	3.0-4.5		
CONELOG® Esthomic® Selection abutment, straight	4.3	1.5-2.5		
CONELOG® Esthomic® Selection abutment, straight	4.3	3.0-4.5		
CONELOG® Esthomic® Selection abutment, straight	5.0	1.5-2.5		
CONELOG® Esthomic® Selection abutment, straight	5.0	3.0-4.5		
CONELOG® Esthomic® Selection abutment, 15° angled, type A	3.8	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 15° angled, type A	4.3	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 15° angled, type A	5.0	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 15° angled, type B	3.8	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 15° angled, type B	4.3	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 15° angled, type B	5.0	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 20° angled, type A	3.8	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 20° angled, type A	4.3	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 20° angled, type A	5.0	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 20° angled, type B	3.8	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 20° angled, type B	4.3	1.5-2.5		
CONELOG® Esthomic® Selection abutment, 20° angled, type B	5.0	1.5-2.5		
CONELOG® Vario SR selection abutment, straight	3.8	1.0		
CONELOG® Vario SR selection abutment, straight	4.3	1.0		
CONELOG® Vario SR selection abutment, straight	5.0	1.0		
CONELOG® Vario SR selection abutment, 20° angled	3.8	3.5-1.9		
CONELOG® Vario SR selection abutment, 20° angled	4.3	3.5-1.9		
CONELOG® Vario SR selection abutment, 20° angled	5.0	4.0-1.8		
CONELOG® Vario SR selection abutment, 30° angled	3.8	3.5-1.1		
CONELOG® Vario SR selection abutment, 30° angled	4.3	3.5-1.1		
CONELOG® Vario SR selection abutment, 30° angled	5.0	4.5-1.3		

GH: gingiva height POM: Polyoxymethylene

Attention, do not use selection abutments on patients!

AUXILIARY ARTICLES

CONTENT

AUXILIARY ARTICLES	65
IMPLANTS FOR PRACTICE	66
DEMONSTRATION MODELS	66
MACRO MODEL	66
CAMLOG LITERATURE	67

AUXILIARY ARTICLES

IMPLANTS FOR PRACTICE

Art. No.	Article	\varnothing in mm	Length in mm	Material
C1069.3813	 CONELOG® SCREW-LINE Implant for practice, incl. insertion post and CONELOG® Cover screw, yellow anodized	3.8	13	Titanium alloy
C1069.4313-01	 CONELOG® SCREW-LINE Implant for practice, incl. insertion post and CONELOG® Cover screw, red anodized (10 units)	4.3	13	Titanium alloy

Attention, do not use implants for practice on patients!

DEMONSTRATION MODELS

C8070.1020		CONELOG® Demonstration model, acrylic glass, upper jaw, 4 CONELOG® SCREW-LINE Implants, 4 x \varnothing 4.3 mm	Acrylic glass/ titanium
C8050.1040		CONELOG® Demonstration model, acrylic glass, lower jaw, 4 CONELOG® SCREW-LINE Implants, 4 x \varnothing 4.3 mm	Acrylic glass/ titanium
J8070.2050		Edentulous mandible incl. mounting plate	Plastic

MACRO MODEL

C8010.1010		CONELOG® Macro model, Scale 3:1 Content: – 1 CONELOG® SCREW-LINE Implant – 1 CONELOG® Esthetic Abutment, straight – 1 CONELOG® Abutment screw, hex – 1 Screwdriver, hex – 1 Premolar, suitable for CONELOG® Esthetic Abutment, straight – 1 Acrylic socket	Plastic/ Stainless steel
------------	---	--	-----------------------------

AUXILIARY ARTICLE

CAMLOG LITERATURE

Art. No.

Article




Patient brochure

Questions and answers to dental implants



Implant pass

Patient-specific documentation of implant restoration
Packaging units: 10 units



Patient advice sheets

Set á 4 sheets, A4



Presentation folder

A4, laminated



Poster

Format: 50 x 70 cm

AUXILIARY ARTICLE

CAMLOG LITERATURE

Art. No.

Article



Appointment pad
50 sheets/pad, A7
Packaging units: 5 units

B2012.0100



IMPLANT PROSTHETICS

DVD compendium

Four teams – their concepts and solutions, Volume 1-4

A. Kirsch, K. L. Ackermann, G. Neuendorff,

A. Happe, A. Nolte, S. Wolfart, V. Weber, F. Beuer,

M. Stimmelmayr, J. Schweiger

2012 Quintessence Publishing Co, Ltd

**INFORMATION
CONTENT**

INFORMATION

69

MATERIALS

70

INFORMATION MATERIALS

Titanium Grade 4

Properties (ASTM F67):

Chemical structure (in %):	O	≤ 0.4
	Fe	≤ 0.5
	C	≤ 0.08
	N	≤ 0.05
	H	≤ 0.015
	Ti	Rest
Mechanical properties:	Tensile strength	$\geq 550 \text{ MPa}$
	Elongation at break	$\geq 15 \%$

Titanium alloy Ti6Al4V ELI

Properties (ASTM F136):

Chemical structure (in %):	Al	5.5 - 6.5
	V	3.5 - 4.5
	Fe	≤ 0.25
	C	≤ 0.08
	N	≤ 0.05
	O	≤ 0.13
	H	≤ 0.012
	Ti	Rest
Mechanical properties:	Tensile strength	$\geq 860 \text{ MPa}$
	Elongation at break	$\geq 10 \%$

Cast-on gold alloy CONELOG® Gold-plastic abutment

Properties:

Chemical structure (in %):	Au	60
	Pd	20
	Pt	19
	Ir	1
Physical properties:	Melting range	1400 - 1490 °C
	Density	17.5 g/cm³
	E-Modul	136 GPa
	Coefficient of thermal expansion (25-500°C)	11.9 µm/m·°C
	Coefficient of thermal expansion (25-600°C)	12.2 µm/m·°C
	Color	white
Mechanical properties:		drawn
	Hardness HV5	> 215
	Tensile strength (Rm)	> 750 MPa
	0.2% Elongation limit (Rp 0.2%)	> 650 MPa
	Elongation at break	> 2%

INFORMATION MATERIALS

**Cast-on gold alloy
Base for CONELOG® Bar
abutment**

Properties:

Chemical structure (in %):	Au	60
	Pt	19
	Pd	20
	Ir	1
Mechanical properties:	Density	17.5 g/cm ³
	Color	white
	Liquidus	1490 °C
	Solidus	1400 °C
Coefficient of thermal expansion (25-500°C)	12.5 µm/m·°C	
Coefficient of thermal expansion (25-600°C)	12.6 µm/m·°C	
E-Modul	136.000 GPa	
	hardened	
	700 °C / 30 min.	
Hardness HV5	210	
0.2 % Elongation limit	450-570 MPa	
Elongation at break	min 10 %	
Tensile strength MPa	530-650	

Solderable gold alloy Base for CONELOG® Bar abutment

Properties:

Chemical structure (in %):	Au	70.00
	Pt	8.50
	Ag	13.40
	Pd	-
	Cu	7.50
	Zn	0.50
	Ir	0.10
	Rh	-
	Ru	-
Mechanical properties:	Color	yellow
	Melting range	895 - 1010 °C
	Hardness	
	annealed HV5	170
	hardened HV5	295
	self hardened HV5	280

INDEX

INDEX

73

ALPHABETICAL

74

ARTICLE NUMBERS

81

INDEX

ALPHABETICAL

Article	Art. No.	Page	Article	Art. No.	Page
A					
Adapter for screw implants, long, for CONELOG® SCREW-LINE Implants			CONELOG® Abutment screws, hex, anodized dark purple		
Ø 3.3 mm	C5302.3310	33	Ø 3.3/3.8/4.3 mm, thread M 1.6 mm	C4015.1601	50
Ø 3.8/4.3 mm	C5302.4310	33	Ø 5.0 mm, thread M 2.0 mm	C4015.2001	50
Adapter ISO shaft for angled hand piece.....	J5002.0010	33	CONELOG® Abutment selection kit	C8011.1000	63
ALTApin applicator, angled 90°	M5100.0030	31	CONELOG® Ball abutments, male part, incl. stabilizing ring		
ALTApin applicator, straight	M5100.0010	31	Ø 3.3 mm, GH 1.5 mm	C2249.3315	57
ALTApin applicator, straight, work element	M5200.0010	31	Ø 3.3 mm, GH 3.0 mm	C2249.3330	57
ALTApin magazine			Ø 3.8 mm, GH 1.5 mm	C2249.3815	57
7 titanium pins, sterile, 1 unit	M1000.0050	31	Ø 3.8 mm, GH 3.0 mm	C2249.3830	57
7 titanium pins, sterile, 3 units.....	M1000.0100	31	Ø 3.8 mm, GH 4.5 mm	C2249.3845	57
ALTApin membrane fixator	M5100.0070	31	Ø 4.3 mm, GH 1.5 mm	C2249.4315	57
ALTApin pricker	M5100.0050	31	Ø 4.3 mm, GH 3.0 mm	C2249.4330	57
ALTApin pricker, insert	M5200.0055	31	Ø 4.3 mm, GH 4.5 mm	C2249.4345	57
ALTApin set	M5600.0110	31	Ø 5.0 mm, GH 1.5 mm	C2249.5015	57
ALTApin single patient drill, ISO shaft	M5500.0050	31	Ø 5.0 mm, GH 3.0 mm	C2249.5030	57
ALTApin surgery mallet	M5100.0100	31	Ø 5.0 mm, GH 4.5 mm	C2249.5045	57
ALTApin tray (without content)	M5600.0210	31	CONELOG® Bar Abutments		
B			Ø 3.3 mm, GH 1.0 mm	C2255.3310	54
Ball abutment analogs, incl. stabilizing ring			Ø 3.3 mm, GH 2.5 mm	C2255.3325	54
Ø 3.3/3.8/4.3 mm	C3015.3300	57	Ø 3.8 mm, GH 1.0 mm	C2255.3810	54
Ø 5.0 mm.....	C3015.5000	57	Ø 3.8 mm, GH 2.5 mm	C2255.3825	54
Base for CONELOG® Bar abutment, burn-out			Ø 3.8 mm, GH 4.0 mm	C2255.3840	54
Ø 3.3/3.8/4.3/5.0 mm	C2257.3300	55	Ø 4.3 mm, GH 1.0 mm	C2255.4310	54
Base for CONELOG® Bar abutment, cast-on			Ø 4.3 mm, GH 2.5 mm	C2255.4325	54
Ø 3.3/3.8/4.3/5.0 mm	C2263.4300	55	Ø 4.3 mm, GH 4.0 mm	C2255.4340	54
Base for CONELOG® Bar abutment, laser-weldable			Ø 5.0 mm, GH 1.0 mm	C2255.5010	54
Ø 3.3/3.8/4.3/5.0 mm	C2262.4300	55	Ø 5.0 mm, GH 2.5 mm	C2255.5025	54
Base for CONELOG® Bar abutment, solderable			Ø 5.0 mm, GH 4.0 mm	C2255.5040	54
Ø 3.3/3.8/4.3/5.0 mm	C2258.4300	55	CONELOG® Ball abutment sets, incl. male part and matrix CM Dalbo®-Plus		
Bite registration caps (5 units)			Ø 3.3 mm, GH 1.5 mm	C2250.3315	56
Ø 3.3 mm.....	J2112.3300	44	Ø 3.3 mm, GH 3.0 mm	C2250.3330	56
Ø 3.8 mm.....	J2112.3800	44	Ø 3.8 mm, GH 1.5 mm	C2250.3815	56
Ø 4.3 mm.....	J2112.4300	44	Ø 3.8 mm, GH 3.0 mm	C2250.3830	56
Ø 5.0 mm.....	J2112.5000	44	Ø 3.8 mm, GH 4.5 mm	C2250.3845	56
C			Ø 4.3 mm, GH 1.5 mm	C2250.4315	56
Cardanic driver (30°).....	J5300.0010	33	Ø 4.3 mm, GH 3.0 mm	C2250.4330	56
Cleaning cannula			Ø 4.3 mm, GH 4.5 mm	C2250.4345	56
for drills with internal irrigation	J5002.0020	33	Ø 5.0 mm, GH 1.5 mm	C2250.5015	56
Cleaning needle			Ø 5.0 mm, GH 3.0 mm	C2250.5030	56
for drills with internal irrigation	J5002.0012	33	Ø 5.0 mm, GH 4.5 mm	C2250.5045	56
CONELOG® Abutment collets, for universal holder			CONELOG® Bite registration posts		
Ø 3.3 mm.....	C3709.3300	62	Ø 3.3 mm.....	C2140.3300	44
Ø 3.8 mm.....	C3709.3800	62	Ø 3.8 mm.....	C2140.3800	44
Ø 4.3 mm.....	C3709.4300	62	Ø 4.3 mm.....	C2140.4300	44
Ø 5.0 mm.....	C3709.5000	62	Ø 5.0 mm.....	C2140.5000	44
CONELOG® Abutments screws, hex			CONELOG® Demonstration model, acrylic glass, lower jaw, 4 CONELOG® SCREW-LINE Implants		
Ø 3.3/3.8/4.3 mm, thread M 1.6 mm	C4005.1601	48, 50, 53, 61	Ø 4.3 mm.....	C8050.1040	66
Ø 5.0 mm, thread M 2.0 mm	C4005.2001	48, 50, 53, 61	CONELOG® Demonstration model, acrylic glass, upper jaw, 4 CONELOG® SCREW-LINE Implants		
CONELOG® Disconnectors for CONELOG® Abutments			Ø 4.3 mm.....	C8070.1020	66
Ø 3.3/3.8/4.3 mm, thread M 1.6 mm	C5300.1601	34	CONELOG® Esthomic® Abutments, 15° angled, type A		
Ø 5.0 mm, thread M 2.0 mm	C5300.2001	34	Ø 3.8 mm, GH 1.5-2.5 mm	C2227.3815	45
CONELOG® Esthomic® Abutments, 15° angled, type A			Ø 3.8 mm, GH 3.0-4.5 mm	C2227.3830	45
Ø 3.8 mm, GH 1.5-2.5 mm	C2227.4315	45	Ø 4.3 mm, GH 1.5-2.5 mm	C2227.4330	45
Ø 3.8 mm, GH 3.0-4.5 mm	C2227.4330	45	Ø 4.3 mm, GH 3.0-4.5 mm	C2227.5015	45
Ø 4.3 mm, GH 1.5-2.5 mm	C2227.4345	45	Ø 5.0 mm, GH 1.5-2.5 mm	C2227.5030	45
Ø 4.3 mm, GH 3.0-4.5 mm	C2227.4345	45	Ø 5.0 mm, GH 3.0-4.5 mm	C2227.5030	45

INDEX

ALPHABETICAL

Article	Art. No.	Page	Article	Art. No.	Page
CONELOG® Esthomic® Abutments, 15° angled, type B			CONELOG® Impression posts, closed tray		
Ø 3.8 mm, GH 1.5-2.5 mm	C2228.3815	46	Ø 3.3 mm.....	C2110.3300	44
Ø 3.8 mm, GH 3.0-4.5 mm	C2228.3830	46	Ø 3.8 mm.....	C2110.3800	44
Ø 4.3 mm, GH 1.5-2.5 mm	C2228.4315	46	Ø 4.3 mm.....	C2110.4300	44
Ø 4.3 mm, GH 3.0-4.5 mm	C2228.4330	46	Ø 5.0 mm.....	C2110.5000	44
Ø 5.0 mm, GH 1.5-2.5 mm	C2228.5015	46			
Ø 5.0 mm, GH 3.0-4.5 mm	C2228.5030	46	CONELOG® Impression posts, open tray		
CONELOG® Esthomic® Abutments, 20° angled, type A			Ø 3.3 mm.....	C2121.3300	44
Ø 3.8 mm, GH 1.5-2.5 mm	C2231.3815	46	Ø 3.8 mm.....	C2121.3800	44
Ø 3.8 mm, GH 3.0-4.5 mm	C2231.3830	46	Ø 4.3 mm.....	C2121.4300	44
Ø 4.3 mm, GH 1.5-2.5 mm	C2231.4315	46	Ø 5.0 mm.....	C2121.5000	44
Ø 4.3 mm, GH 3.0-4.5 mm	C2231.4330	46			
Ø 5.0 mm, GH 1.5-2.5 mm	C2231.5015	46	CONELOG® Lab analogs		
Ø 5.0 mm, GH 3.0-4.5 mm	C2231.5030	46	Ø 3.3 mm.....	C3010.3300	45
CONELOG® Esthomic® Abutments, 20° angled, type B			Ø 3.8 mm.....	C3010.3800	45
Ø 3.8 mm, GH 1.5-2.5 mm	C2232.3815	46	Ø 4.3 mm.....	C3010.4300	45
Ø 3.8 mm, GH 3.0-4.5 mm	C2232.3830	46	Ø 5.0 mm.....	C3010.5000	45
Ø 4.3 mm, GH 1.5-2.5 mm	C2232.4315	46			
Ø 4.3 mm, GH 3.0-4.5 mm	C2232.4330	46	CONELOG® Lab screws, hex, brown anodized		
Ø 5.0 mm, GH 1.5-2.5 mm	C2232.5015	46	Ø 3.3/3.8/4.3 mm, thread M 1.6 mm	C4006.1601	48, 50, 53, 61
Ø 5.0 mm, GH 3.0-4.5 mm	C2232.5030	46	Ø 5.0 mm, thread M 2.0 mm	C4006.2001	48, 50, 53, 61
CONELOG® Esthomic® Abutments, Inset					
Ø 3.3 mm	C2235.3320	46	CONELOG® Lab screws, hex		
Ø 3.8 mm	C2235.3820	46	Ø 3.3/3.8/4.3 mm, thread M 1.6 mm	C4016.1601	50
Ø 4.3 mm	C2235.4320	46	Ø 5.0 mm, thread M 2.0 mm	C4016.2001	50
Ø 5.0 mm	C2235.5020	46			
CONELOG® Esthomic® Abutments, straight			CONELOG® Locator® Abutments		
Ø 3.8 mm, GH 1.5-2.5 mm	C2226.3815	45	Ø 3.3 mm, GH 1.0 mm	C2253.3310	58
Ø 3.8 mm, GH 3.0-4.5 mm	C2226.3830	45	Ø 3.3 mm, GH 2.0 mm	C2253.3320	58
Ø 4.3 mm, GH 1.5-2.5 mm	C2226.4315	45	Ø 3.3 mm, GH 3.0 mm	C2253.3330	58
Ø 4.3 mm, GH 3.0-4.5 mm	C2226.4330	45	Ø 3.3 mm, GH 4.0 mm	C2253.3340	58
Ø 5.0 mm, GH 1.5-2.5 mm	C2226.5015	45	Ø 3.8 mm, GH 1.0 mm	C2253.3810	58
Ø 5.0 mm, GH 3.0-4.5 mm	C2226.5030	45	Ø 3.8 mm, GH 2.0 mm	C2253.3820	58
CONELOG® Gold-plastic abutments, cast-on, incl. CONELOG® Abutment screw			Ø 3.8 mm, GH 3.0 mm	C2253.3830	58
Ø 3.3 mm	C2246.3300	47	Ø 3.8 mm, GH 4.0 mm	C2253.3840	58
Ø 3.8 mm	C2246.3800	47	Ø 3.8 mm, GH 5.0 mm	C2253.3850	58
Ø 4.3 mm	C2246.4300	47	Ø 4.3 mm, GH 1.0 mm	C2253.4310	58
Ø 5.0 mm	C2246.5000	47	Ø 4.3 mm, GH 2.0 mm	C2253.4320	58
CONELOG® Healing cap, bottleneck			Ø 4.3 mm, GH 3.0 mm	C2253.4330	58
Ø 3.3 mm, GH 4.0 mm	C2011.3340	30	Ø 4.3 mm, GH 4.0 mm	C2253.4340	58
Ø 3.8 mm, GH 4.0 mm	C2011.3840	30	Ø 4.3 mm, GH 5.0 mm	C2253.4350	58
Ø 3.8 mm, GH 6.0 mm	C2011.3860	30	Ø 5.0 mm, GH 1.0 mm	C2253.5010	58
Ø 4.3 mm, GH 4.0 mm	C2011.4340	30	Ø 5.0 mm, GH 2.0 mm	C2253.5020	58
Ø 4.3 mm, GH 6.0 mm	C2011.4360	30	Ø 5.0 mm, GH 3.0 mm	C2253.5030	58
Ø 5.0 mm, GH 4.0 mm	C2011.5040	30	Ø 5.0 mm, GH 4.0 mm	C2253.5040	58
Ø 5.0 mm, GH 6.0 mm	C2011.5060	30	Ø 5.0 mm, GH 5.0 mm	C2253.5050	58
CONELOG® Healing cap, cylindrical					
Ø 3.3 mm, GH 2.0 mm	C2015.3320	30	CONELOG® Logfit® Abutments, incl. CONELOG® Abutment screw		
Ø 3.3 mm, GH 4.0 mm	C2015.3340	30	Ø 3.8 mm, GH 1.0 mm	C2550.3810	47
Ø 3.8 mm, GH 2.0 mm	C2015.3820	30	Ø 3.8 mm, GH 2.5 mm	C2550.3825	47
Ø 3.8 mm, GH 4.0 mm	C2015.3840	30	Ø 4.3 mm, GH 1.0 mm	C2550.4310	47
Ø 3.8 mm, GH 6.0 mm	C2015.3860	30	Ø 4.3 mm, GH 2.5 mm	C2550.4325	47
Ø 4.3 mm, GH 2.0 mm	C2015.4320	30	Ø 5.0 mm, GH 1.0 mm	C2550.5010	47
Ø 4.3 mm, GH 4.0 mm	C2015.4340	30	Ø 5.0 mm, GH 2.5 mm	C2550.5025	47
Ø 4.3 mm, GH 6.0 mm	C2015.4360	30			
Ø 5.0 mm, GH 2.0 mm	C2015.5020	30	CONELOG® Macro model		
Ø 5.0 mm, GH 4.0 mm	C2015.5040	30	Scale 3:1	C8010.1010	66
Ø 5.0 mm, GH 6.0 mm	C2015.5060	30			
CONELOG® Healing cap, wide body			CONELOG® Modeling aids, burn-out		
Ø 3.3 mm, GH 4.0 mm	C2014.3340	30	Ø 3.3 mm	C2244.3302	49
Ø 3.8 mm, GH 4.0 mm	C2014.3840	30	Ø 3.8 mm	C2244.3802	49
Ø 3.8 mm, GH 6.0 mm	C2014.3860	30	Ø 4.3 mm	C2244.4302	49
Ø 4.3 mm, GH 4.0 mm	C2014.4340	30	Ø 5.0 mm	C2244.5002	49
Ø 4.3 mm, GH 6.0 mm	C2014.4360	30			
Ø 5.0 mm, GH 4.0 mm	C2014.5040	30	CONELOG® Modeling aids, burn-out		
Ø 5.0 mm, GH 6.0 mm	C2014.5060	30	Ø 3.3 mm	C2242.3302	49
			Ø 3.8 mm	C2242.3802	49
			Ø 4.3 mm	C2242.4302	49
			Ø 5.0 mm	C2242.5002	49
CONELOG® Scanbodies, incl. CONELOG® Abutment screw, sterile			CONELOG® Scanbodies, incl. CONELOG® Abutment screw, sterile		
Ø 3.3 mm	C2600.3310	50	Ø 3.3 mm	C2600.3310	50
Ø 3.8/4.3 mm	C2600.4310	50	Ø 3.8 mm	C2600.4310	50
Ø 5.0 mm	C2600.5010	50			

INDEX

ALPHABETICAL

Article	Art. No.	Page	Article	Art. No.	Page																																																																																																																																																																																																																																																																																																																																																																																																																																																					
CONELOG® ScanPosts for Sirona Scanbody, incl. CONELOG® Abutment screw																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.3 mm.....	C2620.3306	50	Ø 3.3 mm/3.8/4.3/5.0 mm	C3709.0010	62																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 3.8 mm.....	C2620.3806	50	Ø 3.8 mm.....	C2560.3808	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 4.3 mm.....	C2620.4306	50	Ø 4.3 mm.....	C2560.4308	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 5.0 mm.....	C2620.5006	50	Ø 5.0 mm.....	C2560.5008	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
CONELOG® SCREW-LINE Implant for practice																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.8 mm.....	C1069.3813	66	Ø 3.8 mm/4.3 mm, thread M 1.6 mm	C4007.1600	53																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 4.3 mm.....	C1069.4313-01	66	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53																																																																																																																																																																																																																																																																																																																																																																																																																																																					
CONELOG® SCREW-LINE Implant, Promote® plus																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.3 mm, L 9 mm	C1062.3309	15	Ø 3.8 mm.....	C2561.3800	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 3.3 mm, L 11 mm	C1062.3311	15	Ø 4.3 mm.....	C2561.4300	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 3.3 mm, L 13 mm	C1062.3313	15	Ø 5.0 mm.....	C2561.5000	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 3.3 mm, L 16 mm	C1062.3316	15	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw																																																																																																																																																																																																																																																																																																																																																																																																																																																							
Ø 3.8 mm, L 7 mm	C1062.3807	15	Ø 3.8 mm, L 9 mm	C1062.3809	15	Ø 3.8 mm, L 11 mm	C1062.3811	15	Ø 3.8 mm, L 13 mm	C1062.3813	15	Ø 3.8 mm, L 16 mm	C1062.3816	15	Ø 4.3 mm, L 7 mm	C1062.4307	15	Ø 3.8 mm.....	C2562.3800	51	Ø 4.3 mm, L 9 mm	C1062.4309	15	Ø 4.3 mm.....	C2562.4300	51	Ø 4.3 mm, L 11 mm	C1062.4311	15	Ø 5.0 mm.....	C2562.5000	51	Ø 4.3 mm, L 13 mm	C1062.4313	15	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 4.3 mm, L 16 mm	C1062.4316	15	Ø 5.0 mm, L 7 mm	C1062.5007	15	Ø 5.0 mm, L 9 mm	C1062.5009	15	Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35
Ø 3.8 mm, L 9 mm	C1062.3809	15	Ø 3.8 mm, L 11 mm	C1062.3811	15	Ø 3.8 mm, L 13 mm	C1062.3813	15	Ø 3.8 mm, L 16 mm	C1062.3816	15	Ø 4.3 mm, L 7 mm	C1062.4307	15	Ø 3.8 mm.....	C2562.3800	51	Ø 4.3 mm, L 9 mm	C1062.4309	15	Ø 4.3 mm.....	C2562.4300	51	Ø 4.3 mm, L 11 mm	C1062.4311	15	Ø 5.0 mm.....	C2562.5000	51	Ø 4.3 mm, L 13 mm	C1062.4313	15	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 4.3 mm, L 16 mm	C1062.4316	15	Ø 5.0 mm, L 7 mm	C1062.5007	15	Ø 5.0 mm, L 9 mm	C1062.5009	15	Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35			
Ø 3.8 mm, L 11 mm	C1062.3811	15	Ø 3.8 mm, L 13 mm	C1062.3813	15	Ø 3.8 mm, L 16 mm	C1062.3816	15	Ø 4.3 mm, L 7 mm	C1062.4307	15	Ø 3.8 mm.....	C2562.3800	51	Ø 4.3 mm, L 9 mm	C1062.4309	15	Ø 4.3 mm.....	C2562.4300	51	Ø 4.3 mm, L 11 mm	C1062.4311	15	Ø 5.0 mm.....	C2562.5000	51	Ø 4.3 mm, L 13 mm	C1062.4313	15	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 4.3 mm, L 16 mm	C1062.4316	15	Ø 5.0 mm, L 7 mm	C1062.5007	15	Ø 5.0 mm, L 9 mm	C1062.5009	15	Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35						
Ø 3.8 mm, L 13 mm	C1062.3813	15	Ø 3.8 mm, L 16 mm	C1062.3816	15	Ø 4.3 mm, L 7 mm	C1062.4307	15	Ø 3.8 mm.....	C2562.3800	51	Ø 4.3 mm, L 9 mm	C1062.4309	15	Ø 4.3 mm.....	C2562.4300	51	Ø 4.3 mm, L 11 mm	C1062.4311	15	Ø 5.0 mm.....	C2562.5000	51	Ø 4.3 mm, L 13 mm	C1062.4313	15	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 4.3 mm, L 16 mm	C1062.4316	15	Ø 5.0 mm, L 7 mm	C1062.5007	15	Ø 5.0 mm, L 9 mm	C1062.5009	15	Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35									
Ø 3.8 mm, L 16 mm	C1062.3816	15	Ø 4.3 mm, L 7 mm	C1062.4307	15	Ø 3.8 mm.....	C2562.3800	51	Ø 4.3 mm, L 9 mm	C1062.4309	15	Ø 4.3 mm.....	C2562.4300	51	Ø 4.3 mm, L 11 mm	C1062.4311	15	Ø 5.0 mm.....	C2562.5000	51	Ø 4.3 mm, L 13 mm	C1062.4313	15	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 4.3 mm, L 16 mm	C1062.4316	15	Ø 5.0 mm, L 7 mm	C1062.5007	15	Ø 5.0 mm, L 9 mm	C1062.5009	15	Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35												
Ø 4.3 mm, L 7 mm	C1062.4307	15	Ø 3.8 mm.....	C2562.3800	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 4.3 mm, L 9 mm	C1062.4309	15	Ø 4.3 mm.....	C2562.4300	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 4.3 mm, L 11 mm	C1062.4311	15	Ø 5.0 mm.....	C2562.5000	51																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 4.3 mm, L 13 mm	C1062.4313	15	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight																																																																																																																																																																																																																																																																																																																																																																																																																																																							
Ø 4.3 mm, L 16 mm	C1062.4316	15	Ø 5.0 mm, L 7 mm	C1062.5007	15	Ø 5.0 mm, L 9 mm	C1062.5009	15	Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																										
Ø 5.0 mm, L 7 mm	C1062.5007	15	Ø 5.0 mm, L 9 mm	C1062.5009	15	Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																													
Ø 5.0 mm, L 9 mm	C1062.5009	15	Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																
Ø 5.0 mm, L 11 mm	C1062.5011	15	Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																			
Ø 5.0 mm, L 13 mm	C1062.5013	15	Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																						
Ø 5.0 mm, L 16 mm	C1062.5016	15	CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm						Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																									
CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.8 mm.....	C2212.3800	61	CT-tube for drill Ø 2.0 mm																																																																																																																																																																																																																																																																																																																																																																																																																																																							
Ø 4.3 mm.....	C2212.4300	61	Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13	CONELOG® Temporary abutments						Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm						Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																											
Ø 5.0 mm.....	C2212.5000	61	internal diameter 2.1 mm	A2002.2000	13																																																																																																																																																																																																																																																																																																																																																																																																																																																					
CONELOG® Temporary abutments																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.3 mm.....	C2239.3300	45	CT-tube for drill Ø 2.2 mm																																																																																																																																																																																																																																																																																																																																																																																																																																																							
Ø 3.8 mm.....	C2239.3800	45	Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13	Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																			
Ø 4.3 mm.....	C2239.4300	45	internal diameter 2.3 mm	A2222.2200	13																																																																																																																																																																																																																																																																																																																																																																																																																																																					
Ø 5.0 mm.....	C2239.5000	45	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable																																																																																																																																																																																																																																																																																																																																																																																																																																																							
CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.3 mm, GH 0.8 mm	C2244.3308	49	Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																											
Ø 3.3 mm, GH 2.0 mm	C2244.3320	49	Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																														
Ø 3.8 mm, GH 0.8 mm	C2244.3808	49	Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																	
Ø 3.8 mm, GH 2.0 mm	C2244.3820	49	Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																				
Ø 4.3 mm, GH 0.8 mm	C2244.4308	49	Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																							
Ø 4.3 mm, GH 2.0 mm	C2244.4320	49	Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																										
Ø 5.0 mm, GH 0.8 mm	C2244.5008	49	Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18	CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw						Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																													
Ø 5.0 mm, GH 2.0 mm	C2244.5020	49	Ø 3.3 mm.....	J5015.3300	18																																																																																																																																																																																																																																																																																																																																																																																																																																																					
CONELOG® Titanium bases CAD/CAM, incl. CONELOG® Abutment screw																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.3 mm, GH 0.8 mm	C2242.3308	49	Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																												
Ø 3.3 mm, GH 2.0 mm	C2242.3320	49	Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																															
Ø 3.8 mm, GH 0.8 mm	C2242.3808	49	Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																		
Ø 3.8 mm, GH 2.0 mm	C2242.3820	49	Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																					
Ø 4.3 mm, GH 0.8 mm	C2242.4308	49	Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																								
Ø 4.3 mm, GH 2.0 mm	C2242.4320	49	Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																											
Ø 5.0 mm, GH 0.8 mm	C2242.5008	49	Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18	CONELOG® Universal abutments						Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																														
Ø 5.0 mm, GH 2.0 mm	C2242.5020	49	Ø 3.8 mm.....	J5015.3800	18																																																																																																																																																																																																																																																																																																																																																																																																																																																					
CONELOG® Universal abutments																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.3 mm.....	C2211.3300	47	Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																													
Ø 3.8 mm.....	C2211.3800	47, 61	Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																
Ø 4.3 mm.....	C2211.4300	47, 61	Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18	CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets						Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																			
Ø 5.0 mm.....	C2211.5000	47, 61	Ø 4.3 mm.....	J5015.4300	18																																																																																																																																																																																																																																																																																																																																																																																																																																																					
CONELOG® Universal holder, incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.3/3.8/4.3/5.0 mm	C3709.0010	62	Ø 5.0 mm.....	C5015.5000	18	CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw						Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench						Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																		
Ø 5.0 mm.....	C5015.5000	18																																																																																																																																																																																																																																																																																																																																																																																																																																																								
CONELOG® Vario SR abutments, straight, incl. CONELOG® Vario SR abutment screw																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.8 mm.....	C2560.3808	51	Driver for ball abutment, manual/wrench																																																																																																																																																																																																																																																																																																																																																																																																																																																							
Ø 4.3 mm.....	C2560.4308	51	Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																							
Ø 5.0 mm.....	C2560.5008	51	CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																										
CONELOG® Vario SR abutment screw, hex, for CONELOG® Vario SR abutment, straight																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.8/4.3 mm, thread M 1.6 mm	C4007.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																			
Ø 5.0 mm, thread M 2.0 mm	C4007.2000	53	CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench						Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																						
CONELOG® Vario SR abutments, 20° angled, incl. CONELOG® Abutment screw																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.8 mm.....	C2561.3800	51	Driver for cone abutment, manual/wrench																																																																																																																																																																																																																																																																																																																																																																																																																																																							
Ø 4.3 mm.....	C2561.4300	51	Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																								
Ø 5.0 mm.....	C2561.5000	51	CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw						Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench						Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																											
CONELOG® Vario SR abutments, 30° angled, incl. CONELOG® Abutment screw																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.8 mm.....	C2562.3800	51	Driver for core abutment, manual/wrench																																																																																																																																																																																																																																																																																																																																																																																																																																																							
Ø 4.3 mm.....	C2562.4300	51	Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																													
Ø 5.0 mm.....	C2562.5000	51	CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight						Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																
CONELOG® Vario SR Lab screws, hex, for CONELOG® Vario SR abutment, straight																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.8/4.3 mm, thread M 1.6 mm	C4008.1600	53	Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																									
Ø 5.0 mm, thread M 2.0 mm	C4008.2000	53	CT-tube for drill Ø 2.0 mm						internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench						CT-tube for drill Ø 2.2 mm						internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																												
CT-tube for drill Ø 2.0 mm																																																																																																																																																																																																																																																																																																																																																																																																																																																										
internal diameter 2.1 mm	A2002.2000	13	Driver for core abutment, manual/wrench																																																																																																																																																																																																																																																																																																																																																																																																																																																							
CT-tube for drill Ø 2.2 mm																																																																																																																																																																																																																																																																																																																																																																																																																																																										
internal diameter 2.3 mm	A2222.2200	13	Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable						Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench						Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																				
Depth stop for form drills SCREW-LINE, with reduced coil, resterilizable																																																																																																																																																																																																																																																																																																																																																																																																																																																										
Ø 3.3 mm.....	J5015.3300	18	Driver for core abutment, manual/wrench																																																																																																																																																																																																																																																																																																																																																																																																																																																							
Ø 3.8 mm.....	J5015.3800	18	Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																						
Ø 4.3 mm.....	J5015.4300	18	Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																									
Ø 5.0 mm.....	J5015.5000	18	Depth stop SCREW-LINE						for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench						with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																												
Depth stop SCREW-LINE																																																																																																																																																																																																																																																																																																																																																																																																																																																										
for pilot drill (J5051.2000) and pre-drill (J5051.2800)			Driver for ball abutment, manual/wrench																																																																																																																																																																																																																																																																																																																																																																																																																																																							
with reduced coil, resterilizable			L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																														
L 7 mm.....	J5015.0007	17	L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																	
L 9 mm.....	J5015.0009	17	L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																				
L 11 mm.....	J5015.0011	17	L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																							
L 13 mm.....	J5015.0013	17	Drill extension, ISO shaft						L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																										
Drill extension, ISO shaft																																																																																																																																																																																																																																																																																																																																																																																																																																																										
L 26.6 mm.....	J5002.0005	32	Drill for CT-tube						(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																			
Drill for CT-tube																																																																																																																																																																																																																																																																																																																																																																																																																																																										
(for A2002.2000), Ø 2.6 mm	A2050.2600	13	(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																												
(for A2222.2200), Ø 2.8 mm	A2050.2800	13	Driver for ball abutment, manual/wrench						L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																															
Driver for ball abutment, manual/wrench																																																																																																																																																																																																																																																																																																																																																																																																																																																										
L 18.3 mm.....	J5300.0011	35	Driver for CONELOG® Bar abutment						L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																																								
Driver for CONELOG® Bar abutment																																																																																																																																																																																																																																																																																																																																																																																																																																																										
L 18.6 mm.....	C5300.0020	35	Driver for impression post and healing caps						for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																																																	
Driver for impression post and healing caps																																																																																																																																																																																																																																																																																																																																																																																																																																																										
for CONELOG® Bar abutment			L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																																																										
L 19.1 mm.....	C5300.0027	35	Driver for Locator®, manual/wrench, L 24.3 mm						manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																																																													
Driver for Locator®, manual/wrench, L 24.3 mm																																																																																																																																																																																																																																																																																																																																																																																																																																																										
manual/wrench			J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																																																																						
J2253.0001	35																																																																																																																																																																																																																																																																																																																																																																																																																																																									

INDEX

ALPHABETICAL

Article	Art. No.	Page	Article	Art. No.	Page
Driver for screw implants with ISO shaft for angled hand piece.. J5300.0007		Guide System CONELOG® SCREW-LINE Implant, Promote® plus			
Driver, long for screw implants, manual/wrench, L 27.0 mm..... J5300.0009	32	Ø 3.8 mm, L 16 mm C1063.3816 24			
Driver, short for screw implants, manual/wrench, L 18.0 mm..... J5300.0008	32	Ø 4.3 mm, L 7 mm C1063.4307 24			
		Ø 4.3 mm, L 9 mm C1063.4309 24			
		Ø 4.3 mm, L 11 mm C1063.4311 24			
		Ø 4.3 mm, L 13 mm C1063.4313 24			
		Ø 4.3 mm, L 16 mm C1063.4316 24			
E			Guide System Driver, manual/wrench		
Edentulous mandible incl. mounting plate J8070.2050	66	Ø 3.3/3.8/4.3 mm J5303.4300 29			
F			Guide System Driver, with ISO shaft for angled hand piece		
Form drill SCREW-LINE, Cortical bone		Ø 3.3/3.8/4.3 mm J5304.4300 29			
Ø 3.3 mm J5053.3316 19		Guide System Form drill, SCREW-LINE			
Ø 3.8 mm J5053.3816 19		Ø 3.3 mm, L 16 mm J5046.3316 25			
Ø 4.3 mm J5053.4316 19		Ø 3.8 mm, L 16 mm J5046.3816 26			
Ø 5.0 mm J5053.5016 19		Ø 4.3 mm, L 16 mm J5046.4316 26			
Form drill SCREW-LINE, with reduced coil		Guide System Form drill, SCREW-LINE, Cortical bone			
Ø 3.3 mm, L 9 mm J5052.3309 18		Ø 3.3 mm, L 9 mm J5048.3309 27			
Ø 3.3 mm, L 11 mm J5052.3311 18		Ø 3.3 mm, L 11 mm J5048.3311 27			
Ø 3.3 mm, L 13 mm J5052.3313 18		Ø 3.3 mm, L 13 mm J5048.3313 27			
Ø 3.3 mm, L 16 mm J5052.3316 18		Ø 3.3 mm, L 16 mm J5048.3316 27			
Ø 3.8 mm, L 7 mm J5052.3807 18		Ø 3.8 mm, L 7 mm J5048.3807 27			
Ø 3.8 mm, L 9 mm J5052.3809 18		Ø 3.8 mm, L 9 mm J5048.3809 27			
Ø 3.8 mm, L 11 mm J5052.3811 18		Ø 3.8 mm, L 11 mm J5048.3811 27			
Ø 3.8 mm, L 13 mm J5052.3813 18		Ø 3.8 mm, L 13 mm J5048.3813 27			
Ø 3.8 mm, L 16 mm J5052.3816 18		Ø 3.8 mm, L 16 mm J5048.3816 27			
Ø 4.3 mm, L 7 mm J5052.4307 18		Ø 4.3 mm, L 7 mm J5048.4307 27			
Ø 4.3 mm, L 9 mm J5052.4309 18		Ø 4.3 mm, L 9 mm J5048.4309 27			
Ø 4.3 mm, L 11 mm J5052.4311 18		Ø 4.3 mm, L 11 mm J5048.4311 27			
Ø 4.3 mm, L 13 mm J5052.4313 18		Ø 4.3 mm, L 13 mm J5048.4313 27			
Ø 4.3 mm, L 16 mm J5052.4316 18		Ø 4.3 mm, L 16 mm J5048.4316 27			
Ø 5.0 mm, L 7 mm J5052.5007 18		Guide System Gingiva Punches, sterile			
Ø 5.0 mm, L 9 mm J5052.5009 18		Ø 3.3 mm J5041.3300 28			
Ø 5.0 mm, L 11 mm J5052.5011 18		Ø 3.8 mm J5041.3800 28			
Ø 5.0 mm, L 13 mm J5052.5013 18		Ø 4.3 mm J5041.4300 28			
Ø 5.0 mm, L 16 mm J5052.5016 18		Guide System Guiding sleeves			
G		Ø 3.3 mm J3714.3303 28			
Guide System Check-up pin		Ø 3.8 mm J3714.3803 28			
Ø 3.3 mm J5301.3300 29		Ø 4.3 mm J3714.4303 28			
Ø 3.8/4.3 mm J5301.4300 29		Guide System Pilot drill set			
Guide System CONELOG® Insertion posts		Ø 3.3 mm, L 9mm J5043.3309 24			
Ø 3.3 mm C2026.3300 28		Ø 3.3 mm, L 11 mm J5043.3311 24			
Ø 3.8 mm C2026.3800 28		Ø 3.3 mm, L 13 mm J5043.3313 24			
Ø 4.3 mm C2026.4300 28		Ø 3.8/4.3 mm, L 7 mm J5043.4307 24			
Guide System CONELOG® SCREW-LINE Implant, Promote® plus		Ø 3.8/4.3 mm, L 9 mm J5043.4309 24			
Ø 3.3 mm, L 9 mm C1063.3309 24		Ø 3.8/4.3 mm, L 11 mm J5043.4311 25			
Ø 3.3 mm, L 11 mm C1063.3311 24		Ø 3.8/4.3 mm, L 13 mm J5043.4313 25			
Ø 3.3 mm, L 13 mm C1063.3313 24		Guide System Pilot drill			
Ø 3.3 mm, L 16 mm C1063.3316 24		Ø 3.3 mm, L 16 mm J5044.3316 24			
Ø 3.8 mm, L 7 mm C1063.3807 24		Ø 3.8/4.3 mm, L 16 mm J5044.4316 25			
Ø 3.8 mm, L 9 mm C1063.3809 24		Guide System Seating tool			
Ø 3.8 mm, L 11 mm C1063.3811 24		Ø 3.3 mm J3716.3300 29			
Ø 3.8 mm, L 13 mm C1063.3813 24		Ø 3.8/4.3 mm J3716.4300 29			
		Guide System Surgery set, SCREW-LINE			
		Ø 3.3 mm, L 9 mm J5045.3309 25			
		Ø 3.3 mm, L 11 mm J5045.3311 25			
		Ø 3.3 mm, L 13 mm J5045.3313 25			
		Ø 3.8 mm, L 7 mm J5045.3807 26			
		Ø 3.8 mm, L 9 mm J5045.3809 26			
		Ø 3.8 mm, L 11 mm J5045.3811 26			
		Ø 3.8 mm, L 13 mm J5045.3813 26			
		Ø 4.3 mm, L 7 mm J5045.4307 26			
		Ø 4.3 mm, L 9 mm J5045.4309 26			
		Ø 4.3 mm, L 11 mm J5045.4311 26			
		Ø 4.3 mm, L 13 mm J5045.4313 26			

INDEX

ALPHABETICAL

Article	Art. No.	Page
Guide System Template drill		
Ø 3.3 mm	J3713.3300	29
Ø 3.8/4.3 mm	J3713.4300	29

H

Healing cap for CONELOG® Bar abutment		
Ø 3.3/3.8/4.3/5.0 mm	C2030.4300	54
Holding key for insertion post	J5302.0010	33
Holding pin for CONELOG® Bar abutment/ CONELOG® Bar lab analog		
Ø 3.3/3.8/4.3/5.0 mm, thread M 1.6 mm	C2256.1600	56
Holding sleeve for screw implants		
Ø 3.3 mm	J5302.3300	33
Ø 3.8 mm	J5302.3800	33
Ø 4.3 mm	J5302.4300	33

IMPLANT PROSTHETICS. DVD compendium B2012.0100 68

Impression caps for Impression posts, (5 units)	
Ø 3.3 mm.....	J2111.3300 44
Ø 3.8 mm.....	J2111.3800 44
Ø 4.3 mm.....	J2111.4300 44
Ø 5.0 mm.....	J2111.5000 44

Impression post for CONELOG® Bar abutment

\emptyset 3.3/3.8/4.3/5.0 mm C2130.4300 54

Lamella retention insert

Lumena Retention Insert, for matrix CM Dalbo®-Plus		
Ø 3.3/3.8/4.3/5.0 mm	J2250.0007	57
Locator® Analog		
Ø 3.3/3.8/4.3/5.0 mm	J2253.0340	58
Locator® Angle measurement guide	J2253.0003	35
Locator® Block out spacer		
Ø 3.3/3.8/4.3/5.0 mm	J2253.0401	59
Locator® Impression cap		
Ø 3.3/3.8/4.3/5.0 mm	J2253.0200	58
Locator® Instrument, threepart		
L 83 mm.....	J2253.0002	35
Locator® Male processing packages		
Ø 3.3/3.8/4.3/5.0 mm	J2253.0102	59
Locator® Male processing packages for extended range		
Ø 3.8/4.3/5.0 mm	J2253.0112	59
Locator® Parallel post	J2253.0004	35
Locator® Processing replacement male		
Ø 3.3/3.8/4.3/5.0 mm	J2253.0402	59
Locator® Replacement male for extended range, gray		
Ø 3.8/4.3/5.0 mm	J2253.2000	60

Article	Art. No.	Page
Locator® Replacement male for extended range, green		
Ø 3.8/4.3/5.0 mm	J2253.2004	60
Locator® Replacement male for extended range, orange		
Ø 3.8/4.3/5.0 mm	J2253.2003	60
Locator® Replacement male for extended range, red		
Ø 3.8/4.3/5.0 mm	J2253.2002	60
Locator® Replacement male, blue		
Ø 3.3/3.8/4.3/5.0 mm	J2253.1002	60
Locator® Replacement male, clear		
Ø 3.3/3.8/4.3/5.0 mm	J2253.1005	60
Locator® Replacement male, pink		
Ø 3.3/3.8/4.3/5.0 mm	J2253.1003	60
Logfit® Analogs		
Ø 3.8/4.3 mm	J2552.4300	47
Ø 5.0 mm	J2552.6000	47
Logfit® Impression caps		
Ø 3.8/4.3 mm	J2551.4300	47
Ø 5.0 mm	J2551.6000	47
Logfit® Plastic copings, for bridges, burn-out		
Ø 3.8/4.3 mm	J2553.4301	48
Ø 5.0 mm	J2553.6001	48
Logfit® Plastic copings, for crowns, burn-out		
Ø 3.8/4.3 mm	J2553.4302	48
Ø 5.0 mm	J2553.6002	48

M

Matrix CM Dalbo®-Plus,
for ball abutment, incl. lamella retention insert
 \varnothing 3.3/3.8/4.3/5.0 mm J2250.0005 57

0

Osteotome SCREW-LINE, angled concave		
Ø 3.3 mm.....	J5420.3310	23
Ø 3.8 mm.....	J5420.3810	23
Ø 4.3 mm.....	J5420.4310	23
Ø 5.0 mm.....	J5420.5010	23
Osteotome SCREW-LINE, angled convex		
Ø 3.3 mm.....	J5418.3310	21
Ø 3.8 mm.....	J5418.3810	21
Ø 4.3 mm.....	J5418.4310	21
Ø 5.0 mm.....	J5418.5010	21
Osteotome SCREW-LINE, straight concave		
Ø 3.3 mm.....	J5420.3300	22
Ø 3.8 mm.....	J5420.3800	22
Ø 4.3 mm.....	J5420.4300	22
Ø 5.0 mm.....	J5420.5000	22
Osteotome SCREW-LINE, straight convex		
Ø 3.3 mm.....	J5418.3300	20
Ø 3.8 mm.....	J5418.3800	20
Ø 4.3 mm.....	J5418.4300	20
Ø 5.0 mm.....	J5418.5000	20
Osteotomy Set CAMLOG®/CONELOG® SCREW-LINE		
angled concave.....	J5420.0030	23
angled convex.....	J5418.0030	21
straight concave.....	J5420.0020	22
straight convex.....	J5418.0020	20

INDEX

ALPHABETICAL

Article	Art. No.	Page	Article	Art. No.	Page
P					
Paralleling pin					
SCREW-LINE					
Ø 1.7-2.8/2.0 mm	J5300.2028	32			
PickUp instrument	J5300.0030	33			
Pilot drill SCREW-LINE					
reduced coil, resterilizable					
Ø 2.0 mm	J5051.2000	17			
Pilot drill without coil, resterilizable					
Ø 2.0 mm	J5051.2003	17			
Pre-drill SCREW-LINE					
Ø 1.7-2.8 mm	J5051.2800	17			
Pre-Osteotome SCREW-LINE , straight concave					
Ø 1.7-2.8 mm	J5419.2800	22, 23			
Pre-Osteotome SCREW-LINE , straight convex					
Ø 1.7-2.8 mm	J5417.2800	20, 21			
Prosthetic screw, for bar abutments, hex, (also for Passive-Fit)					
Ø 3.3/3.8/4.3/5.0 mm, thread M 1.6 mm	J4005.1602	56			
Prosthetic set	J5330.8600	34			
Prosthetic tray (without content)	J5330.8500	34			
R					
Reworking reamer, for base for bar abutment, plane surface/cone seat					
Ø 3.3/3.8/4.3/5.0 mm	J3711.0010	62			
Reworking reamer, for base for bar abutment, screw seat					
Ø 3.3/3.8/4.3/5.0 mm	J3711.0020	62			
Round bur , resterilizable					
Ø 2.3 mm	J5050.2300	17			
Ø 3.5 mm	J5050.3500	17			
S					
Screw, hex, for bar abutment , screw for bases for bar abutments for impression taking, open tray and soldering					
Ø 3.3/3.8/4.3/5.0 mm, thread M 1.6 mm	J4005.1610	54			
Screwdriver Activator ,					
for ball abutment matrix CM Dalbo®-Plus	J5315.0005	35			
Screwdriver hex , long, ISO shaft.....	J5316.0503	34			
Screwdriver hex , long, manual/wrench.....	J5316.0502	34			
Screwdriver hex , extra short, manual/wrench	J5316.0510	34			
Screwdriver hex , short, ISO shaft.....	J5316.0504	34			
Screwdriver hex , short, manual/wrench.....	J5316.0501	34			
Sleeve for titanium bonding base , burn-out, Passive-Fit					
Ø 3.3/3.8/4.3/5.0 mm	C2261.3300	56			
Soldering aid/bar-lab analog					
for CONELOG® Bar abutment					
Ø 3.3/3.8/4.3/5.0 mm	C3020.4300	54			
Surgery Set CAMLOG®/CONELOG® SCREW-LINE	J5300.0057	16			
T					
Tap adapter					
for tap SCREW-LINE					
short L 18.0 mm	J5322.0010	32			
long L 23.0 mm	J5322.0011	32			
Tap SCREW-LINE , with hexagon					
Ø 3.3 mm	J5054.3309	19			
Ø 3.8 mm	J5054.3809	19			
Ø 4.3 mm	J5054.4309	19			
Ø 5.0 mm	J5054.5009	19			
Titanium bonding base for CONELOG® Bar abutment , Passive-Fit					
Ø 3.3/3.8/4.3/5.0 mm	C2260.3300	55			
Torque wrench	J5320.1030	34			
U					
Universal holder	J3709.0015	62			
V					
Vario SR analog					
Ø 3.8/4.3 mm	J2567.4300	51			
Ø 5.0 mm	J2567.6000	51			
Vario SR aligning tool 20° , for insertion post					
Ø 3.8/4.3/5.0 mm	J2569.0001	53			
Vario SR aligning tool 30° , for insertion post					
Ø 3.8/4.3/5.0 mm	J2569.0002	53			
Vario SR base for bar , laser-weldable, incl. Vario SR prosthetic screw					
Ø 3.8/4.3 mm	J2570.4300	52			
Ø 5.0 mm	J2570.6000	52			
Vario SR impression cap, closed tray					
Ø 3.8/4.3 mm	J2565.4300	51			
Ø 5.0 mm	J2565.6000	51			
Vario SR impression cap, open tray , incl. fixing screw					
Ø 3.8/4.3 mm	J2566.4300	51			
Ø 5.0 mm	J2566.6000	51			
Vario SR plastic coping, bridge , burn-out					
Ø 3.8/4.3 mm	J2563.4301	52			
Ø 5.0 mm	J2563.6001	52			
Vario SR plastic coping, crown , with triple-surface antirotation, burn-out					
Ø 3.8/4.3 mm	J2563.4302	52			
Ø 5.0 mm	J2563.6002	52			
Vario SR plastic coping, crown , with single-surface antirotation, for angled Vario SR abutments, burn-out					
Ø 3.8/4.3 mm	J2563.4303	52			
Ø 5.0 mm	J2563.6003	52			
Vario SR protection cap					
Ø 3.8/4.3 mm	J2568.4300	52			
Ø 5.0 mm	J2568.6000	52			

INDEX

ALPHABETICAL

Article	Art. No.	Page
Vario SR prosthetic screw, hex Ø 3.8/4.3/5.0 mm, thread M 2.0 mm	J4005.2004	53
Vario SR titanium cap, bridge, incl. Vario SR prosthetic screw Ø 3.8/4.3 mm	J2564.4301	52
Ø 5.0 mm.....	J2564.6001	52

X

X-Ray Planning foil 1.25:1, for CONELOG® SCREW-LINE Implants, Magnification 25%	C5300.9010	12
X-Ray Planning foil 1.4:1, for CONELOG® SCREW-LINE Implants, Magnification 40%	C5300.9011	12
X-Ray Transfer pictures 1.25:1, for CONELOG® SCREW-LINE Implants, Planning foils, self-adhesive, Magnification 25% Ø 3.3 mm.....	C5300.9080	12
Ø 3.8 mm.....	C5300.9081	12
Ø 4.3 mm.....	C5300.9082	12
Ø 5.0 mm.....	C5300.9083	12

INDEX

ARTICLE NUMBERS

Art. No.	Article	Page	Art. No.	Article	Page
	CT-tube for drill Ø 2.0 mm			CONELOG® Healing cap, cylindrical	
A2002.2000	internal diameter 2.1 mm	13	C2015.3320	Ø 3.3 mm, GH 2.0 mm	30
	Drill for CT-tube		C2015.3340	Ø 3.3 mm, GH 4.0 mm	30
A2050.2600	(for A2002.2000), Ø 2.6 mm.....	13	C2015.3820	Ø 3.8 mm, GH 2.0 mm	30
A2050.2800	(for A2222.2200), Ø 2.8 mm	13	C2015.3840	Ø 3.8 mm, GH 4.0 mm	30
	CT-tube for drill Ø 2.2 mm		C2015.3860	Ø 3.8 mm, GH 6.0 mm	30
A2222.2200	internal diameter 2.3 mm	13	C2015.4320	Ø 4.3 mm, GH 2.0 mm	30
B2012.0100	IMPLANT PROSTHETICS, DVD compendium	68	C2015.4340	Ø 4.3 mm, GH 4.0 mm	30
	Guide System CONELOG® SCREW-LINE Implant, Promote® plus		C2015.4360	Ø 4.3 mm, GH 6.0 mm	30
C1063.3309	Ø 3.3 mm, L 9 mm	24	C2015.5020	Ø 5.0 mm, GH 2.0 mm	30
C1063.3311	Ø 3.3 mm, L 11 mm	24	C2015.5040	Ø 5.0 mm, GH 4.0 mm	30
C1063.3313	Ø 3.3 mm, L 13 mm	24	C2015.5060	Ø 5.0 mm, GH 6.0 mm	30
C1063.3316	Ø 3.3 mm, L 16 mm	24		Guide System CONELOG® Insertion posts	
C1063.3807	Ø 3.8 mm, L 7 mm	24	C2026.3300	Ø 3.3 mm	28
C1063.3809	Ø 3.8 mm, L 9 mm	24	C2026.3800	Ø 3.8 mm	28
C1063.3811	Ø 3.8 mm, L 11 mm	24	C2026.4300	Ø 4.3 mm	28
C1063.3813	Ø 3.8 mm, L 13 mm	24		Healing cap for CONELOG® Bar abutment	
C1063.3816	Ø 3.8 mm, L 16 mm	24	C2030.4300	Ø 3.3/3.8/4.3/5.0 mm	54
C1063.4307	Ø 4.3 mm, L 7 mm	24		CONELOG® Impression posts, closed tray	
C1063.4309	Ø 4.3 mm, L 9 mm	24	C2110.3300	Ø 3.3 mm	44
C1063.4311	Ø 4.3 mm, L 11 mm	24	C2110.3800	Ø 3.8 mm	44
C1063.4313	Ø 4.3 mm, L 13 mm	24	C2110.4300	Ø 4.3 mm	44
C1063.4316	Ø 4.3 mm, L 16 mm	24	C2110.5000	Ø 5.0 mm	44
	CONELOG® SCREW-LINE Implant, Promote® plus			CONELOG® Impression posts, open tray	
C1062.3309	Ø 3.3 mm, L 9 mm	15	C2121.3300	Ø 3.3 mm	44
C1062.3311	Ø 3.3 mm, L 11 mm	15	C2121.3800	Ø 3.8 mm	44
C1062.3313	Ø 3.3 mm, L 13 mm	15	C2121.4300	Ø 4.3 mm	44
C1062.3316	Ø 3.3 mm, L 16 mm	15	C2121.5000	Ø 5.0 mm	44
C1062.3807	Ø 3.8 mm, L 7 mm	15		Impression post for CONELOG® Bar abutment	
C1062.3809	Ø 3.8 mm, L 9 mm	15	C2130.4300	Ø 3.3/3.8/4.3/5.0 mm	54
C1062.3811	Ø 3.8 mm, L 11 mm	15		CONELOG® Bite registration posts	
C1062.3813	Ø 3.8 mm, L 13 mm	15	C2140.3300	Ø 3.3 mm	44
C1062.3816	Ø 3.8 mm, L 16 mm	15	C2140.3800	Ø 3.8 mm	44
C1062.4307	Ø 4.3 mm, L 7 mm	15	C2140.4300	Ø 4.3 mm	44
C1062.4309	Ø 4.3 mm, L 9 mm	15	C2140.5000	Ø 5.0 mm	44
C1062.4311	Ø 4.3 mm, L 11 mm	15		CONELOG® Universal abutments, preparable, incl. CONELOG® Abutment screw	
C1062.4313	Ø 4.3 mm, L 13 mm	15	C2211.3300	Ø 3.3 mm	47
C1062.4316	Ø 4.3 mm, L 16 mm	15	C2211.3800	Ø 3.8 mm	47, 61
C1062.5007	Ø 5.0 mm, L 7 mm	15	C2211.4300	Ø 4.3 mm	47, 61
C1062.5009	Ø 5.0 mm, L 9 mm	15	C2211.5000	Ø 5.0 mm	47, 61
C1062.5011	Ø 5.0 mm, L 11 mm	15		CONELOG® Telescope abutments, preparable, incl. CONELOG® Abutment screw	
C1062.5013	Ø 5.0 mm, L 13 mm	15	C2212.3800	Ø 3.8 mm	61
C1062.5016	Ø 5.0 mm, L 16 mm	15	C2212.4300	Ø 4.3 mm	61
	CONELOG® SCREW-LINE Implant for practice		C2212.5000	Ø 5.0 mm	61
C1069.3813	Ø 3.8 mm	66		CONELOG® Esthomic® Abutments, straight	
C1069.4313-01	Ø 4.3 mm	66	C2226.3815	Ø 3.8 mm, GH 1.5-2.5 mm	45
	CONELOG® Healing cap, bottleneck		C2226.3830	Ø 3.8 mm, GH 3.0-4.5 mm	45
C2011.3340	Ø 3.3 mm, GH 4.0 mm	30	C2226.4315	Ø 4.3 mm, GH 1.5-2.5 mm	45
C2011.3840	Ø 3.8 mm, GH 4.0 mm	30	C2226.4330	Ø 4.3 mm, GH 3.0-4.5 mm	45
C2011.3860	Ø 3.8 mm, GH 6.0 mm	30	C2226.5015	Ø 5.0 mm, GH 1.5-2.5 mm	45
C2011.4340	Ø 4.3 mm, GH 4.0 mm	30	C2226.5030	Ø 5.0 mm, GH 3.0-4.5 mm	45
C2011.4360	Ø 4.3 mm, GH 6.0 mm	30		CONELOG® Esthomic® Abutments, 15° angled, type A	
C2011.5040	Ø 5.0 mm, GH 4.0 mm	30	C2227.3815	Ø 3.8 mm, GH 1.5-2.5 mm	45
C2011.5060	Ø 5.0 mm, GH 6.0 mm	30	C2227.3830	Ø 3.8 mm, GH 3.0-4.5 mm	45
	CONELOG® Healing cap, wide body		C2227.4315	Ø 4.3 mm, GH 1.5-2.5 mm	45
C2014.3340	Ø 3.3 mm, GH 4.0 mm	30	C2227.4330	Ø 4.3 mm, GH 3.0-4.5 mm	45
C2014.3840	Ø 3.8 mm, GH 4.0 mm	30	C2227.5015	Ø 5.0 mm, GH 1.5-2.5 mm	45
C2014.3860	Ø 3.8 mm, GH 6.0 mm	30	C2227.5030	Ø 5.0 mm, GH 3.0-4.5 mm	45
C2014.4340	Ø 4.3 mm, GH 4.0 mm	30		CONELOG® Esthomic® Abutments, 15° angled, type B	
C2014.4360	Ø 4.3 mm, GH 6.0 mm	30	C2228.3815	Ø 3.8 mm, GH 1.5-2.5 mm	46
C2014.5040	Ø 5.0 mm, GH 4.0 mm	30	C2228.3830	Ø 3.8 mm, GH 3.0-4.5 mm	46
C2014.5060	Ø 5.0 mm, GH 6.0 mm	30			

INDEX

ARTICLE NUMBERS

Art. No.	Article	Page	Art. No.	Article	Page
	CONELOG® Esthomic® Abutments, 15° angled, type B			CONELOG® Ball abutments, male part, incl. stabilizing ring	
C2228.4315	Ø 4.3 mm, GH 1.5-2.5 mm	46	C2249.3315	Ø 3.3 mm, GH 1.5 mm	57
C2228.4330	Ø 4.3 mm, GH 3.0-4.5 mm	46	C2249.3330	Ø 3.3 mm, GH 3.0 mm	57
C2228.5015	Ø 5.0 mm, GH 1.5-2.5 mm	46	C2249.3815	Ø 3.8 mm, GH 1.5 mm	57
C2228.5030	Ø 5.0 mm, GH 3.0-4.5 mm	46	C2249.3830	Ø 3.8 mm, GH 3.0 mm	57
	CONELOG® Esthomic® Abutments, 20° angled, type A		C2249.3845	Ø 3.8 mm, GH 4.5 mm	57
C2231.3815	Ø 3.8 mm, GH 1.5-2.5 mm	46	C2249.4315	Ø 4.3 mm, GH 1.5 mm	57
C2231.3830	Ø 3.8 mm, GH 3.0-4.5 mm	46	C2249.4330	Ø 4.3 mm, GH 3.0 mm	57
C2231.4315	Ø 4.3 mm, GH 1.5-2.5 mm	46	C2249.4345	Ø 4.3 mm, GH 4.5 mm	57
C2231.4330	Ø 4.3 mm, GH 3.0-4.5 mm	46	C2249.5015	Ø 5.0 mm, GH 1.5 mm	57
C2231.5015	Ø 5.0 mm, GH 1.5-2.5 mm	46	C2249.5030	Ø 5.0 mm, GH 3.0 mm	57
C2231.5030	Ø 5.0 mm, GH 3.0-4.5 mm	46	C2249.5045	Ø 5.0 mm, GH 4.5 mm	57
	CONELOG® Esthomic® Abutments, 20° angled, type B			CONELOG® Ball abutment sets, incl. male part and matrix CM Dalbo®-Plus	
C2232.3815	Ø 3.8 mm, GH 1.5-2.5 mm	46	C2250.3315	Ø 3.3 mm, GH 1.5 mm	56
C2232.3830	Ø 3.8 mm, GH 3.0-4.5 mm	46	C2250.3330	Ø 3.3 mm, GH 3.0 mm	56
C2232.4315	Ø 4.3 mm, GH 1.5-2.5 mm	46	C2250.3815	Ø 3.8 mm, GH 1.5 mm	56
C2232.4330	Ø 4.3 mm, GH 3.0-4.5 mm	46	C2250.3830	Ø 3.8 mm, GH 3.0 mm	56
C2232.5015	Ø 5.0 mm, GH 1.5-2.5 mm	46	C2250.3845	Ø 3.8 mm, GH 4.5 mm	56
C2232.5030	Ø 5.0 mm, GH 3.0-4.5 mm	46	C2250.4315	Ø 4.3 mm, GH 1.5 mm	56
	CONELOG® Esthomic® Abutments, Inset		C2250.4330	Ø 4.3 mm, GH 3.0 mm	56
C2235.3320	Ø 3.3 mm	46	C2250.4345	Ø 4.3 mm, GH 4.5 mm	56
C2235.3820	Ø 3.8 mm	46	C2250.5015	Ø 5.0 mm, GH 1.5 mm	56
C2235.4320	Ø 4.3 mm	46	C2250.5030	Ø 5.0 mm, GH 3.0 mm	56
C2235.5020	Ø 5.0 mm	46	C2250.5045	Ø 5.0 mm, GH 4.5 mm	56
	CONELOG® Temporary abutments			CONELOG® Locator® Abutments	
C2239.3300	Ø 3.3 mm	45	C2253.3310	Ø 3.3 mm, GH 1.0 mm	58
C2239.3800	Ø 3.8 mm	45	C2253.3320	Ø 3.3 mm, GH 2.0 mm	58
C2239.4300	Ø 4.3 mm	45	C2253.3330	Ø 3.3 mm, GH 3.0 mm	58
C2239.5000	Ø 5.0 mm	45	C2253.3340	Ø 3.3 mm, GH 4.0 mm	58
	CONELOG® Modeling aids for CONELOG® Titanium base CAD/CAM, burn-out		C2253.3810	Ø 3.8 mm, GH 1.0 mm	58
C2242.3302	Ø 3.3 mm	49	C2253.3820	Ø 3.8 mm, GH 2.0 mm	58
C2242.3802	Ø 3.8 mm	49	C2253.3830	Ø 3.8 mm, GH 3.0 mm	58
C2242.4302	Ø 4.3 mm	49	C2253.3840	Ø 3.8 mm, GH 4.0 mm	58
C2242.5002	Ø 5.0 mm	49	C2253.3850	Ø 3.8 mm, GH 5.0 mm	58
	CONELOG® Titanium bases CAD/CAM		C2253.4310	Ø 4.3 mm, GH 1.0 mm	58
C2242.3308	Ø 3.3 mm, GH 0.8 mm	49	C2253.4320	Ø 4.3 mm, GH 2.0 mm	58
C2242.3320	Ø 3.3 mm, GH 2.0 mm	49	C2253.4330	Ø 4.3 mm, GH 3.0 mm	58
C2242.3808	Ø 3.8 mm, GH 0.8 mm	49	C2253.4340	Ø 4.3 mm, GH 4.0 mm	58
C2242.3820	Ø 3.8 mm, GH 2.0 mm	49	C2253.4350	Ø 4.3 mm, GH 5.0 mm	58
C2242.4308	Ø 4.3 mm, GH 0.8 mm	49	C2253.5010	Ø 5.0 mm, GH 1.0 mm	58
C2242.4320	Ø 4.3 mm, GH 2.0 mm	49	C2253.5020	Ø 5.0 mm, GH 2.0 mm	58
C2242.5008	Ø 5.0 mm, GH 0.8 mm	49	C2253.5030	Ø 5.0 mm, GH 3.0 mm	58
C2242.5020	Ø 5.0 mm, GH 2.0 mm	49	C2253.5040	Ø 5.0 mm, GH 4.0 mm	58
	CONELOG® Modeling aids for CONELOG® Titanium base CAD/CAM, burn-out		C2253.5050	Ø 5.0 mm, GH 5.0 mm	58
C2244.3302	Ø 3.3 mm	49		CONELOG® Bar Abutments	
C2244.3802	Ø 3.8 mm	49	C2255.3310	Ø 3.3 mm, GH 1.0 mm	54
C2244.4302	Ø 4.3 mm	49	C2255.3325	Ø 3.3 mm, GH 2.5 mm	54
C2244.5002	Ø 5.0 mm	49	C2255.3810	Ø 3.8 mm, GH 1.0 mm	54
	CONELOG® Titanium bases CAD/CAM		C2255.3825	Ø 3.8 mm, GH 2.5 mm	54
C2244.3308	Ø 3.3 mm, GH 0.8 mm	49	C2255.3840	Ø 3.8 mm, GH 4.0 mm	54
C2244.3320	Ø 3.3 mm, GH 2.0 mm	49	C2255.4310	Ø 4.3 mm, GH 1.0 mm	54
C2244.3808	Ø 3.8 mm, GH 0.8 mm	49	C2255.4325	Ø 4.3 mm, GH 2.5 mm	54
C2244.3820	Ø 3.8 mm, GH 2.0 mm	49	C2255.4340	Ø 4.3 mm, GH 4.0 mm	54
C2244.4308	Ø 4.3 mm, GH 0.8 mm	49	C2255.5010	Ø 5.0 mm, GH 1.0 mm	54
C2244.4320	Ø 4.3 mm, GH 2.0 mm	49	C2255.5025	Ø 5.0 mm, GH 2.5 mm	54
C2244.5008	Ø 5.0 mm, GH 0.8 mm	49	C2255.5040	Ø 5.0 mm, GH 4.0 mm	54
C2244.5020	Ø 5.0 mm, GH 2.0 mm	49		Holding pin for CONELOG® Bar abutment/ CONELOG® Bar lab analog	
	CONELOG® Gold-plastic abutments		C2256.1600	Ø 3.3/3.8/4.3/5.0 mm, thread M 1.6 mm	56
C2246.3300	Ø 3.3 mm	47		Base for CONELOG® Bar abutment, burn-out	
C2246.3800	Ø 3.8 mm	47	C2257.3300	Ø 3.3/3.8/4.3/5.0 mm	55
C2246.4300	Ø 4.3 mm	47		Base for CONELOG® Bar abutment, solderable	
C2246.5000	Ø 5.0 mm	47	C2258.4300	Ø 3.3/3.8/4.3/5.0 mm	55

INDEX

ARTICLE NUMBERS

Art. No.	Article	Page
C2260.3300	Titanium bonding base for CONELOG® Bar abutment , Passive-Fit Ø 3.3/3.8/4.3/5.0 mm	55
C2261.3300	Sleeve for titanium bonding base , burn-out, Passive-Fit Ø 3.3/3.8/4.3/5.0 mm	56
C2262.4300	Base for CONELOG® Bar abutment , laser-weldable Ø 3.3/3.8/4.3/5.0 mm	55
C2263.4300	Base for CONELOG® Bar abutment , cast-on Ø 3.3/3.8/4.3/5.0 mm	55
C2550.3810	CONELOG® Logfit® Abutments , incl. CONELOG® Abutment screw Ø 3.8 mm, GH 1.0 mm	47
C2550.3825	Ø 3.8 mm, GH 2.5 mm	47
C2550.4310	Ø 4.3 mm, GH 1.0 mm	47
C2550.4325	Ø 4.3 mm, GH 2.5 mm	47
C2550.5010	Ø 5.0 mm, GH 1.0 mm	47
C2550.5025	Ø 5.0 mm, GH 2.5 mm	47
C2560.3808	CONELOG® Vario SR abutments , straight, incl. CONELOG® Vario SR abutment screw Ø 3.8 mm	51
C2560.4308	Ø 4.3 mm	51
C2560.5008	Ø 5.0 mm	51
C2561.3800	CONELOG® Vario SR abutments , 20° angled, incl. CONELOG® Abutment screw Ø 3.8 mm	51
C2561.4300	Ø 4.3 mm	51
C2561.5000	Ø 5.0 mm	51
C2562.3800	CONELOG® Vario SR abutments , 30° angled, incl. CONELOG® Abutment screw Ø 3.8 mm	51
C2562.4300	Ø 4.3 mm	51
C2562.5000	Ø 5.0 mm	51
C2600.3310	CONELOG® Scanbodies , incl. CONELOG® Abutment screw, sterile Ø 3.3 mm	50
C2600.4310	Ø 3.8/4.3 mm	50
C2600.5010	Ø 5.0 mm	50
C2620.3306	CONELOG® ScanPosts for Sirona Scanbody , incl. CONELOG® Abutment screw Ø 3.3 mm	50
C2620.3806	Ø 3.8 mm	50
C2620.4306	Ø 4.3 mm	50
C2620.5006	Ø 5.0 mm	50
C3010.3300	CONELOG® Lab analogs Ø 3.3 mm	45
C3010.3800	Ø 3.8 mm	45
C3010.4300	Ø 4.3 mm	45
C3010.5000	Ø 5.0 mm	45
C3015.3300	Ball abutment analogs , incl. stabilizing ring Ø 3.3/3.8/4.3 mm	57
C3015.5000	Ø 5.0 mm	57
C3020.4300	Soldering aid/bar-lab analog for CONELOG® Bar abutment Ø 3.3/3.8/4.3/5.0 mm	54
C3709.0010	CONELOG® Universal holder , incl. 2 CONELOG® Lab screws, hex, and CONELOG® Abutment collets Ø 3.3/3.8/4.3/5.0 mm	62
C3709.3300	CONELOG® Abutment collets , for universal holder Ø 3.3 mm	62
C3709.3800	Ø 3.8 mm	62
C3709.4300	Ø 4.3 mm	62
C3709.5000	Ø 5.0 mm	62

Art. No.	Article	Page
C4005.1601	CONELOG® Abutment screws, hex Ø 3.3/3.8/4.3 mm, thread M 1.6 mm.....	48, 50, 53, 61
C4005.2001	Ø 5.0 mm, thread M 2.0 mm	48, 50, 53, 61
C4006.1601	CONELOG® Lab screws, hex , brown anodized Ø 3.3/3.8/4.3 mm, thread M 1.6 mm.....	48, 50, 53, 61
C4006.2001	Ø 5.0 mm, thread M 2.0 mm	48, 50, 53, 61
C4007.1600	CONELOG® Vario SR abutment screws , hex, for CONELOG® Vario SR abutment, straight Ø 3.8/4.3 mm, thread M 1.6 mm	53
C4007.2000	Ø 5.0 mm, thread M 2.0 mm	53
C4008.1600	CONELOG® Vario SR lab screws , hex, for CONELOG® Vario SR abutment, straight Ø 3.8/4.3 mm, thread M 1.6 mm	53
C4008.2000	Ø 5.0 mm, thread M 2.0 mm	53
C4015.1601	CONELOG® Abutment screws for Titanium base CAD/CAM , hex, dark purple anodized Ø 3.3/3.8/4.3 mm, thread M 1.6 mm.....	50
C4015.2001	Ø 5.0 mm, thread M 2.0 mm	50
C5300.0020	CONELOG® Lab screws for Titanium base CAD/CAM , hex L 18.6 mm	35
C5300.0027	Driver for CONELOG® Bar abutment L 19.1 mm	35
C5300.1601	Driver for impression post and healing caps for CONELOG® Bar abutment Ø 3.3/3.8/4.3 mm, thread M 1.6 mm.....	34
C5300.2001	Ø 5.0 mm, thread M 2.0 mm	34
C5300.9010	X-Ray Planning foil 1.25:1 , for CONELOG® SCREW-LINE Implants, Magnification 25%	12
C5300.9011	X-Ray Planungsfolie 1.4:1 , for CONELOG® SCREW-LINE Implants, Magnification 40%	12
C5300.9080	X-Ray Transfer pictures 1.25:1 , Planning foils, self-adhesive, for CONELOG® SCREW-LINE Implants, Magnification 25%	12
C5300.9081	Ø 3.3 mm	12
C5300.9082	Ø 3.8 mm	12
C5300.9083	Ø 4.3 mm	12
C5300.9083	Ø 5.0 mm	12
C5302.3310	Adapter for screw implants, long , for CONELOG® SCREW-LINE Implants Ø 3.3 mm	33
C5302.4310	Ø 3.8/4.3 mm	33
C8010.1010	CONELOG® Macro model Scale 3:1	66
C8011.1000	CONELOG® Abutment selection kit	63
C8050.1040	CONELOG® Demonstration model , acrylic glass, lower jaw, 4 CONELOG® SCREW-LINE Implants Ø 4.3 mm	66
C8070.1020	CONELOG® Demonstration model , acrylic glass, upper jaw, 4 CONELOG® SCREW-LINE Implants Ø 4.3 mm	66

INDEX

ARTICLE NUMBERS

Art. No.	Article	Page	Art. No.	Article	Page
	Impression caps for Impression posts, (5 units)			Vario SR plastic coping, bridge, burn-out	
J2111.3300	Ø 3.3 mm	44	J2563.4301	Ø 3.8/4.3 mm	52
J2111.3800	Ø 3.8 mm	44		Vario SR plastic coping, crown, with triple-surface antirotation, burn-out	
J2111.4300	Ø 4.3 mm	44	J2563.4302	Ø 3.8/4.3 mm	52
J2111.5000	Ø 5.0 mm	44		Vario SR plastic coping, crown, with single-surface antirotation, for angled Vario SR abutments, burn-out	
	Bite registration caps (5 units)		J2563.4303	Ø 3.8/4.3 mm	52
J2112.3300	Ø 3.3 mm	44		Vario SR plastic coping, bridge, burn-out	
J2112.3800	Ø 3.8 mm	44	J2563.6001	Ø 5.0 mm	52
J2112.4300	Ø 4.3 mm	44		Vario SR plastic coping, crown, with triple-surface antirotation, burn-out	
J2112.5000	Ø 5.0 mm	44	J2563.6002	Ø 5.0 mm	52
	Matrix CM Dalbo®-Plus, for ball abutment, incl. lamella retention insert			Vario SR plastic coping, crown, with single-surface antirotation, for angled Vario SR abutments, burn-out	
J2250.0005	Ø 3.3/3.8/4.3/5.0 mm	57	J2563.6003	Ø 5.0 mm	52
	Lamella retention insert, for matrix CM Dalbo®-Plus			Vario SR titanium cap, bridge, incl. Vario SR prosthetic screw	
J2250.0007	Ø 3.3/3.8/4.3/5.0 mm	57	J2564.4301	Ø 3.8/4.3 mm	52
	Driver for Locator®, manual/wrench, L 24.3 mm	35	J2564.6001	Ø 5.0 mm	52
J2253.0001	Locator® Instrument , threepart	35		Vario SR impression cap, closed tray	
J2253.0002	L 83 mm	35	J2565.4300	Ø 3.8/4.3 mm	51
J2253.0003	Locator® Angle measurement guide	35	J2565.6000	Ø 5.0 mm	51
J2253.0004	Locator® Parallel post	35		Vario SR impression cap, open tray, incl. fixing screw	
	Locator® Impression cap		J2566.4300	Ø 3.8/4.3 mm	51
J2253.0200	Ø 3.3/3.8/4.3/5.0 mm	58	J2566.6000	Ø 5.0 mm	51
	Locator® Analog			Vario SR analog	
J2253.0340	Ø 3.3/3.8/4.3/5.0 mm	58	J2567.4300	Ø 3.8/4.3 mm	51
	Locator® Male processing packages		J2567.6000	Ø 5.0 mm	51
J2253.0102	Ø 3.3/3.8/4.3/5.0 mm	59		Vario SR protection cap	
	Locator® Male processing packages for extended range		J2568.4300	Ø 3.8/4.3 mm	52
J2253.0112	Ø 3.8/4.3/5.0 mm	59	J2568.6000	Ø 5.0 mm	52
	Locator® Block out spacer			Vario SR aligning tool 20° , for insertion post, sterile	
J2253.0401	Ø 3.3/3.8/4.3/5.0 mm	59	J2569.0001	Ø 3.8/4.3/5.0 mm	53
	Locator® Processing replacement male		J2569.0002	Ø 3.8/4.3/5.0 mm	53
J2253.0402	Ø 3.3/3.8/4.3/5.0 mm	59		Vario SR base for bar , laser-weldable, incl. Vario SR prosthetic screw	
	Locator® Replacement male , blue		J2570.4300	Ø 3.8/4.3 mm	52
J2253.1002	Ø 3.3/3.8/4.3/5.0 mm	60	J2570.6000	Ø 5.0 mm	52
	Locator® Replacement male , pink			J3709.0015 Universal holder	62
J2253.1003	Ø 3.3/3.8/4.3/5.0 mm	60		Reworking reamer, for base for bar abutment, burn-out, plane surface/cone seat	
	Locator® Replacement male , clear		J3711.0010	Ø 3.3/3.8/4.3/5.0 mm	62
J2253.1005	Ø 3.3/3.8/4.3/5.0 mm	60		Reworking reamer, for base for bar abutment, burn-out, screw seat	
	Locator® Replacement male for extended range, red		J3711.0020	Ø 3.3/3.8/4.3/5.0 mm	62
J2253.2002	Ø 3.8/4.3/5.0 mm	60		Guide System Template drill	
	Locator® Replacement male for extended range, orange		J3713.3300	Ø 3.3 mm	29
J2253.2003	Ø 3.8/4.3/5.0 mm	60	J3713.4300	Ø 3.8/4.3 mm	29
	Locator® Replacement male for extended range, green			Guide System Guiding sleeves	
J2253.2004	Ø 3.8/4.3/5.0 mm	60	J3714.3303	Ø 3.3 mm	28
	Locator® Replacement male for extended range, gray		J3714.3803	Ø 3.8 mm	28
J2253.2000	Ø 3.8/4.3/5.0 mm	60	J3714.4303	Ø 4.3 mm	28
	Logfit® Impression caps			Guide System Seating tool	
J2551.4300	Ø 3.8/4.3 mm	47	J3716.3300	Ø 3.3 mm	29
J2551.6000	Ø 5.0 mm	47	J3716.4300	Ø 3.8/4.3 mm	29
	Logfit® Analogs			Prosthetic screw, for bar abutments, hex, (also for Passive-Fit)	
J2552.4300	Ø 3.8/4.3 mm	47	J4005.1602	Ø 3.3/3.8/4.3/5.0 mm, thread M 1.6 mm	56
J2552.6000	Ø 5.0 mm	47			
	Logfit® Plastic copings, for crowns , burn-out				
J2553.4302	Ø 3.8/4.3 mm	48			
J2553.6002	Ø 5.0 mm	48			
	Logfit® Plastic copings, for bridges , burn-out				
J2553.4301	Ø 3.8/4.3 mm	48			
J2553.6001	Ø 5.0 mm	48			

INDEX

ARTICLE NUMBERS

Art. No.	Article	Page	Art. No.	Article	Page
	Screw, hex, for bar abutment, screw for bases for bar abutments for impression taking, open tray and soldering			Guide System Form drill, SCREW-LINE, Cortical bone	
J4005.1610	Ø 3.3/3.8/4.3/5.0 mm, thread M 1.6 mm	54	J5048.3813	Ø 3.8 mm, L 13 mm	27
	Vario SR prosthetic screw, hex		J5048.3816	Ø 3.8 mm, L 16 mm	27
J4005.2004	Ø 3.8/4.3/5.0 mm, thread M 2.0 mm	53	J5048.4307	Ø 4.3 mm, L 7 mm	27
	Drill extension, ISO shaft		J5048.4309	Ø 4.3 mm, L 9 mm	27
J5002.0005	L 26.6 mm	32	J5048.4311	Ø 4.3 mm, L 11 mm	27
J5002.0010	Adapter, ISO shaft for angled hand piece	33	J5048.4313	Ø 4.3 mm, L 13 mm	27
J5002.0012	Cleaning needle for drills with internal irrigation ...	33	J5048.4316	Ø 4.3 mm, L 16 mm	27
J5002.0020	Cleaning cannula for drills with internal irrigation..	33		Round bur , resterilizable	
	Depth stop SCREW-LINE for pilot drill (J5051.2000) and pre-drill (J5051.2800) with reduced coil		J5050.2300	Ø 2.3 mm	17
J5015.0007	L 7 mm	17	J5050.3500	Ø 3.5 mm	17
J5015.0009	L 9 mm	17		Pilot drill SCREW-LINE , reduced coil, resterilizable	
J5015.0011	L 11 mm	17	J5051.2000	Ø 2.0 mm	17
J5015.0013	L 13 mm	17		Pilot drill without coil, resterilizable	
	Depth stop for form drills SCREW-LINE, with reduced coil		J5051.2003	Ø 2.0 mm	17
J5015.3300	Ø 3.3 mm	18		Pre-drill SCREW-LINE , resterilizable	
J5015.3800	Ø 3.8 mm	18	J5051.2800	Ø 1.7-2.8 mm	17
J5015.4300	Ø 4.3 mm	18		Form drill SCREW-LINE, with reduced coil, resterilizable	
J5015.5000	Ø 5.0 mm	18	J5052.3309	Ø 3.3 mm, L 9 mm	18
	Guide System Gingiva Punches		J5052.3311	Ø 3.3 mm, L 11 mm	18
J5041.3300	Ø 3.3 mm	28	J5052.3313	Ø 3.3 mm, L 13 mm	18
J5041.3800	Ø 3.8 mm	28	J5052.3316	Ø 3.3 mm, L 16 mm	18
J5041.4300	Ø 4.3 mm	28	J5052.3807	Ø 3.8 mm, L 7 mm	18
	Guide System Pilot drill set		J5052.3809	Ø 3.8 mm, L 9 mm	18
J5043.3309	Ø 3.3 mm, L 9 mm	24	J5052.3811	Ø 3.8 mm, L 11 mm	18
J5043.3311	Ø 3.3 mm, L 11 mm	24	J5052.3813	Ø 3.8 mm, L 13 mm	18
J5043.3313	Ø 3.3 mm, L 13 mm	24	J5052.3816	Ø 3.8 mm, L 16 mm	18
J5043.4307	Ø 3.8/4.3 mm, L 7 mm	24	J5052.4307	Ø 4.3 mm, L 7 mm	18
J5043.4309	Ø 3.8/4.3 mm, L 9 mm	24	J5052.4309	Ø 4.3 mm, L 9 mm	18
J5043.4311	Ø 3.8/4.3 mm, L 11 mm	25	J5052.4311	Ø 4.3 mm, L 11 mm	18
J5043.4313	Ø 3.8/4.3 mm, L 13 mm	25	J5052.4313	Ø 4.3 mm, L 13 mm	18
	Guide System Pilot drill		J5052.4316	Ø 4.3 mm, L 16 mm	18
J5044.3316	Ø 3.3 mm, L 16 mm	24	J5052.5007	Ø 5.0 mm, L 7 mm	18
J5044.4316	Ø 3.8/4.3 mm, L 16 mm	25	J5052.5009	Ø 5.0 mm, L 9 mm	18
	Guide System Surgery set, SCREW-LINE		J5052.5011	Ø 5.0 mm, L 11 mm	18
J5045.3309	Ø 3.3 mm, L 9 mm	25	J5052.5013	Ø 5.0 mm, L 13 mm	18
J5045.3311	Ø 3.3 mm, L 11 mm	25	J5052.5016	Ø 5.0 mm, L 16 mm	18
J5045.3313	Ø 3.3 mm, L 13 mm	25		Form drill SCREW-LINE, Cortical bone, resterilizable	
J5045.3807	Ø 3.8 mm, L 7 mm	26	J5053.3316	Ø 3.3 mm	19
J5045.3809	Ø 3.8 mm, L 9 mm	26	J5053.3816	Ø 3.8 mm	19
J5045.3811	Ø 3.8 mm, L 11 mm	26	J5053.4316	Ø 4.3 mm	19
J5045.3813	Ø 3.8 mm, L 13 mm	26	J5053.5016	Ø 5.0 mm	19
J5045.4307	Ø 4.3 mm, L 7 mm	26		Tap SCREW-LINE, with hexagon, resterilizable	
J5045.4309	Ø 4.3 mm, L 9 mm	26	J5054.3309	Ø 3.3 mm	19
J5045.4311	Ø 4.3 mm, L 11 mm	26	J5054.3809	Ø 3.8 mm	19
J5045.4313	Ø 4.3 mm, L 13 mm	26	J5054.4309	Ø 4.3 mm	19
	Guide System Form drill, SCREW-LINE		J5054.5009	Ø 5.0 mm	19
J5046.3316	Ø 3.3 mm, L 16 mm	25		Driver for screw implants, with ISO shaft for angled hand piece	
J5046.3816	Ø 3.8 mm, L 16 mm	26	J5300.0007	L 27.5 mm	32
J5046.4316	Ø 4.3 mm, L 16 mm	26		Driver, short for screw implants, manual/wrench	
	Guide System Form drill, SCREW-LINE, Cortical bone		J5300.0008	L 18.0 mm	32
J5048.3309	Ø 3.3 mm, L 9 mm	27		Driver, long for screw implants, manual/wrench	
J5048.3311	Ø 3.3 mm, L 11 mm	27	J5300.0009	L 27.0 mm	32
J5048.3313	Ø 3.3 mm, L 13 mm	27	J5300.0010	Cardanic driver (30°)	33
J5048.3316	Ø 3.3 mm, L 16 mm	27		Driver for ball abutment, manual/wrench	
J5048.3807	Ø 3.8 mm, L 7 mm	27	J5300.0011	L 18.3 mm	35
J5048.3809	Ø 3.8 mm, L 9 mm	27			
J5048.3811	Ø 3.8 mm, L 11 mm	27			

INDEX

ARTICLE NUMBERS

Art. No.	Article	Page
J5300.0030	PickUp instrument	33
J5300.0057	Surgery Set CAMLOG®/CONELOG® SCREW-LINE	16
	Paralleling pin SCREW-LINE	
J5300.2028	Ø 1.7-2.8/2.0 mm	32
	Guide System Check-up pin	
J5301.3300	Ø 3.3 mm	29
J5301.4300	Ø 3.8/4.3 mm	29
J5302.0010	Holding key for insertion post	33
	Holding sleeve for screw implants	
J5302.3300	Ø 3.3 mm	33
J5302.3800	Ø 3.8 mm	33
J5302.4300	Ø 4.3 mm	33
	Guide System Driver , manual/wrench	
J5303.4300	3.3/3.8/4.3 mm	29
	Guide System Driver , with ISO shaft for angled hand piece	
J5304.4300	3.3/3.8/4.3 mm	29
	Screwdriver Activator , for ball abutment matrix CM Dalbo®-Plus	
J5315.0005	J5316.0501	35
J5316.0501	Screwdriver hex , short, manual/wrench	34
J5316.0502	Screwdriver hex , long, manual/wrench	34
J5316.0503	Screwdriver hex , long, ISO shaft	34
J5316.0504	Screwdriver hex , short, ISO shaft	34
J5316.0510	Screwdriver hex , extra short, manual/wrench.....	34
J5320.1030	Torque wrench	34
	Tap adapter for tap SCREW-LINE	
J5322.0010	L 18.0 mm short.....	32
J5322.0011	L 23.0 mm long.....	32
J5330.8500	Prosthetic tray (without content).....	34
J5330.8600	Prosthetic set	34
	Pre-Osteotome SCREW-LINE , straight convex	
J5417.2800	Ø 1.7-2.8 mm	20, 21
	Osteotome Set CAMLOG®/CONELOG® SCREW-LINE	
J5418.0020	straight convex.....	20
J5420.0020	straight concave.....	22
J5418.0030	angled convex.....	21
J5420.0030	angled concave	23
	Osteotome SCREW-LINE , straight convex	
J5418.3300	Ø 3.3 mm	20
J5418.3800	Ø 3.8 mm	20
J5418.4300	Ø 4.3 mm	20
J5418.5000	Ø 5.0 mm	20
	Osteotome SCREW-LINE , angled convex	
J5418.3310	Ø 3.3 mm	21
J5418.3810	Ø 3.8 mm	21
J5418.4310	Ø 4.3 mm	21
J5418.5010	Ø 5.0 mm	21
	Pre-Osteotome SCREW-LINE , straight concave	
J5419.2800	Ø 1.7-2.8 mm	22, 23
	Osteotome SCREW-LINE , straight concave	
J5420.3300	Ø 3.3 mm	22
J5420.3800	Ø 3.8 mm	22
J5420.4300	Ø 4.3 mm	22
J5420.5000	Ø 5.0 mm	22
	Osteotome SCREW-LINE , angled concave	
J5420.3310	Ø 3.3 mm	23
J5420.3810	Ø 3.8 mm	23
J5420.4310	Ø 4.3 mm	23
J5420.5010	Ø 5.0 mm	23
J8070.2050	Edentulous mandible incl. mounting plate	66

Art. No.	Article	Page
	ALTApin magazine	
M1000.0050	7 titanium pins, sterile, 1 unit	31
M1000.0100	7 titanium pins, sterile, 3 units	31
M5100.0010	ALTApin applicator, straight	31
M5100.0030	ALTApin applicator, angled 90°	31
M5100.0050	ALTApin pricker	31
M5100.0070	ALTApin membrane fixator	31
M5100.0100	ALTApin surgery mallet	31
M5200.0010	ALTApin applicator, straight, work element	31
M5200.0055	ALTApin pricker, insert	31
M5500.0050	ALTApin single patient drill, ISO shaft	31
M5600.0110	ALTApin set	31
M5600.0210	ALTApin tray (without content)	31

Adress/Stamp

Date/Signature

Order Form

Your Customer no.


Delivery

- normal express

We thank you for your order.

CE0123

Art. No. J80010009 Rev. 2 05/2014



HEADQUARTERS

CAMLOG Biotechnologies AG | Margarethenstrasse 38 | CH - 4053 Basel
Telefon +41 61 565 41 00 | Fax +41 61 565 41 01 | info@camlog.com | www.camlog.com

Manufacturer of CAMLOG®/CONELOG® Products: ALTATEC GmbH, Maybachstraße 5, D-71299 Wimsheim

camlog