



CAMLOG® IMPLANT SYSTEM
INFORMATION ABOUT THE
PREVIOUS CAMLOG®
IMPLANT PACKAGING



Supplement for the CAMLOG® Product Catalog 2017

INFORMATION ABOUT THE PREVIOUS CAMLOG® IMPLANT PACKAGING

EXTERIOR PACKAGING

Sealed, color-coded, folding carton with product label

SECONDARY PACKAGING

Sealed

IMPLANT INTERNAL PACKAGING (PRIMARY PACKAGING)

Sealed, color-coded

EXAMPLE PRODUCT LABEL OF EXTERIOR IMPLANT PACKAGING

Article number with diameter and length in mm

Product family

Ø

Length

Do not reuse

Implant type:
SL: SCREW-LINE
RL: ROOT-LINE 2

Implant/surface

camlog K1052.4313 Ø 4,3 L 13 1 St./pc.

SCREW-LINE Implantat, Promote® plus
SCREW-LINE Implant, Promote® plus

UDI

STERILE R

LOT XXXXXXXXXX

Do not reuse

HIBC – Barcode

caution consult accompanying documents !

Caution: US Federal law restricts this device to sale by or on the order of a dentist or physician.

2016-09-05
Use by
2021-09-05

Inhalt ist steril bei unbeschädigter Verpackung
Contents sterile unless packaging is opened or damaged

+E219K105243131/\$\$3210905XXXXXXXXXX/16D20160905D

Rx only

Hersteller / Manufacturer
ALTATEC GmbH
Maybachstr. 5
71299 Wimsheim
GERMANY

CE 0123

Made in Germany

Symbol for sterilization with gamma-radiation

Use-by date

Lot number

Caution, consult accompanying documents

The green dot

Identification of registry in accordance with guideline 93/42 EEC: TÜV SÜD Product Service GmbH

Color-coding according to Implant diameter

INFORMATION ABOUT THE PREVIOUS CAMLOG® IMPLANT PACKAGING

STILL AVAILABLE
TILL JULY
2017

CAMLOG® SCREW-LINE IMPLANT PROMOTE®

	Article	Art.-No.	Ø	L	A Ø
	CAMLOG® SCREW-LINE Implant, Promote® incl. insertion post and cover screw, sterile Material Titanium Grade 4	K1042.3311	3.3 mm	11 mm	2.7 mm
		K1042.3313		13 mm	
		K1042.3316		16 mm	
		K1042.3809	3.8 mm	9 mm	3.5 mm
		K1042.3811		11 mm	
		K1042.3813		13 mm	
		K1042.3816		16 mm	
		K1042.4309	4.3 mm	9 mm	3.9 mm
		K1042.4311		11 mm	
		K1042.4313		13 mm	
		K1042.4316	5.0 mm	16 mm	4.6 mm
		K1042.5009		9 mm	
		K1042.5011		11 mm	
		K1042.5013	6.0 mm	13 mm	5.5 mm
		K1042.5016		16 mm	
		K1042.6009		9 mm	
		K1042.6011	6.0 mm	11 mm	5.5 mm
		K1042.6013		13 mm	
K1042.6016	16 mm				

STILL AVAILABLE
TILL OCTOBER
2017

CAMLOG® SCREW-LINE IMPLANT PROMOTE® PLUS

	Article	Art.-No.	Ø	L	A Ø
	CAMLOG® SCREW-LINE Implant, Promote® plus incl. insertion post and cover screw, sterile Material Titanium Grade 4	K1052.3311	3.3 mm	11 mm	2.7 mm
		K1052.3313		13 mm	
		K1052.3316		16 mm	
		K1052.3809	3.8 mm	9 mm	3.5 mm
		K1052.3811		11 mm	
		K1052.3813		13 mm	
		K1052.3816		16 mm	
		K1052.4309	4.3 mm	9 mm	3.9 mm
		K1052.4311		11 mm	
		K1052.4313		13 mm	
		K1052.4316	5.0 mm	16 mm	4.6 mm
		K1052.5009		9 mm	
		K1052.5011		11 mm	
		K1052.5013	6.0 mm	13 mm	5.5 mm
		K1052.5016		16 mm	
		K1052.6009		9 mm	
		K1052.6011	6.0 mm	11 mm	5.5 mm
		K1052.6013		13 mm	
K1052.6016	16 mm				

CE0123

Art.-Nr. J8001.0016 Rev.03/2017 EN

HEADQUARTERS

CAMLOG Biotechnologies AG | Margarethenstr. 38 | 4053 Basel | Switzerland
Phone +41 61 565 41 00 | Fax +41 61 565 41 01 | info@camlog.com | www.camlog.com

Manufacturer CAMLOG® and CONELOG® Products: ALTATEC GmbH | Maybachstr. 5 | 71299 Wimsheim | Germany

camlog