

A cooperation between Camlog and Amann Girrbach





# CAM blanks type AG

The original type AG CAM blanks were developed jointly by CAMLOG Biotechnologies GmbH and Amann Girrbach AG. They feature the blank holder interface patented by Amann Girrbach and are compatible with the Ceramill<sup>®</sup> CAD/CAM system. The new blanks enable laboratories to fabricate customized one-piece abutments and healing caps made of titanium alloy on in-house milling machines.

# Benefits at a glance

- Original implant connection geometries for maximum precision
- Best possible customer support during the entire workflow
- A well-matched system: compatible design combined with high-precision in-house fabrication on the milling machines **Ceramill® Motion 2, Motion 3 and Matik**
- Color coding of the CAMLOG<sup>®</sup> and CONELOG<sup>®</sup> CAM blanks for easy assignment of the implant diameter



# **Product overview**

#### CAMLOG<sup>®</sup> Implant System

	Art. No. (2 units)	Description	Implant Ø
CAMLOG Ø 43	K2471.3327 K2471.3827 K2471.4327 K2471.5027 K2471.6027	CAMLOG <sup>®</sup> CAM titanium blank, type AG, for the Ceramill <sup>®</sup> CAD/CAM system from Amann Girrbach, Ø 12 mm, titanium alloy. Delivery includes 2 separately packaged abutment screws.	3.3 mm 3.8 mm 4.3 mm 5.0 mm 6.0 mm



CONELOG® Implant System						
	Art. No. (2 units)	Description	Implant Ø			
CONFLOG 8 4 3	C2471.3327 C2471.3827 C2471.4327 C2471.5027	CONELOG® CAM titanium blank, type AG, for the Ceramill® CAD/CAM system from Amann Girrbach, Ø 12 mm, titanium alloy. Delivery includes 2 separately packaged abutment screws.	3.3 mm 3.8 mm 4.3 mm 5.0 mm			



iSy® Implant System					
	Art. No. (2 units)	Description	Implant Ø		
54 154 154	P2471.5027	iSy® CAM titanium blank, type AG, for the Ceramill® CAD/CAM system from Amann Girrbach, Ø 12 mm, titanium alloy. Delivery includes 2 separately packaged abutment screws.	3.8/4.4/5.0 mm		



## Headquarters

CAMLOG Biotechnologies GmbH | Margarethenstr. 38 | 4053 Basel | Switzerland Phone +41 61 565 41 00 | Fax +41 61 565 41 01 | info@camlog.com | www.biohorizonscamlog.com

### **Cooperation partner**

Amann Girrbach AG | Herrschaftswiesen 1 | 6842 Koblach | Austria Telephone +43 5523 62333-200 | austria@amanngirrbach.com | www.amanngirrbach.com

 ${\sf CAMLOG}^{\scriptscriptstyle \oplus}, {\sf CONELOG}^{\scriptscriptstyle \oplus} \text{ and } iSy^{\scriptscriptstyle \oplus} \text{ are registered trademarks of CAMLOG Biotechnologies GmbH.} \\ {\sf Ceramill}^{\scriptscriptstyle \oplus} \text{ and Amann Girrbach}^{\scriptscriptstyle \oplus} \text{ are registered trademarks of Amann Girrbach AG.} \\$ 

