














## SCREW OVERVIEW

For the CONOLOG<sup>®</sup> Implant system



# ABUTMENT AND PROSTHETIC SCREWS INTRAORAL USE






## Implant-Abutment connection

	Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm		
	M 1.6			M 2.0	Tightening torque	
CONELOG® Abutment screw						
 <p>Scanbody ScanPost for Sirona Scanbody</p>					tightening by hand*	
 <p>Temporary Abutments titanium, crown and bridge</p>						
 <p>Esthomic® Abutments</p>		8.9 mm  C4005.1601		8.9 mm  C4005.2001	20 Ncm	
 <p>Universal Abutment Telescope Abutment Gold-plastic Abutment Logfit® Abutment</p>						
 <p>Vario SR Abutments, 20° and 30° angled</p>						
CONELOG® Abutment screw for titanium bases CAD/CAM, dark purple anodized						
 <p>Titanium bases CAD/CAM, crown and bridge</p>		8.9 mm  C4015.1601		8.9 mm  C4015.2001		20 Ncm
CONELOG® Vario SR Abutment screws						
 <p>Vario SR Abutment, straight</p>		10.6 mm  C4007.1600		10.6 mm  C4007.2000	20 Ncm	
CONELOG® Abutment screws with reduced head, light blue anodized						
 <p><b>COMFOUR™</b> Bar Abutments, 17° and 30° angled</p>		7.8 mm  C4004.1601		7.8 mm  C4004.2001	20 Ncm	

\*Optionally für temporary abutments titanium: tightening torque after successful healing phase 20 Ncm













# ABUTMENT AND PROSTHETIC SCREWS INTRAORAL USE

## Abutment-Prosthetic connection

		Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm		
		M 1.6			M 2.0		Tightening torque
Prosthetic screws for bar abutments, light blue anodized							
 <p><b>COMFOUR™</b> Bar Abutments, 17° and 30° angled</p>	3.6 mm	 J4012.1601		3.8 mm	 J4012.2001		15 Ncm
	Vario SR Prosthetic screw, yellow anodized						
 <p>Vario SR Abutments, straight, 20° and 30° angled</p>	4 mm	 J4005.2004					15 Ncm




















# AUXILIARY SCREWS INTRA- AND EXTRAORAL USE

## Abutment-Prosthetic connection

		Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm			
		M 1.6			M 2.0		Tightening torque	
Prosthetic screws for bar abutments, light blue anodized								
 <p>Scanning cap for bar abutments</p>	3.6 mm	 J4012.1601		3.8 mm	 J4012.2001		tightening by hand	
Screws for bar abutments, for impression taking open tray and for soldering, light blue anodized								
 <p><b>COMFOUR™</b> Bar abutments, straight, 17° and 30° angled</p>	12 mm	 J4012.1610		12.2 mm	 J4012.2010		tightening by hand	
	17 mm	 J4012.1615		17.2 mm	 J4012.2015			
	22 mm	 J4012.1620		22.2 mm	 J4012.2020			
	Plastic screws for bar abutment, as fixation and bonding aid, beige							
	29 mm					29.2 mm		tightening by hand
	J4009.1627					J4009.2027		

# LAB SCREWS EXTRAORAL

## Lab analog-Abutment connection












	Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	
	M 1.6		M 2.0		Tightening torque
<b>CONELOG® Lab screws*, brown anodized</b>					
 <p>Scanbody ScanPost for Sirona Scanbody</p>					
 <p>Temporary Abutments titanium, crown and bridge</p>					
 <p>Esthomic® Abutments</p>		8.9 mm  C4006.1601		8.9 mm  C4006.2001	tightening by hand
 <p>Universal Abutment Telescope Abutment Gold-plastic Abutment</p>					
 <p>Vario SR Abutments, 20° and 30° angled</p>					
<b>CONELOG® Lab screws for titanium bases CAD/CAM*, brown partially anodized</b>					
 <p>Titanium bases CAD/CAM, crown and bridge</p>		8.9 mm  C4016.1601		8.9 mm  C4016.2001	tightening by hand
<b>CONELOG® Bonding aids**</b>					
 <p>Titanium bases CAD/CAM, crown and bridge</p>		26 mm 		26 mm 	tightening by hand
<b>CONELOG® Vario SR Lab screws*, brown anodized</b>					
 <p>Vario SR Abutment, straight</p>		10.6 mm  C4008.1600		10.6 mm  C4008.2000	tightening by hand
<b>CONELOG® Lab screws with reduced head*, light blue partially anodized</b>					
 <p><b>COMFOUR™</b> Bar Abutments, 17° and 30° angled</p>		7.8 mm  C4004.1600		7.8 mm  C4004.2000	tightening by hand

\*Note: Lab screws may not be used on patients.

\*\*not separate available, enclosed with the titanium bases CAD/CAM.

# LAB SCREWS EXTRAORAL

## Abutment-Prosthetic connection

		Ø 3.3 mm	Ø 3.8 mm	Ø 4.3 mm	Ø 5.0 mm	
		M 1.6			M 2.0	Tightening torque
Lab prosthetic screws for bar abutments*, brown anodized						
	Scanning cap for bar abutments					
	<b>COMFOUR™</b> Bar Abutments, straight, 17° and 30° angled		3.6 mm  J4013.1601		3.8 mm  J4013.2001	tightening by hand
	Bar lab analog for bar abutments					
Vario SR Prosthetic screw, yellow anodized						
	Vario SR Abutments, straight, 20° and 30° angled			4 mm  J4005.2004		tightening by hand
	Vario SR Analog					
Prosthetic screw for bar abutments*, for fabrication of the wax up on the bar sleeve for titanium bonding base, Passive-Fit, on the bar lab analog						
	Titanium bonding base for bar abutments and bar sleeve for titanium bonding base, burn-out, Passive-Fit		5.5 mm  J4005.1602		5.5 mm  J4005.2002	tightening by hand

\*Note: Lab screws may not be used on patients.



**HEADQUARTERS**

CAMLOG Biotechnologies AG | Margarethenstr. 38 | 4053 Basel | Switzerland  
Phone +41 61 565 41 00 | Fax +41 61 565 41 01 | [info@camlog.com](mailto:info@camlog.com) | [www.camlog.com](http://www.camlog.com)

Manufacturer CAMLOG®/CONOLOG® Products: ALTATEC GmbH, Maybachstr. 5, 71299 Wimsheim, Germany

