


CONELOG[®]
SYSTEM



SYSTEM OVERVIEW



a perfect fit™






camlog

SYSTEM OVERVIEW CONELOG® IMPLANT SYSTEM











SURGERY

			∅ 3.3 mm	∅ 3.8 mm	∅ 4.3 mm	∅ 5.0 mm
IMPLANTS	CONELOG® SCREW-LINE Implants, Promote® plus	Length in mm				
		 7 mm	–	C1062.3807	C1062.4307	C1062.5007
		9 mm	C1062.3309	C1062.3809	C1062.4309	C1062.5009
		11 mm	C1062.3311	C1062.3811	C1062.4311	C1062.5011
		13 mm	C1062.3313	C1062.3813	C1062.4313	C1062.5013
		16 mm	C1062.3316	C1062.3816	C1062.4316	C1062.5016
	Guide System CONELOG® SCREW-LINE Implants, Promote® plus	 7 mm	–	C1063.3807	C1063.4307	–
		9 mm	C1063.3309	C1063.3809	C1063.4309	–
		11 mm	C1063.3311	C1063.3811	C1063.4311	–
		13 mm	C1063.3313	C1063.3813	C1063.4313	–
		16 mm	C1063.3316	C1063.3816	C1063.4316	–
HEALING CAPS	CONELOG® Healing cap, cylindrical	 GH 2.0	C2015.3320	C2015.3820	C2015.4320	C2015.5020
		 GH 4.0	C2015.3340	C2015.3840	C2015.4340	C2015.5040
		GH 6.0	–	C2015.3860	C2015.4360	C2015.5060
	CONELOG® Healing cap, wide body	 GH 4.0	C2014.3340	C2014.3840	C2014.4340	C2014.5040
		 GH 6.0	–	C2014.3860	C2014.4360	C2014.5060
	CONELOG® Healing cap bottleneck	 GH 4.0	C2011.3340	C2011.3840	C2011.4340	C2011.5040
		 GH 6.0	–	C2011.3860	C2011.4360	C2011.5060

PROSTHETICS









IMPLANT IMPRESSION TAKING	CONELOG® Impression post, open tray		C2121.3300	C2121.3800	C2121.4300	C2121.5000
			C2110.3300	C2110.3800	C2110.4300	C2110.5000
	Impression cap for impression post, closed tray		J2111.3300	J2111.3800	J2111.4300	J2111.5000
	CONELOG® Bite registration posts		C2140.3300	C2140.3800	C2140.4300	C2140.5000
FABRICATION OF THE PLASTER MODEL	CONELOG® Lab analogs		C3010.3300	C3010.3800	C3010.4300	C3010.5000

SYSTEM OVERVIEW CONELOG® IMPLANT SYSTEM







			∅ 3.3 mm	∅ 3.8 mm	∅ 4.3 mm	∅ 5.0 mm
CONELOG® ABUTMENTS FOR CROWN AND BRIDGE RESTORATIONS	CONELOG® Temporary abutments, titanium alloy		C2239.3300	C2239.3800	C2239.4300	C2239.5000
	CONELOG® Esthomic® Abutments, straight	 GH 1.5-2.5 3.0-4.5	–	C2226.3815 C2226.3830	C2226.4315 C2226.4330	C2226.5015 C2226.5030
	CONELOG® Esthomic® Abutments, 15° angled, type A	 GH 1.5-2.5 3.0-4.5	–	C2227.3815 C2227.3830	C2227.4315 C2227.4330	C2227.5015 C2227.5030
	CONELOG® Esthomic® Abutments, 15° angled, type B	 GH 1.5-2.5 3.0-4.5	–	C2228.3815 C2228.3830	C2228.4315 C2228.4330	C2228.5015 C2228.5030
	CONELOG® Esthomic® Abutments, 20° angled, type A	 GH 1.5-2.5 3.0-4.5	–	C2231.3815 C2231.3830	C2231.4315 C2231.4330	C2231.5015 C2231.5030
	CONELOG® Esthomic® Abutments, 20° angled, type B	 GH 1.5-2.5 3.0-4.5	–	C2232.3815 C2232.3830	C2232.4315 C2232.4330	C2232.5015 C2232.5030
	CONELOG® Esthomic® Abutments, Inset	 GH 2.0-3.3	C2235.3320	C2235.3820	C2235.4320	C2235.5020
	CONELOG® Universal abutments		C2211.3300	C2211.3800	C2211.4300	C2211.5000
	CONELOG® Gold-plastic abutments		C2246.3300	C2246.3800	C2246.4300	C2246.5000
	CONELOG® Logfit® Abutments	 GH 1.0 GH 2.5	–	C2550.3810 C2550.3825	C2550.4310 C2550.4325	C2550.5010 C2550.5025

GH = gingiva height

SYSTEM OVERVIEW CONELOG® IMPLANT SYSTEM

		∅ 3.3 mm	∅ 3.8 mm	∅ 4.3 mm	∅ 5.0 mm		
CONELOG® ABUTMENTS FOR CROWN AND BRIDGE RESTORATIONS	CONELOG® Telescope abutments, for double crown restorations		–	C2212.3800	C2212.4300	C2212.5000	
CONELOG® CAD/CAM PROSTHETICS FOR CROWNS, BRIDGES AND HYBRID RESTORATIONS	CONELOG® Titanium bases CAD/CAM		GH 0.8	C2242.3308	C2242.3808	C2242.4308	C2242.5008
			GH 2.0	C2242.3320	C2242.3820	C2242.4320	C2242.5020
CONELOG® VARIO SR ABUTMENTS FOR CROWNS, BRIDGES AND HYBRID RESTORATIONS OCCLUSALLY SCREW-RETAINED	CONELOG® Vario SR Abutments, straight		GH 1.0	–	C2560.3808	C2560.4308	C2560.5008
	CONELOG® Vario SR Abutments, 20° angled		GH				
3.5 - 1.9			–	C2561.3800	C2561.4300	–	
4.0 - 1.8			–	–	–	C2561.5000	
	CONELOG® Vario SR Abutments, 30° angled		GH				
3.5 - 1.1			–	C2562.3800	C2562.4300	–	
4.5 - 1.3			–	–	–	C2562.5000	
	Vario SR Impression caps, closed tray (without rotation securing device)		–		J2565.4300	J2565.6000	
	Vario SR Impression caps, open tray (without rotation securing device)		–		J2566.4300	J2566.6000	
	Vario SR analogs		–		J2567.4300	J2567.6000	

SYSTEM OVERVIEW CONELOG® IMPLANT SYSTEM

		∅ 3.3 mm	∅ 3.8 mm	∅ 4.3 mm	∅ 5.0 mm
CONELOG® ABUTMENTS FOR HYBRID RESTORATIONS	CONELOG® Bar abutments				
					
	GH 1.0	C2255.3310	C2255.3810	C2255.4310	C2255.5010
	GH 2.5	C2255.3325	C2255.3825	C2255.4325	C2255.5025
	GH 4.0	–	C2255.3840	C2255.4340	C2255.5040
	CONELOG® Ball abutment sets, incl. male part and matrix				
	CM Dalbo®-Plus, red				
	dublication aid/spacer,				
	stabilizing ring and ball				
	abutment analog				
					
	GH 1.5	C2250.3315	C2250.3815	C2250.4315	C2250.5015
	GH 3.0	C2250.3330	C2250.3830	C2250.4330	C2250.5030
	GH 4.5	–	C2250.3845	C2250.4345	C2250.5045
	CONELOG® Ball abutments, male part				
					
GH 1.5	C2249.3315	C2249.3815	C2249.4315	C2249.5015	
GH 3.0	C2249.3330	C2249.3830	C2249.4330	C2249.5030	
GH 4.5	–	C2249.3845	C2249.4345	C2249.5045	
CONELOG® Locator® Abutments					
					
GH 1.0	C2253.3310	C2253.3810	C2253.4310	C2253.5010	
GH 2.0	C2253.3320	C2253.3820	C2253.4320	C2253.5020	
GH 3.0	C2253.3330	C2253.3830	C2253.4330	C2253.5030	
GH 4.0	C2253.3340	C2253.3840	C2253.4340	C2253.5040	
GH 5.0	–	C2253.3850	C2253.4350	C2253.5050	
CAD/CAM PROSTHETICS	CONELOG® Scanbodies				
		C2600.3310	C2600.4310	C2600.5010	
	CONELOG® ScanPosts for Sirona Scanbody				
		C2620.3306	C2620.3806	C2620.4306	C2620.5006

TORQUE WRENCH SETTINGS

FOR FINAL FIXATION OF SURGICAL COMPONENTS AND SUPRASTRUCTURES



Article

J5320.1030

				Torque			
CONELOG® Implant cover screws				tightened by hand			
CONELOG® Healing caps cylindrical wide body bottleneck							
CONELOG® Impression posts CONELOG® Bite registration posts							
CONELOG® Lab screws		C4006.1601 C4016.1601 (∅ 3.3/3.8/4.3 mm)		C4006.2001 C4016.2001 (∅ 5.0 mm)			
CONELOG® Vario SR lab screws		C4008.1600 (∅ 3.8/4.3 mm)		C4008.2000 (∅ 5.0 mm)			
CONELOG® Temporary abutments (Optional: Torque after finished healing phase 20 Ncm)				20 Ncm			
CONELOG® Abutment screws		C4005.1601 (∅ 3.3/3.8/4.3 mm)				C4005.2001 (∅ 5.0 mm)	
		C4015.1601 (∅ 3.3/3.8/4.3 mm)				C4015.2001 (∅ 5.0 mm)	
CONELOG® Vario SR abutment screws		C4007.1600 (∅ 3.8/4.3 mm)		C4007.2000 (∅ 5.0 mm)			
CONELOG® Esthomic® Abutments, straight CONELOG® Esthomic® Abutments, angled 15°/20°, type A and B CONELOG® Esthomic® Abutments, Inset				tightened by hand			
CONELOG® Gold-plastic abutments CONELOG® Universal abutments CONELOG® Telescope abutments							
CONELOG® Logfit® Abutments CONELOG® Titanium bases CAD/CAM							
(from left to right)							
CONELOG® Vario SR abutments, straight, 20°/30° angled							
Vario SR protection caps							
Vario SR prosthetic screw		J4005.2004 (∅ 3.8/4.3/5.0 mm)					
Vario SR plastic copings/ titanium caps/bases for bar				15 Ncm			

All screws have to be tightened after at least 5 minutes with the corresponding torque!

TORQUE WRENCH SETTINGS

FOR FINAL FIXATION OF SURGICAL COMPONENTS AND SUPRASTRUCTURES



Article

J5320.1030

		Torque	in mm	
CONELOG® Bar abutments		C5300.0020 (∅ 3.3/3.8/4.3/5.0 mm)	20 Ncm	3.3
			30 Ncm	3.8 4.3 5.0
Healing caps for bar abutments		C5300.0027 (∅ 3.3/3.8/4.3/5.0 mm)	tightened by hand	
Prosthetic screw	J4005.1602 (∅ 3.3/3.8/4.3/5.0 mm)	J5316.0510 (hex, extra short) J5316.0501 (hex, short) J5316.0502 (hex, long)	15 Ncm	
Bases for bar abutments Titanium bonding bases for bar abutments				
CONELOG® Ball abutments		J5300.0011	20 Ncm	3.3
			30 Ncm	3.8 4.3 5.0
CONELOG® Locator® Abutments		J2253.0001	20 Ncm	3.3
			30 Ncm	3.8 4.3 5.0
CONELOG® Scanbodies (∅ 3.3/3.8/4.3/5.0 mm)		J5316.0501 (hex, short) J5316.0502 (hex, long)	tightened by hand	
CONELOG® ScanPosts for Sirona Scanbody (∅ 3.3/3.8/4.3/5.0 mm)				

All screws have to be tightened after at least 5 minutes with the corresponding torque!

Not all products are available in all countries.

HEADQUARTERS

CAMLOG Biotechnologies AG | Margarethenstrasse 38 | CH - 4053 Basel
Telefon +41 61 565 41 00 | Fax +41 61 565 41 01 | info@camlog.com | www.camlog.com

Manufacturer of CAMLOG®/CONELOG® products: ALTATEC GmbH, Maybachstraße 5, D-71299 Wimsheim

camlog